SensationTM Short Throw
CaptivatorTM
Captivator IITM
CaptiflexTM

Single-Use Polypectomy Snares

en	Instructions for Use	2
es	Instrucciones de uso	5
fr	Mode d'emploi	8
de	Gebrauchsanweisung	11
it	Istruzioni per l'uso	15
nl	Instructies voor gebruik	18
pt-eu	Instruções de utilização	21

Sensation™ Short Throw Captivator™ Captivator II™ Captiflex™ Single-Use Polypectomy Snares

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

REUSE WARNING

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

DEVICE DESCRIPTION

Contents

(1) Single-Use Polypectomy Snare

Operating Principle

The Single-use Polypectomy Snare consists of a flexible wire cable and loop which can be extended and retracted from the snare's flexible outer sheath using a three-ring handle. The inner diameter of the sheath is lubricant coated to provide minimal friction during extension and retraction of the loop from the sheath. The snare loop is fully opened and placed over the tissue. The snare loop is slowly closed capturing the tissue. A monopolar electrosurgical generator may be used to apply current to aid resection and cauterization. The snare loop is subsequently tightened, while current may be applied, until the tissue is resected. The snare loop can be fully actuated without current to resect tissue as well, known as cold snaring.

The Single-Use Polypectomy Snares adapt to a variety of monopolar generators using the appropriate active cord as indicated below.

The Single-Use Polypectomy Snares are available with the following active cord configuration:

ACTIVE CORD ORDER NUMBER	GENERATOR COMPATIBILITY
M00561250	Olympus Brand Generator
M00561270	Bovie, Valleylab, Aspen Labs, ERBE Brand Generators, EndoStat, EndoStat II and EndoStat III

User Information

A thorough understanding of the technical principles, clinical applications, and risks associated with monopolar and/or cold snaring (non-electrical) polypectomy and tissue resection is necessary before using this product. The Single-use Polypectomy Snare should only be used by or under the supervision of physicians thoroughly trained in endoscopic polypectomy and tissue resection.

Materials

The Single-Use Polypectomy Snare contains the following substance(s) defined as CMR 1B (Carcinogenic, Mutagenic, or Reprotoxic) (Presumed human CMR) in a concentration above 0.1% weight by weight:

Cobalt; CAS No. 7440-48-4; EC No. 231-158-0. This device is made with stainless steel which contains cobalt. Current scientific evidence supports that stainless steels used in medical devices, do not cause an increased risk of cancer or adverse reproductive effects.

INTENDED USE/INDICATIONS FOR USE

The Single-Use Polypectomy Snares are indicated for use endoscopically for the removal and or cauterization of diminutive polyps, sessile polyps, pedunculated polyps and tissue from within the gastrointestinal tract.

Clinical Benefit Statement

Endoscopic resection of pre-cancerous and/or cancerous tissue from within the gastrointestinal tract.

CONTRAINDICATIONS

Contraindications for these devices are those specific to endoscopic Polypectomy and tissue resection.

WARNINGS

Always follow the generator manufacturer's operation suggestions to prevent unnecessary hazard to the operator and/or the patient.

- Modification to this device could result in damage to device and/or patient/user injury.
- When simultaneously using an electrosurgical unit and physiological monitoring equipment you should position any monitoring electrodes as far as possible from the surgical electrodes.
- Consult the neutral electrode manufacturer about proper grounding of the patient. It is recommended that a monitoring neutral electrode be used, if a contact quality monitor is available, or built into the generator. The entire area of the neutral electrode should be attached reliably to the patient's body, and as close to the operating field as possible. The patient should not come into contact with metal parts or objects that may be grounded to earth in order to avoid patient injury. The use of antistatic drapes is recommended in order to insulate the patient.
- Skin-to-skin contact should be avoided (for example between the patient's arms and body) by way of dry cloth or gauze in order to prevent a possible thermal/electrical injury.
- The lines to the surgical electrodes should be fitted so that they do not touch the patient or other lines. Any active electrodes not currently in use should be put to one side so that they cannot touch the patient in order to avoid patient injury.
- It is suggested that the operator and the assistant wear protective gloves to prevent accidental burns. Universal precautions should be used in all cases.

- Flammable agents used for cleaning or disinfecting, or as solvents of adhesives, should be allowed to evaporate before the procedure in order to avoid thermal injuries to patients and/or users.
- Fluids or flammable agents that may pool under the patient or in body depressions or cavities should be mopped prior to using the Snare Device.
- While operating the device, avoid contact with the patient in order to avoid user injuries such as electrical shock or burn.
- This device is not intended to be used in the presence of flammable liquid, in an oxygen enriched atmosphere or in the presence of explosive gases. This may lead to patient or user injuries such as burns.
- Gas embolism may occur as a result of over-insufflation of air during the procedure or during the presence of inert gas prior to high frequency surgery. Endogenous gases should be sucked away, if possible, prior to procedure.
- Leakage current to patient from endoscope, as well as energized polypectomy snare, are additive. Consult the endoscope manufacturer about the proper grounding of the endoscope.
- Monopolar diathermy or electrosurgical cautery in patients with pacemakers or implantable cardiac defibrillators can result in electrical reset of the cardiac device, inappropriate sensing and/or therapy, tissue damage around the implanted electrodes, or permanent damage to the pulse generator. A cardiologist should be consulted prior to using Snare in these patients.
- Should be used with caution and only after careful consideration in patients who are at risk for bleeding complications.

PRECAUTIONS

- Consult the ESU manufacturer manual to ensure proper settings and recommendations by the manufacturer are followed.
- The Single-use Polypectomy Snare must be used in conjunction with a Type BF or CF generator; see Generator Compatibility section. The generator connector (Active Cord sold separately) is connected to the snare handle by a plug pushed onto the connector as far as possible so none of the connecting pin is visible. The opposite end of the active cord is inserted into the generator.
- The output power setting selected should remain under 50 Watts, and be as low as possible for the intended purpose.
- It is very important that if the proper setting of the generator is not known, one should set the unit at a power setting lower than the recommended range and cautiously increase the power until the desired effect is achieved. It is recommended that the operator not use the device with any generator setting which may output a voltage exceeding the Maximum Accessory Voltage Rating. Maximum Accessory Voltage Rating: 3800 V peak (7600 V peak-to-peak). Active accessories should be selected that have an Accessory Voltage Rating equal to or greater than 3800 V peak. Consult the ESU manual for recommended monopolar settings.

Precaution: The active length information provided is intended to be used with the guidance from the Electrosurgical Generator manufacturer to determine if the overall length of an Electrosurgical Accessory remains within the maximum allowable length.

The overall active length of each Electrosurgical Accessory is the following:

- Captiflex (M00562470) and Captivator II (M00561220): 251 cm.
- Captiflex (M00562420), Captivator (M00562300 and M00562450) and Captivator II (M00561230): 252 cm.

- Captivator II (M00561240 and M00561190): 253 cm.
- Captiflex (M00562400), Captivator (M00561310, M00562320, M00562340 and M00562370) and Captivator II (M00561290): 255 cm.
- Captivator (M00562390): 256 cm.
- Sensation (M00560320 and M00562690): 262.1 cm.
- Sensation (M00560300, M00560310, M00560330, M00562650 and M00562670): 264.3 cm.

The overall active length includes the connector, cord handle and loop. The length of other connecting cords must be considered when used with the snare.

ADVERSE EVENTS

Complications include, but may not be limited to:

- Fulgurations
- Hemorrhage
- Infection
- Pain
- Perforation
- Stricture Formation
- Tissue Damage
- Transmural burn
- Transient ileus

HOW SUPPLIED

Single-Use Polypectomy Snares are supplied sterile using an ethylene oxide(EO) process and are packaged in boxes containing multiple units with one (1) unit per sealed pouch.

Do not use if package is damaged or unintentionally opened before use.

Do not use if labeling is incomplete or illegible.

Handling and Storage

This product has no special handling or storage requirements. Check device label for expiration date.

OPERATIONAL INSTRUCTIONS

Preparation for Use

1. Open the snare package and remove it from the pouch.

Precaution: To avoid damage, uncoil the snare before proceeding with preparation.

2. Visually inspect the snare for kinks, loop fraying, and overall product integrity. If any damage is found, do not use the snare and return to Boston Scientific for replacement.

Precaution: Please review the operations and service manuals of the electrosurgical generator for proper setup and operation prior to using the Snare.

Warning: To avoid patient injury, ensure that the patient is properly grounded prior to use of the monopolar electrosurgical generator and snare.

- 3. Test the snare handle assembly, prior to passing it through the endoscope, by sliding the finger grip forward and backward several times. The snare loop should completely extend from and retract into the tip of the sheath.
- 4. The snare is now ready for use.

Directions

- 1. Before inserting the snare into endoscope, fully retract the loop into the sheath.
- 2. Prior to polypectomy or tissue resection, identify and position the endoscope proximal to the intended site. Repeat irrigation or suctioning may be necessary to clear visual field facilitating accurate snare loop placement.
- 3. If intending to apply monopolar electrical current: Connect the appropriate active cord to the connector on the snare handle by completely inserting it until no pin is visible. Insert the other end of the active cord into the appropriate monopolar electrosurgical generator. If not intending to apply monopolar electrical current: Connection to an active cord is not required.

Warning: The electrosurgical generator should be placed in the OFF position prior to advancing the snare through the endoscope to avoid injury to patient or equipment resulting from improper electrical circuit grounding.

Precaution: Advance the snare through the endoscope using short, deliberate 2-3 cm (20-30 mm) movements to prevent inadvertent damage to the snare or endoscope.

- 4. Advance the retracted snare through the endoscope until visible endoscopically and follow standard techniques for performing polypectomy or tissue resection. If intending to apply monopolar electrical current: It is recommended that the operator not use the device with any generator setting which may output a voltage exceeding the Maximum Accessory Voltage Rating. Maximum Accessory Voltage Rating: 3800 V peak (7600 V peak-to-peak). Active accessories should be selected that have an Accessory Voltage Rating equal to or greater than 3800 V peak.

Warning: Check for the proper position of the snare loop using direct vision. Positioning the snare loop in an improper location may lead to patient injury.

Warning: Polypectomy and tissue resection should be performed under direct endoscopic vision at all times.

Improper use could result in patient injury. Determining the approximate retraction of the snare loop. Polypectomy and tissue resection should be performed under direct endoscopic vision at all times.

Sheath Withdrawal

- 1. If using monopolar electrical current: Prior to removal of the snare from the endoscope, ensure the electrosurgical generator is turned OFF.

Warning: Before removing the snare from the endoscope, fully retract the loop into the sheath. Remove the snare from the endoscope slowly and deliberately to prevent inadvertent damage to the endoscope and or patient injury.

Disposal

To minimize the risk of infection or microbial hazards after use, dispose device and packaging as follows:

After use, device may contain biohazardous substances. The device and packaging should be treated and disposed of as biohazardous waste or have them treated and

disposed of in accordance with any applicable hospital, administrative, and/or local government regulations. Use of a biohazardous container with biological hazard symbol is recommended. Untreated biohazardous waste should not be disposed of in the municipal waste system.

Post-Procedure

Any serious incident that occurs in relation to this device should be reported to the manufacturer and relevant local regulatory authority.

PATIENT COUNSELING INFORMATION

Brief the patient on any relevant post-procedure instructions, contraindications, warnings, precautions and/or adverse events found within this Instructions for Use (IFU), pertaining to the patient.

WARRANTY

For device warranty information, visit (www.bostonscientific.com/warranty).

Captivator, Captiflex and Sensation are trademarks of Boston Scientific Corporation or its affiliates. All other trademarks are the property of their respective owners.

SYMBOL DEFINITIONS

Commonly used medical device symbols that appear on the labeling are defined at www.bostonscientific.com/SymbolsGlossary.

Table 1. Symbol Definitions

	Contents		Minimum Required Working Channel
---	----------	---	----------------------------------

Sensation™ Short Throw Captivator™ Captivator II™ Captiflex™

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All rights reserved.
Boston Scientific (Australia) Pty Ltd
Boston Scientific Limited
© 2024 Boston Scientific Corporation or its affiliates.