Alliance™ II

Integrated Inflation/Litho Handle

Directions for Use	2
Instrucciones de uso	5
Mode d'emploi	8
Gebrauchsanweisung	11
Istruzioni per l'uso	14
Gebruiksaanwijzing	17
Instruções de Utilização	20

51105965-01 Rev. A

2020-03

Alliance™ II

Integrated Inflation/Litho Handle

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

DEVICE DESCRIPTION

The Alliance II Integrated Inflation/Litho Handle is a nonsterile, reusable tool used with the sterile, single use Alliance II Syringe/Gauge Assembly (sold separately), designed to simplify inflation and deflation of balloon dilatation catheters. The Alliance II Integrated Inflation/Litho Handle is also used with the sterile, single-use Trapezoid™ RX Wireguided Retrieval Basket, which is supplied separately. When assembled, the system is intended to aid in the crushing of gallstones during endoscopic retrograde cholangiopancreatography.

The diagram shown below (Figure 1) describes the key components of the Alliance II Integrated Inflation/Litho Handle when used with the Syringe/Gauge Assembly. On the left side of the handle is a blue directional switch that controls the movement of the Syringe/Gauge Assembly plunger or the movement of the basket. Pumping the handle lever inflates and deflates the balloon catheter or extends and retracts the basket to crush the captured stone.


Figure 1.

Contents

(1) Alliance II Integrated Inflation/ Litho Handle

SINGLE-USE SYRINGE/GAUGE ASSEMBLY

The Alliance Syringe/Gauge Assembly is a disposable syringe pre-assembled with a disposable pressure gauge. The gauge is graduated in both atmospheres (12 ATM maximum) and kilopascals (1216 kPa maximum). The Syringe/Gauge Assembly is to be used with the reusable Alliance II Integrated Inflation/Litho Handle when performing balloon dilatation.

INTENDED USE / INDICATIONS FOR USE

The Alliance II Integrated Inflation/Litho Handle is indicated for inflation and deflation of balloon dilatation catheters while simultaneously monitoring balloon catheter pressures and for the crushing of biliary calculi when used with the Trapezoid RX Wireguided Retrieval Basket.

Clinical Benefit Statement

To provide inflation and deflation of balloon catheters and to aid in the crushing of calculi during ERCP procedures.

ADVERSE EVENTS

None known

HOW SUPPLIED

The Alliance II Integrated Inflation/Litho Handle is packaged with one non-sterile Alliance II Integrated Inflation/Litho Handle.

Do not use if package is damaged or unintentionally opened before use. Do not use if labeling is incomplete or illegible.

Handling and Storage

This product has no special handling or storage requirements.

Precaution: A thorough understanding of the technical principles, clinical applications, and risks associated with balloon dilatation and/or ERCP and mechanical lithotripsy is necessary before using this product.

Inspect the Alliance II Integrated Inflation / Litho Handle for damage and utility. Pump the handle lever to ensure the device's ability to move in the direction indicated by the arrows. The Syringe/Gauge Assembly component of the Alliance II Integrated Inflation/Litho System is supplied ETO sterilized in a sealed tray. Carefully examine the unit to verify that neither the contents nor the sterile package

has been damaged during shipment. **DO NOT USE if damaged.**
Immediately return the damaged product to Boston Scientific.

INSTRUCTIONS FOR USE FOR DILATATION

1. To insert the Alliance II Syringe/Gauge Assembly into the Alliance II Integrated Inflation/Litho Handle, turn the blue directional switch to the upward position. Align the syringe flanges with the slot in the handle and align the back of the syringe plunger with the slot in the heel plate. Push down on the Syringe/Gauge Assembly firmly to insure proper fit. Snap the heel clip over the syringe plunger. Move the heel plate and syringe plunger all the way forward.

Note: Improper Syringe/Gauge Assembly alignment may cause inability to properly inflate the balloon.

Optional: Extension tubing is provided for your added convenience. If using the extension tubing, connect to the Syringe/Gauge Assembly and the balloon catheter hub using the Luer-Lok®.

2. Turn the blue directional switch backward toward the red arrow. Pump the handle lever to fill the Syringe/Gauge Assembly with 35 ml of sterile water for balloons up to a 20 mm diameter.

Note: NEVER use air or any gas medium to inflate the balloon catheters. Maximum dilatation results are achieved using fluids.

3. Purge excess air from the Syringe/Gauge Assembly by turning the blue directional switch forward toward the green arrow, tilting the assembly so that the air bubbles float towards the Luer-Lok end and pumping the handle lever until all air is purged from the syringe.
4. Attach the balloon catheter hub to the Syringe/Gauge Assembly or to the extension tubing via the Luer-Lok.
5. Pull a vacuum on the balloon catheter by moving the blue directional switch backward toward the red arrow and pump the handle lever until it stops.
6. Remove the protective sleeve from the balloon.

Optional: Apply silicone spray to the balloon to facilitate passage through the endoscope.

Note: DO NOT pre-inflate or pre-test the balloon.

7. To inflate the balloon, move the blue directional switch forward toward the green arrow, holding the Alliance II Integrated Inflation System at a downward angle to prevent air from entering the balloon.
8. Pump the handle lever until the gauge on the Syringe/Gauge Assembly registers pressure.

Note: Refer to the balloon package label for the maximum rated balloon pressure. The gauge will read pressure in both atmospheres (ATMs) and kilopascals (kPa).

9. Continue inflation until the desired pressure is achieved.

PRECAUTION: DO NOT exceed the maximum balloon pressure specified on the balloon.

10. To deflate, turn the blue directional switch backward toward the red arrow and pump the handle lever until it stops.

Note: The blue directional switch may be harder to turn when the Syringe/Gauge Assembly is inflated to high pressures.

PRECAUTION: Slide arm may retract rapidly when used with balloons rated above 5 ATM (507 kPa).

11. Return the blue directional switch to the upward position to detach the balloon catheter.

DIRECTIONS FOR USE WITH LITHOTRIPSY

1. To attach the Trapezoid RX Wireguided Retrieval Basket to the Alliance II Integrated Inflation/Litho Handle, turn the blue directional switch to the upward position.
2. Align the thumb ring and the handle ring with the heel plate and body portion of the Alliance II Integrated Inflation/Litho Handle respectively.
3. Push the basket handle down into the Alliance II Integrated Inflation/Litho Handle until secure.

Note: The heel clip is not intended for use with the Trapezoid RX Wireguided Retrieval Basket.

4. To retract the basket, turn the blue directional switch forward toward the green arrow. Firmly pump the handle lever until the calculus is fragmented.

Note: Once the calculus is fragmented, stop squeezing the handle of the Alliance II device to avoid disengaging the Trapezoid basket tip.

Note: Basket retraction and extension may be viewed under fluoroscopy.

5. Stop pumping the handle lever when the basket fully retracts into the catheter tip.
6. To extend the basket, turn the blue directional switch backward toward the red arrow. Firmly pump the handle lever until the basket is fully extended.
7. To detach the lithotripsy basket from the Alliance II Integrated Inflation/Litho Handle, return the blue directional switch to the upward position, and pull upward on the lithotripsy basket handle.

Disposal

After each use, clean and disinfect all external and accessible surfaces using the cleaning and disinfection instructions included in the user manual.

When the device reaches the end of its useful life, to minimize the risk of infection or microbial hazards, dispose of device and packaging as follows:

Device may contain biohazardous substances. The device and packaging should be treated and disposed of as biohazardous waste. They should be treated and disposed of in accordance with any applicable hospital, administrative, and/or local government regulations. Use of a biohazard container with biological hazard symbol is recommended. Untreated biohazardous waste should not be disposed of in a municipal waste system.

Post-Procedure

Any serious incident that occurs in relation to this device should be reported to the manufacturer and relevant local regulatory authority.

CLEANING INSTRUCTIONS

DO NOT immerse the unit in disinfecting solutions. DO NOT steam autoclave.

The Alliance II Integrated Inflation/Litho Handle can be wiped with clean water and disinfecting solution such as 2% glutaraldehyde. Follow the solution manufacturer's recommended instructions for safe handling.

PRECAUTION: Ensure that the Alliance II Integrated Inflation/Litho Handle is dry prior to use, otherwise slippage may occur.

WARRANTY

For device warranty information, visit (www.bostonscientific.com/warranty).

EU Importer: Boston Scientific International B.V., Veststraat 6, 6468 EX Kerkrade, The Netherlands

Alliance and Trapezoid are trademarks of Boston Scientific Corporation or its affiliates.

All other trademarks are the property of their respective owners.

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All rights reserved.
© 2020 Boston Scientific Corporation or its affiliates.