HurricaneTM RX

BILINARY

Balloon Dilatation Catheter

en	Instructions for Use	2
es	Instrucciones de uso	8
fr	Mode d'emploi	14
de	Gebrauchsanweisung	20
it	Istruzioni per l'uso	27
nl	Instructies voor gebruik	34
pt-eu	Instruções de Utilização	41

51173927-01

2024-12

TABLE OF CONTENTS

REUSE WARNING3

DEVICE DESCRIPTION3

 Contents.....3

 Operating Principle and Device Detail3

 Balloon Construction.....3

 Catheter Construction3

 Materials4

 User Information4

INTENDED USE/INDICATIONS FOR USE4

 Clinical Benefit Statement.....4

CONTRAINDICATIONS.....4

WARNINGS.....4

PRECAUTIONS.....4

ADVERSE EVENTS4

HOW SUPPLIED.....5

 Device Details.....5

 Handling and Storage5

OPERATIONAL INSTRUCTIONS.....5

 Inspection and Preparation5

 Catheter Insertion5

 Balloon Inflation6

 Catheter Withdrawal6

 Disposal7

 Post-Procedure7

PATIENT COUNSELING INFORMATION7

WARRANTY.....7

SYMBOL DEFINITIONS.....7

 Table 1. Symbol Definitions7

Hurricane™ RX

BILIARY

Balloon Dilatation Catheter

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

REUSE WARNING

For single use only. Do not reuse, reprocess or sterilize. Reuse, reprocessing or sterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or sterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

DEVICE DESCRIPTION

Contents

Balloon dilatation catheter

(1) Stopcock

(1) FloSwitch

(1) Preloaded stylet

Operating Principle and Device Detail

The Hurricane RX Biliary Balloon Dilatation Catheter is a double lumen catheter with a balloon mounted at the distal tip. Dilatation balloon catheters are used to exert radial force to dilate narrow duct segments, as well as the Sphincter of Oddi.

Balloon Construction

Balloons are constructed with specially treated polymers and processes, which provide maximum strength. Each balloon inflates to its stated diameter and length at a specific pressure. Maximum pressure is dependent on the balloon size and is identified on the label.

Warning: Do not exceed the maximum pressure printed on the label as this can lead to balloon burst and/or patient injuries, such as perforation, hemorrhage and/or tissue damage.

Catheter Construction

The catheter is multi-lumen. The guidewire lumen begins at the distal tip of the catheter, travels under the black sleeve where it exits at the proximal end of the black sleeve, which is 25 cm from the catheter tip. This lumen is used to pass the catheter over a 0.035 in (0.89 mm) guidewire to provide for the rapid exchange feature. The lumen marked “**Balloon**” is the balloon inflation lumen. The catheter shaft tapers beneath the balloon segment to achieve the lowest possible deflated profile. The lumen marked “**Injection**”, at the proximal end of the catheter, can be used for contrast medium injection for fluoroscopic visualization of the bile duct. There is a FloSwitch Adapter connected to the hub of the “**Injection**” lumen. The FloSwitch Adapter is to remain closed until distal dye contrast injections are required.

Two radiopaque markers are placed under the balloon segment of the catheter to provide visual reference points for balloon positioning within the stricture or duct. A stainless steel stiffening stylet is inserted in the proximal portion of the catheter shaft to augment catheter pushability through the endoscope and should be left in place. This stylet ends under the black tunnel.

Materials

Contains cobalt: CAS No. 7440-48-4; EN No. 231-158-0.
Defined as a 1B carcinogen and reproductive toxicant according to the European Commission in a concentration above 0.1% weight by weight.

Note: This device is made with stainless steel which contains cobalt. Current scientific evidence supports that stainless steels used in medical devices, do not cause an increased risk of cancer or adverse reproductive effects.

User Information

A thorough understanding of the technical principles, clinical applications and risks associated with balloon dilatation of the biliary tract is necessary before using this product.

INTENDED USE/INDICATIONS FOR USE

The Hurricane RX Biliary Balloon Dilatation Catheter is recommended for endoscopic dilatation of strictures of the biliary tree and the Sphincter of Oddi.

The Hurricane RX Biliary Balloon Dilatation Catheter may be used for injection of contrast medium for fluoroscopic visualization of the bile ducts.

Clinical Benefit Statement

The clinical benefit of the Hurricane RX Biliary Balloon Dilatation Catheter is to facilitate endoscopic dilatation of strictures of the biliary tree and sphincter of Oddi and enables distal contrast injections for fluoroscopic visualization of bile ducts.

CONTRAINDICATIONS

None known.

WARNINGS

- To prevent balloon burst, do not exceed the Maximum Pressure printed on the label as this can lead to balloon burst and/or patient injuries, such as perforation, hemorrhage and/or tissue damage.
- To ensure minimum balloon profile, do not pre-inflate, pretest balloon, or attempt to refold balloon into protective sleeve as this could lead to patient injury.

PRECAUTIONS

Any use of this product for procedures other than those indicated in these instructions is not recommended.

- Monitor pressure using an inflation device and pressure gauge. As dilatation takes place, the pressure reading may fluctuate. Adjust the balloon inflation pressure, as necessary to maintain desired pressure.
- If the balloon does rupture or a significant loss of pressure within the balloon occurs, deflate the balloon and carefully remove the balloon and endoscope together as a unit. Do not attempt to withdraw a ruptured balloon through the endoscope.

ADVERSE EVENTS

The adverse events that may result from a balloon dilatation procedure include but may not be limited to:

- Allergic reaction
- Hemorrhage
- Infection/sepsis
- Inflammation
- Pain and tenderness
- Perforation
- Tissue damage

HOW SUPPLIED

Device Details

The Hurricane RX Biliary Balloon Dilatation Catheter is supplied sterile using an ethylene oxide (EO) process.

Do not use if package is damaged or unintentionally opened before use.

Do not use if labeling is incomplete or illegible. If damage is found, call your Boston Scientific representative. See product label for expiration date.

Handling and Storage

This product has no special handling or storage requirements.

OPERATIONAL INSTRUCTIONS

Inspection and Preparation

Precaution: The Hurricane RX Biliary Balloon Dilatation Catheter is packaged sterile in a pouch. Before using, inspect the pouch for any breach of the package to ensure a sterile product.

1. Prior to use, carefully examine the unit to verify that the catheter or sterile package has not been damaged in shipment.
2. Remove plastic mandrel from distal end of catheter.

Precaution: Do not remove stylet from proximal injection hub as the stylet is not intended to be removed from the device. Attempting to remove the stylet may cause device damage.

3. Close FloSwitch Adapter attached to injection hub.
4. **Prepare balloon inflation syringe with contrast medium*.**
5. Attach balloon inflation hub to inflation syringe with a gauge to monitor pressure.
6. To maximize endoscope passage, apply a vacuum to the catheter before removing protective sleeve from the balloon section of the catheter. Close stopcock, disconnect syringe from stopcock, remove air from syringe and reconnect to balloon hub.
7. Remove protective sleeve from balloon prior to insertion.
8. **For Contrast Injection:** Prime second syringe (10 ml (cc) recommended) with contrast medium* and attach to FloSwitch Adapter, open FloSwitch Adapter and flush if needed.

*** Note:** Commercially available contrast fluids have different viscosities that may affect balloon inflation and deflation performance as well as distal contrast injection performance. The higher the viscosity (thicker) the slower the fluid will move through the catheter. Dilute the contrast fluid with saline as required to improve inflate-deflate time as well as distal contrast injection performance.

Warning: To ensure minimum balloon profile, do not pre-inflate, pretest balloon, or attempt to refold balloon into protective sleeve as this could lead to device damage and/or patient injury.

Catheter Insertion

The Hurricane RX Biliary Balloon Dilatation Catheter is designed to pass through a duodenoscope with a minimum 3.2 mm working channel. Each catheter's working length is 180 cm. The balloon should be advanced over the guidewire and through the scope channel using short, deliberate 2-3 cm movements to prevent inadvertent kinking of the catheter shaft. Fluoroscopy can be used to confirm proper placement.

1. The Hurricane RX Biliary Balloon Dilatation Catheters are to be used with a Boston Scientific RX locking device.
2. Ensure the RX locking device is firmly seated onto the endoscope against the side of the working channel port.

3. With the 0.035 in (0.89 mm) guidewire positioned and locked into place, backload the Hurricane RX Biliary Balloon Dilatation Catheter over the guidewire. Slide the Hurricane RX Biliary Balloon Dilatation Catheter onto the guidewire and through the first 25 cm until the guidewire exits the proximal end of the black sleeve. Grasp the guidewire and slide the catheter over the guidewire until it reaches the RX locking device.
4. Unlock the guidewire from guidewire locking device. Pass the distal 25 cm section of the catheter past the guidewire locking device while maintaining tension on the guidewire. Lock the guidewire into place and continue advancing the Hurricane RX Biliary Balloon Dilatation Catheter through the scope until cannulation is achieved.
5. Two radiopaque markers located under the dilatation section of the balloon will aid in positioning the catheter while using fluoroscopy.
6. Contrast medium can be injected through the FloSwitch Injection Port to visualize the bile duct if needed. The FloSwitch Adapter must be opened prior to injection. The guidewire tip may be pulled back into the catheter to help facilitate contrast injection.

Note: Use caution if retracting wire for aid in contrast injection. Retracting the guidewire tip more than 20 cm into the catheter tip may facilitate loss of guidewire cannulation through the Hurricane RX Biliary Balloon Dilatation Catheter.

Balloon Inflation

Under fluoroscopic guidance, the balloon should be advanced and positioned across the stricture to be dilated. The balloon can be inflated by means of a hand-held syringe.

Warning: To prevent balloon burst and patient injury, do not exceed the inflation pressure given for the diameter on the catheter's package label. If the balloon does rupture or a significant loss of pressure within the balloon occurs, deflate the balloon completely and carefully remove the balloon and endoscope together as a unit. Do not attempt to withdraw a ruptured balloon through the endoscope. Continue procedure with a new catheter.

It is recommended that a 10 ml (cc) or larger syringe with a gauge to monitor balloon pressure be used when inflating a Hurricane RX Biliary Balloon Dilatation Catheter. The typical inflation medium is a 1:1 mixture of contrast medium and normal saline.

Precaution: Do not use air or a gas medium to inflate the balloon as this will result in reduced balloon effectiveness.

Catheter Withdrawal

1. When dilatation has been completed, deflate the balloon by applying suction to the balloon lumen.
2. Do not pull back on the catheter until the balloon is completely deflated.

Warning: The balloon must be thoroughly deflated, and all fluid removed prior to withdrawal in order to prevent patient injury and/or damage to the device (up to 30 seconds depending on inflation medium).

3. Slowly remove the catheter from the endoscope.

Warning: If excessive resistance is felt, remove the endoscope and balloon catheter together as a complete unit to prevent damage to body tissue, the catheter, or the endoscope.

4. The catheter should be disposed of in accordance with accepted hospital guidelines.

Note: The larger the syringe diameter, the greater the suction that is applied. For maximum deflation, a 50 ml (cc) syringe is recommended.

Disposal

To minimize the risk of infection or microbial hazards after use, dispose device and packaging as follows:

After use, device may contain biohazardous substances. The device and packaging should be treated and disposed of as biohazardous waste or have them treated and disposed of in accordance with any applicable hospital, administrative, and/or local government regulations. Use of a biohazardous container with biological hazard symbol is recommended. Untreated biohazardous waste should not to be disposed of in the municipal waste system.

Post-Procedure

Any serious incident that occurs in relation to this device should be reported to the manufacturer and relevant local regulatory authority.

PATIENT COUNSELING INFORMATION

Brief the patient on any relevant post-procedure instructions, contraindications, warnings, precautions and/or adverse events found within this Instructions for Use (IFU), pertaining to the patient.

WARRANTY

For device warranty information, visit (www.bostonscientific.com/warranty).

Hurricane is a trademark of Boston Scientific Corporation or its affiliates. All other trademarks are the property of their respective owners.

SYMBOL DEFINITIONS

Commonly used medical device symbols that appear on the labeling are defined at www.bostonscientific.com/SymbolsGlossary.

Table 1. Symbol Definitions

	Contents		Minimum Required Working Channel
	Maximum Pressure		Recommended Guidewire
	Includes Stopcock		Remove & Discard

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or its affiliates. All rights reserved.
Boston Scientific Limited
Boston Scientific (Australia) Pty Ltd
© 2024 Boston Scientific Corporation