# Rotablator™ or ROTAPRO™

Rotational Atherectomy System

# ROTAWIRE™ Drive and wireClip™ Torquer

Guidewire and Guidewire Manipulation Device

|       |                          |    |
|-------|--------------------------|----|
| en    | Instructions for Use     | 2  |
| es    | Instrucciones de uso     | 8  |
| fr    | Mode d'emploi            | 14 |
| de    | Gebrauchsanweisung       | 20 |
| it    | Istruzioni per l'uso     | 26 |
| nl    | Instructies voor gebruik | 32 |
| pt-EU | Instruções de Utilização | 38 |

TABLE OF CONTENTS

**REUSE WARNING .....3**

**DEVICE DESCRIPTION .....3**

    Contents .....3

    Guidewire .....3

        Table 1.....3

    wireClip Torquer.....3

        Figure 1. wireClip Torquer.....3

    Patient Population.....3

    Non-Pyrogenic.....3

    User Information.....3

**INTENDED USE/INDICATIONS FOR USE .....3**

**CLINICAL BENEFIT STATEMENT.....4**

**SUMMARY OF SAFETY AND CLINICAL PERFORMANCE.....4**

**CONTRAINDICATIONS.....4**

**WARNINGS.....4**

**PRECAUTIONS.....4**

    Figure 2.....5

    Figure 3.....5

    Figure 4.....5

    Figure 5.....5

**ADVERSE EVENTS .....5**

**CLINICAL STUDIES .....6**

    Rotablator System Multi-Center Registry Safety Results.....6

        Table 2. Rotablator System Multi-Center Registry - Safety Results at 1 Year

        Follow-up.....6

**HOW SUPPLIED .....6**

    Device Details .....6

    Handling and Storage.....6

**OPERATIONAL INSTRUCTIONS .....6**

    Instructions for use.....6

        Figure 6.....6

    Disposal .....7

    Post-Procedure .....7

**INFORMATION TO BRIEF THE PATIENT .....7**

**WARRANTY.....7**

**SYMBOL DEFINITIONS.....7**

# Rotablator™ or ROTAPRO™

Rotational Atherectomy System

## ROTAWIRE™ Drive and wireClip™ Torquer

Guidewire and Guidewire Manipulation Device

### Rx ONLY

**Caution:** Federal Law (USA) restricts this device to sale by or on the order of a physician.

#### REUSE WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, contact your Boston Scientific Representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death.

Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

#### DEVICE DESCRIPTION

##### Contents

##### Quantity    Material

One (1)    ROTAWIRE Drive Guidewire

One (1)    wireClip Torquer

##### Guidewire

The ROTAWIRE Drive Guidewires are available for use with the Rotablator and ROTAPRO Systems. Models available are Floppy and Extra Support ROTAWIRE Drive Guidewires.

These guidewires are all 0.009 in (0.23 mm) in diameter with an enlarged distal spring tip of 0.014 in (0.36 mm) diameter. The overall length of the guidewire is 330 cm. The guidewires are differentiated on the basis of the spring tip length and the stiffness of the shaft proximal to the spring. A comparison of the guidewires is shown in Table 1.

Table 1.

| Wire Type                    | Shaft Characteristic                                       | Spring Tip Length |
|------------------------------|------------------------------------------------------------|-------------------|
| ROTAWIRE Drive Floppy        | Long taper, minimum guidewire bias (most flexible)         | 2.2 cm            |
| ROTAWIRE Drive Extra Support | Short taper, stiffer, and increased guidewire bias (stiff) | 2.8 cm            |

The spring tip configuration is atraumatic, radiopaque, and the straight tip can be shaped to form a steerable system. The wire shaft is constructed of stainless steel with a smooth finish. These guidewires, designed exclusively for the Rotablator and ROTAPRO Systems, can be independently advanced and steered. The ROTAWIRE Drive is lubricated with Hystrene, which is applied to the entire length of the wire other than the spring tip. The purpose of Hystrene is to assist with the performance of the Rotablator or ROTAPRO Systems.

#### wireClip Torquer

The wireClip Torquer is a plastic device that attaches to guidewires that have shaft diameters from 0.009 in (0.23 mm) to 0.018 in (0.46 mm). The wireClip Torquer provides a convenient gripping surface for manipulating guidewires (reference Figure 1).


Figure 1. wireClip Torquer

#### Patient Population

Rotational atherectomy is indicated in patients who have:

- Single vessel atherosclerotic coronary artery disease with a stenosis.
- Multiple vessel coronary artery disease that in the physician's judgement does not pose undue risk to the patient.
- Certain patients who have had prior percutaneous transluminal coronary angioplasty (PTCA), and who have a restenosis of the native vessel; or,
- Native vessel atherosclerotic coronary artery disease that is less than 25 mm in length.

#### Non-Pyrogenic

The ROTAWIRE Drive Guidewire is non-pyrogenic.

#### User Information

Intended users of the ROTAWIRE Drive Guidewire are primarily interventional cardiologists experienced in performing percutaneous coronary intervention in the sterile field. Cath lab technicians and/or nurses experienced in assisting percutaneous rotational coronary angioplasty are also intended users for preparation and assistance.

#### INTENDED USE/INDICATIONS FOR USE

These guidewires are intended for use with the Rotablator or ROTAPRO Rotational Atherectomy System.

## CLINICAL BENEFIT STATEMENT

The ROTAWIRE Drive Guidewire facilitates rotational atherectomy procedures. Rotational atherectomy is a technique used during percutaneous coronary interventions, which have the clinical benefit of improving myocardial perfusion.

## SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

For customers in the European Union, use the device name found in the labeling to search for the device's Summary of Safety and Clinical Performance, which is available on the European database on medical devices (EUDAMED) website: (<https://ec.europa.eu/tools/eudamed>).

## CONTRAINDICATIONS

Carefully read this document and refer to the Rotablator and ROTAPRO System Console "Operator's Manual", RotaLink™ Exchangeable Burr Catheter "Instructions for Use", and ROTAPRO "Instructions for Use" observing all Contraindications, Restrictions, Warnings, and Precautions for specific information on the use of these components.

## WARNINGS

When advancing or removing the guidewire, always use fluoroscopic guidance with radiographic equipment that provides high resolution images.

Guidewire manipulation within the vasculature should be done under fluoroscopic visualization.

Since the guidewire functions as a monorail that the advancer/catheter/burr tracks over, it is imperative that you initially place the guidewire in the stenotic lumen or the virtual lumen of the vessel and not in a false channel.

Do not allow the individual burr run time to exceed 30 seconds as this may lead to wire fracture/tip separation that may result in perforation, dissection, embolism, myocardial infarction and in rare cases, death. The ROTAWIRE Drive Guidewire has an expected functional life of 5 minutes (total of individual burr run times).

The use of a ROTAWIRE Drive Guidewire may straighten a tortuous vessel such that it places the point of burr-vessel interaction on the lesser curvature of the vessel (burr bias) potentially inducing vasospasm and pseudostenoses that may result in perforation and/or dissection during rotational atherectomy treatment.

Care needs to be taken to maintain coaxial alignment of the guide catheter and ROTAWIRE Drive Guidewire/burr assembly during ablation. Failure to do so may cause a transected ROTAWIRE Drive Guidewire that may result in embolism, myocardial infarction, dissection, and /or surgical intervention and in rare cases, death.

Do not torque, advance or withdraw guidewire if significant resistance is felt. Torqueing, advancing, or withdrawing guidewire under resistance may cause guidewire and/or vessel damage.

Exercise care in handling of the ROTAWIRE Drive Guidewire during the procedure to reduce the possibility of accidental breakage, bending, kinking, coil separation, or loop making in the aorta. A tight loop, kink or sharp bend (greater than 90 degrees in the guidewire) may cause fracture during use. Resulting wire fracture may require additional percutaneous intervention or surgery.

Never advance the rotating burr by advancing the sheath. Guidewire buckling may occur, and perforation or vascular trauma may result. Always advance the rotating burr by using the advancer knob.

Never advance the rotating burr to the point of contact with the guidewire spring tip, as this may result in distal detachment and embolization of the tip. The distal spring tip is designed to not pass through the lumen of the burr catheter. Do not pass through the lumen of the burr catheter.

Never operate the guidewire brake release unless you have a firm grip on the guidewire using the wireClip Torquer. Releasing the brake without first securing the guidewire may result in rotation and entanglement of the guidewire.

Do not allow the burr to remain in one location while rotating at high speeds, as this may lead to wear of the guidewire. Gently advance or retract the burr while it is in high-speed rotary motion. In instances when long ablation runs are required, particularly in calcified, angulated lesions, reposition the guidewire to expose a previously unused segment or exchange the guidewire to prevent damage.

## PRECAUTIONS

The guidewire should not be advanced within vessels or stenoses with a narrower diameter than the spring tip to avoid guidewire and vessel damage.

Withdrawing the device with too much force could cause the spring tip to stretch or separate. The distal spring tip is designed to not pass through the lumen of the burr catheter. The ROTAWIRE Drive cannot be withdrawn through the advancer. The advancer must be removed prior to ROTAWIRE Drive exchange.

Handle the guidewire carefully to prevent a tight loop, kink, or sharp bend ( $> 90^\circ$ ) from forming in the guidewire, which may cause it to fracture during use. Treating certain types and/or locations of lesions or patients with certain conditions is inherently riskier, regardless of the therapeutic device being used.

Physicians should be aware of the higher risk when treating patients, such as:

1. Patients who are not candidates for coronary artery bypass.
2. Patients with severe, diffuse three-vessel disease (multiple diseased vessels should be treated in separate sessions).
3. Patients with unprotected left main coronary artery disease.
4. Patients with ejection fraction less than 30 %.
5. Lesions longer than 25 mm.
6. Angulated ( $\geq 45^\circ$ ) lesions.
7. Patients for whom anticoagulation is not indicated.


When performing percutaneous rotational atherectomy with the ROTAPRO System on-site surgical backup should be included as a clinical consideration.

Do not proximally remove guidewire from carrier tube as it may result in kinked or bent guidewires or tips.

Do not operate the Rotablator or ROTAPRO System if there is a bend, kink, or loop in the guidewire or if the spring tip is prolapsed.


If a loop forms in the guidewire, never pull to straighten it. To remove loop, follow these steps:

**Step 1.** Pull back slightly on the guide catheter until the catheter is out of the ostium (reference Figure 2).


**Figure 2.**

**Step 2.** Using the wireClip Torquer, rotate the guidewire one-half turn clockwise (reference Figure 3).


**Figure 3.**

**Step 3.** If the loop still exists, rotate the wire one-half turn in the opposite direction (counterclockwise) and reevaluate (reference Figure 4).

**Figure 4.**

**Step 4.** With the loop removed, reseat the guide catheter while pulling back slightly on the guidewire (reference Figure 5).

**Figure 5.**

If patient defibrillation is necessary, the burr should be pulled back from the treatment area while the guidewire stays in position.

Before starting defibrillation, be sure that the guidewire is sufficiently insulated from the ground.

Contact of the guidewire with the ground must be avoided to ensure a type CF (cardiac floating) insulation for the patient.

Be sure that the proximal end of the guidewire is at all times located on the sterile covering cloth or on similar insulating protective material.

#### **ADVERSE EVENTS**

Potential adverse events which may result from the use of ROTAWIRE Drive Guidewires, Rotablator, and ROTAPRO Systems include but are not limited to:

- Additional intervention or surgery
- Allergic reaction (contrast, device or other)
- Angina
- Arrhythmia
- Cerebral vascular accident/Transient ischemic attack
- Death
- Electric shock
- Embolism (air, device, plaque, thrombus, tissue, or other)
- Hemorrhage or hematoma
- Hypotension/hypertension (including cardiogenic shock)
- Infection/sepsis
- Myocardial ischemia or infarction
- Pericardial effusion/cardiac tamponade
- Pulmonary edema
- Radiation injury
- Thrombosis/thrombus
- Vasospasm
- Vessel occlusion (slow flow, no reflow, abrupt vessel closure)
- Vessel trauma (arteriovenous fistula, dissection, perforation, pseudoaneurysm, rupture or injury)

There may also be complications associated with distortion, kinks, and fracture of the guidewire and physical deterioration or malfunction of the device, which can lead to patient injury or death.

CLINICAL STUDIES

Rotablator System Multi-Center Registry Safety Results

The Rotablator System Multi-Center Registry evaluates the safety and efficacy of the Rotablator System when used as a stand-alone device or with adjunctive PTCA.

The following data was compiled at the 1-year follow-up from the Multi-Center Registry (2,736 procedures and 3,424 lesions at 22 centers). (Reference Table 2).

Table 2. Rotablator System Multi-Center Registry - Safety Results at 1 Year Follow-up

| Major Complications             | Frequency | Angiographic Complications | Frequency |
|---------------------------------|-----------|----------------------------|-----------|
| Death                           | 1.0 %     | Intimal Dissection         | 13.7 %    |
| Q-Wave Myocardial Infarction    | 1.1 %     | Acute Vessel Closure       | 5.1 %     |
| Coronary Artery Bypass Grafting | 2.5 %     | Vessel Perforation or Tear | 0.7 %     |
|                                 |           | Arrhythmia                 | 2.7 %     |

HOW SUPPLIED

Contents supplied STERILE using an ethylene oxide (EO) process in a peelable pouch. The wireClip Torquer is packaged on a holding card in the same sterile pouch. Do not use if sterile barrier is damaged. If damage is found, contact your Boston Scientific Representative.

Device Details

Do not use if package is damaged or unintentionally opened before use.

Do not use if labeling is incomplete or illegible.

Handling and Storage

This product has no special handling or storage requirements.

Use the device prior to the "Use by" date noted on the product label.

OPERATIONAL INSTRUCTIONS

Instructions for use

Select the appropriate guidewire for the procedure; then follow these steps:

- Using sterile technique, open the pouch and extract the packaging coil that contains the guidewire.
- Unload the guidewire from the packaging coil as follows:
  - Locate the proximal wire retainer on the inside diameter of the coil. Carefully remove the wire from the retainer. This will expose the proximal end of the guidewire.
  - Locate the guard tubing on the outside of the coil. Grasp the tubing and pull back gently to release the end from the small distal anchor tube.
  - Remove the guard tubing, sliding it forward off the coil. This will expose the distal end of the guidewire and the spring tip.
  - Grasp the exposed distal guidewire near the end of the coil tube and gently pull it out of the package. Care should be taken to avoid grasping the spring tip.

**Note:** The guidewire may be unloaded directly from the package into the guide catheter (reference Figure 6). Inspect the guidewire for damage. If damaged, do not use.


Figure 6.

- The guidewire is lubricated with a thin film of lubricant which may appear as a white powder. Do not wipe off lubricant. The lubricant can cause the guidewire to adhere to the inside of the tube, making unloading difficult. If this happens, gently tap the outside of the coil to loosen the wire. Use care not to stretch or damage the spring tip.
- Inspect the spring tip region for damage. Gently shape the spring tip of the guidewire. If damaged, do not use. Handle the guidewire carefully to prevent a tight loop, kink, or sharp bend (> 90°) from forming. These will make passage of the wire through the advancer difficult and may cause it to fracture during use.
- Using standard sterile angioplasty procedure, and under fluoroscopic guidance; gently advance the guidewire past the lesion, using a bare-wire or free-wire technique. The spring tip must be placed distal to the lesion.
- Grasp the proximal tip of the guidewire and thread this end into the hole in the tip of the burr. Continue feeding the wire into the Rotablator or ROTAPRO System catheter until it appears at the back of the advancer, then grasp and gently pull the exposed wire until the burr is a few centimeters from the hemostasis valve.
- If you have difficulty guiding the wire through the advancer, slide the advancer knob back and forth while gently pushing the wire. This usually helps to ease the wire through the advancer. Remove any lubricant that may have built up on the burr by gently wiping it with a gloved fingertip.
- The advancer has an internal guidewire brake that is automatically applied when compressed gas is supplied from the console. This brake prevents the guidewire from rotating during operation of the advancer.
- When using the RotaLink Burr or ROTAPRO Advancer, a wireClip Torquer should always be placed on the guidewire. Attach a wireClip Torquer to the guidewire at a point a few centimeters behind the end of the advancer. To attach the wireClip Torquer: squeeze the two handles to open the clip jaws, place the clip adjacent to the guidewire, then move the clip until the wire is entirely engaged in the groove. Release the handles to allow the wireClip Torquer to securely grasp the wire. This will facilitate steering and advancement of the wire, if the wire is not already placed distal to the lesion.

10. Manipulate the guidewire with the wireClip Torquer by firmly gripping only the cylindrical portion of the wireClip Torquer. The clip can be repositioned as often as necessary. Always ensure that the wire is entirely engaged in the groove of the clip.
11. To test the internal guidewire brake, grasp the wireClip Torquer while the Rotablator or ROTAPRO System is running and attempt to retract the guidewire at the point where it exits from the back of the advancer. During normal operation, the wire is securely gripped by the internal automatic brake and resists any attempts at rotation or advancement. In some cases, however, you may want to release the brake in order to better steer the guidewire or to exchange the advancer. When using the brake release, always ensure that the wireClip Torquer is attached properly and hold the wireClip Torquer firmly to prevent the guidewire from rotating. Secure the guidewire in place while tracking the interventional device over the guidewire and across the lesion.
12. Once the procedure is completed, first remove the wireClip Torquer from the ROTAWIRE Drive Guidewire before removing the RotaLink Advancer or ROTAPRO Advancer and ROTAWIRE Drive Guidewire assembly.  
If the guidewire spring appears to be unraveling during removal, stop the removal procedure. Carefully place a balloon catheter or exchange catheter over the guidewire, positioning the device as distal as possible and remove as a unit.

---

**Note:** This wire cannot be exchanged through the ROTAPRO or Rotablator systems. Secure the guidewire in place while tracking the advancer over the guidewire and across the lesion.

---

If you use a balloon catheter, briefly inflate the balloon as necessary to relieve any spasm. When the spasm stops, continue to gently remove the guidewire.

## Disposal

To minimize the risk of infection or microbial hazards after use, dispose device and packaging as follows:

After use, device and packaging may contain biohazardous substances. Any device and packaging that came into contact with biohazardous substances should be treated and disposed of as biohazardous waste or be treated and disposed of in accordance with any applicable hospital, administrative, and/or local government regulations. Use of a biohazardous container with biological hazard symbol is recommended. Untreated biohazardous waste should not be disposed of in the municipal waste system.

## Post-Procedure

Assess patient for hematoma and/or other signs of bleeding at the access site.

- Any serious incident that occurs in relation to this device should be reported to the manufacturer and relevant local regulatory authority.

## INFORMATION TO BRIEF THE PATIENT

The physician should consider the following points while counseling the patient on the use of the ROTAWIRE Drive in association with the interventional procedure:

- Discuss the risks and benefits, including review of potential adverse events listed in this document, both for the ROTAWIRE Drive along with the Rotablator or ROTAPRO System, and for other interventional treatments likely to be employed.
- Discuss post-procedure instructions, including any lifestyle changes, medications, and home care and rehabilitation guidelines.

No separate patient information is provided because the ROTAWIRE Drive along with the Rotablator or ROTAPRO System, is used in association with other interventional procedures. Risks and benefits to the patient are similar to and, from the patient's perspective, a part of those interventional procedures.

## WARRANTY

For device warranty information, visit ([www.bostonscientific.com/warranty](http://www.bostonscientific.com/warranty)).

## SYMBOL DEFINITIONS

Commonly used medical device symbols that appear on the labeling are defined at [www.bostonscientific.com/SymbolsGlossary](http://www.bostonscientific.com/SymbolsGlossary).

Additional symbols are defined at the end of this document.

ROTAPRO, RotaLink, Rotablator, Dynaglide, wireClip, and ROTAWIRE are trademarks of Boston Scientific Corporation or its affiliates.

All other trademarks are the property of their respective owners.

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All rights reserved.
Boston Scientific Limited
Boston Scientific (Australia) Pty Ltd
© 2022 Boston Scientific Corporation or its affiliates.