Coyote™ ES

OVER - THE - WIRE

PTA Balloon Dilatation Catheter

en	Instructions for Use	2
es	Instrucciones de uso	6
fr	Mode d'emploi	10
de	Gebrauchsanweisung	14
it	Istruzioni per l'uso	19
nl	Instructies voor gebruik	24
pt-EU	Instruções de Utilização	29

TABLE OF CONTENTS

REUSE WARNING STATEMENT 3

DEVICE DESCRIPTION 3

 Contents 3

 Operating Principle..... 3

 Materials..... 3

 Non-Pyrogenic..... 3

 User Information..... 3

INTENDED USE/INDICATIONS FOR USE 3

 Clinical Benefit Statement 3

CONTRAINDICATIONS..... 3

WARNINGS..... 3

PRECAUTIONS..... 3

ADVERSE EVENTS 4

HOW SUPPLIED..... 4

 Device Details 4

 Handling and Storage..... 4

OPERATIONAL INSTRUCTIONS..... 4

 Additional Items for Safe Use 4

 Preparation..... 4

 Inspection Prior to Use 4

 Patient Preparation..... 4

 Inflation Device Preparation..... 4

 Dilatation Catheter Selection..... 4

 Dilatation Catheter Preparation..... 4

 Inflation Device Connection to Catheter 4

 Procedure 4

 Use of the Coyote ES OTW PTA Balloon Dilatation Catheter..... 4

 Table 1. Typical Coyote ES OTW PTA Balloon Dilatation Catheter Compliance 5

 Dilatation Catheter Exchange Procedure..... 5

 Disposal..... 5

 Post-Procedure 5

INFORMATION TO BRIEF THE PATIENT..... 5

WARRANTY..... 5

SYMBOL DEFINITIONS..... 5

Coyote™ ES

OVER-THE-WIRE

PTA Balloon Dilatation Catheter

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

REUSE WARNING STATEMENT

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

DEVICE DESCRIPTION

The Coyote ES Over-The-Wire Percutaneous Transluminal Angioplasty (PTA) Balloon Dilatation Catheter is a high performance balloon catheter for peripheral indications. The device features an ultra low profile, semi-compliant balloon combined with a low profile tip. The catheter is compatible with 0.014 in (0.39 mm) guidewires.

The proximal portion of the catheter includes one female luer-lock port connected to the inflation lumen, and one female luer-lock port for guidewire lumen.

Hydrophilic (Bioslide) coating is applied from the distal tip to just proximal of the balloon and hydrophobic (Xtra) coating is applied from the distal tip to just proximal of the balloon.

Contents

Quantity Material

1 Coyote ES Over-The-Wire PTA Balloon Dilatation Catheter

Operating Principle

The Coyote ES Over-The-Wire PTA Balloon Dilatation Catheter is an Over-The-Wire (OTW) catheter with a semi-compliant balloon fixed at the distal tip. The balloon catheter has a coaxial shaft design. The outer lumen is used for inflation of the balloon, and the wire lumen permits the use of guidewires 0.014 in (0.39 mm) to facilitate advancement of the catheter to and through the stenosis to be dilated. The balloon is designed to provide an inflatable segment of known diameter and length at recommended pressures. The catheter includes a tapered tip to facilitate advancement of the catheter to and through the stenosis. The working length of the balloon catheter is approximately 142 cm. The catheters are available in two versions (1.5 mm version and 2.0 mm through 4.0 mm) that provide different performance characteristics for treating different types of stenoses. Marks on the proximal portion of the catheter shaft at 90 cm and 100 cm indicate the exit of the dilatation catheter tip out of the guide catheter.

Radiopaque marker bands, in conjunction with fluoroscopy, aid in the placement of the catheter's balloon segment. The 1.5 mm model has one radiopaque marker band located at the mid point of the balloon, while all other models have two radiopaque marker bands (one proximal and one distal).

Materials

The Coyote ES OTW PTA Balloon Dilatation Catheter consists of a polyamide outer shaft, and a polyamide semicompliant balloon with a Hydrophilic (Bioslide) coating applied from the distal tip to just proximal of the balloon to reduce friction when tracking the catheter through tortuous vasculature and a hydrophobic (Xtra) coating applied from the distal tip to just proximal of the balloon to minimize the potential for auto-adhesion. The proximal shaft is a polyamide sheathed stainless steel hypotube.

Non-Pyrogenic

This device meets pyrogen limit specifications.

User Information

Only physicians, technicians, and nurses experienced in preparing for and performing peripheral endovascular procedures should use this device.

INTENDED USE/INDICATIONS FOR USE

The Coyote ES OTW PTA Balloon Dilatation Catheter is indicated for Percutaneous Transluminal Angioplasty (PTA) in the peripheral vasculature, including iliac, femoral, popliteal, infra-popliteal and renal arteries, and for the treatment of obstructive lesions of native or synthetic arteriovenous dialysis fistulae. The 2.00 mm - 4.00 mm balloon devices are also indicated for post-dilatation of balloon expandable and self-expanding stents in the peripheral vasculature.

Clinical Benefit Statement

The PTA Balloon Dilatation Catheters are intended to be used to treat obstructive lesions in the peripheral vasculature. The clinical benefit can be measured by overall clinical outcomes as demonstrated by technical success and a safety profile comparable to other existing therapies.

CONTRAINDICATIONS

None known.

WARNINGS

- To reduce the potential for vessel trauma, the inflated diameter of the balloon should approximate, or be less than, the diameter of the vessel just proximal and distal to the stenosis.
- When the catheter is exposed to the vascular system, it should be manipulated while under high-quality fluoroscopic observation. Do not advance or retract the catheter unless the balloon is fully deflated under vacuum. If resistance is felt during manipulation, determine the cause of the resistance before proceeding.
- Inflation in excess of the rated burst pressure may cause the balloon to rupture. In-vivo balloon pressures shall never exceed the rated burst pressure (refer to Table 1 for rated burst pressures).
- Use only the appropriate balloon inflation medium (typically a 50/50 mixture by volume of contrast medium and sterile saline). Never use air or any gaseous medium to inflate the balloon.
- Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

PRECAUTIONS

- Care should be taken to control the position of the guide catheter tip during manipulation of the balloon catheter.
- If a distal protection device is used, follow the manufacturer's instructions for use. The maximum distance from the distal edge of the most distal radiopaque marker to the distal end of this catheter is 9 mm.
- Coyote ES OTW PTA Balloon Dilatation Catheters have soft tips that can become engaged with the distal protection device if advanced beyond the distance recommended by the manufacturer.
- Carefully inspect the catheter prior to use to verify that catheter has not been damaged during shipment and that its size, shape, and condition are suitable for the procedure for which it is to be used.
- Use the catheter prior to the "Use By" date specified on the package.
- Bench testing was conducted with the Coyote ES OTW PTA Dilatation Catheters and marketed Boston Scientific balloon expandable and self-expanding stents. Consideration should be taken when this device is used with different manufacturers' stents due to differences in stent design.
- The Coyote ES OTW PTA Balloon Dilatation Catheter should be used with caution for procedures involving calcified lesions or synthetic vascular grafts due to the abrasive nature of these lesions.
- The Coyote ES OTW PTA Balloon Dilatation Catheters are not intended for injection of contrast medium.
- Non dilatable lesions are infrequently encountered in the population, do not exceed RATED BURST PRESSURE (RBP). To prevent over pressurization, use of a pressure monitoring device is recommended.
- If resistance is felt during post procedure withdrawal of the catheter, it is recommended to extract the entire system with the guide catheter.
- Short-term and long-term biological effects at pressures above the nominal pressure are not known.
- Balloon dilatation procedures should be performed only under fluoroscopic observation with radiographic equipment that provides high resolution images.
- Do not advance any portion of the dilatation catheter system against significant resistance. The cause of resistance should be determined via fluoroscopy before proceeding.
- Never advance the angioplasty catheter without the guidewire extending from the tip.
- Do not pull the Balloon Protector proximally on to the catheter shaft.
- If unusual resistance is felt during removal of the balloon protector and mandrel, do not use the catheter and replace with another. Follow the product returns procedure for unused product.
- To minimize the possible introduction of air into the system, it is imperative that prior to proceeding, careful attention is paid to the maintenance of tight catheter connections and thorough aspiration and flushing of the system.
- Any use for procedures other than those indicated in these instructions is not recommended.
- Precautions to prevent or reduce thrombosis/thrombus should be taken when any catheter is used:
 - Consider systemic anticoagulation.

- Flush or rinse all products entering the vascular system with sterile isotonic saline or a similar solution prior to use.
- Do not store catheters where they are directly exposed to organic solvents or ionizing radiation.

ADVERSE EVENTS

Potential adverse events which may be associated with the use of a peripheral balloon when used as indicated include, but are not limited to:

- Allergic reaction (to contrast, device, or other)
- Death
- Embolism (air, device, plaque, thrombus, tissue, or other)
- Hematoma
- Hemorrhage
- Infection/Sepsis
- Ischemia/Necrosis
- Need for additional intervention or surgery
- Renal Insufficiency or Renal Failure
- Restenosis
- Thrombus/Thrombosis
- Vasospasm
- Vessel occlusion
- Vessel Trauma (AV fistula, dissection, perforation, pseudoaneurysm, rupture, or injury)

HOW SUPPLIED

The Coyote ES OTW PTA Balloon Dilatation Catheters are supplied STERILE using an ethylene oxide (EO) process.

Device Details

One Coyote ES OTW PTA Balloon Dilatation Catheter is loaded into a coiled carrier tube which is placed inside a sterile barrier pouch and packed within a shelf carton.

Do not use if package is damaged or unintentionally opened before use.

Do not use if labeling is incomplete or illegible.

Handling and Storage

This product has no special handling or storage requirements.

OPERATIONAL INSTRUCTIONS

Additional Items for Safe Use

Materials Required for PTA with the Coyote ES OTW PTA Balloon Dilatation Catheter:

Quantity	Material
1 or 2	Guidewire(s) of appropriate size for advancement of guide sheath or guide catheter
1	Appropriate arterial sheath and dilator set (for femoral approach only) (per Coyote ES OTW PTA Balloon Dilatation Catheter device labeling)
1 or 2	Guiding catheters in the appropriate size and configuration to select the targeted artery (per Coyote ES OTW PTA Balloon Dilatation Catheter device labeling)
1	Vial of contrast medium
1	Vial of sterile saline solution
1	Inflation device with manometer
1 or 2	Coyote ES OTW PTA Balloon Dilatation Catheter(s)
1 or 2	0.014 in (0.39 mm) guidewire(s) that are of appropriate length relative to the procedure requirements
1	Luer-lock syringe
1	Hemostatic adapter
1	Three-way stopcock

Preparation

Inspection Prior to Use

Prior to use, carefully examine all equipment to be used during the procedure, including the catheter, to verify proper function. Verify that the catheter and sterile packaging have not been damaged. Verify that the catheter size is suitable for the specific procedure for which it is intended.

Check the pouch for "Use By" date. Carefully inspect the sterile pouch before opening.

Do not use the product after the "Use By" date. If the integrity of the sterile package has been compromised prior to the product "Use By" date (e.g., damage of the package), contact your local Boston Scientific representative for return information. Do not use if any defects are noted.

Note: Do not continue to use the catheter if damage occurs or sterility is compromised during use.

Patient Preparation

Patient preparation and sterile precautions should be the same as for any endovascular procedure.

Inflation Device Preparation

1. Prepare the inflation device according to the manufacturer's instructions.
2. Purge the system of air.

Dilatation Catheter Selection

The inflation diameter of the balloon must not exceed the diameter of the targeted artery proximal and distal to the stenosis. Similarly, the inflated length of the balloon (shoulder to shoulder) should not exceed the length of the stenosis within the targeted artery. Careful consideration should also be given to select an appropriate catheter shaft length that considers the relationship between the arterial access point and the location of the targeted lesion or stenosis.

Dilatation Catheter Preparation

1. The catheter is packaged in a protective coil. Remove catheter from carrier tube. Use care when removing the catheter to avoid damage.
2. Remove the balloon protector and mandrel by grasping the catheter just proximal to the balloon (at the proximal balloon bond site). With the other hand, gently grasp the proximal section of the balloon protector and remove distally. For Coyote ES Over-The-Wire Catheters, the mandrel will slide off with the balloon protector.

Precaution: If unusual resistance is felt during removal of the balloon protector or mandrel, do not use the catheter and replace with another. Follow the product returns procedure for unused product.

3. Prepare the balloon dilatation catheter for purging. Fill a 10 ml, 12 ml, or 20 ml luer-lock syringe with 3 ml of contrast medium. Use only the appropriate balloon inflation medium (e.g., the equivalent of a 50 % contrast medium/50 % sterile saline solution).

Precaution: Use only the appropriate balloon inflation medium (typically a 50/50 mixture by volume of contrast medium and sterile saline). Never use air or any gaseous medium to inflate the balloon.

4. Connect a three-way stopcock to the port fitting on the dilatation catheter. Flush through the stopcock.
5. Connect the syringe to the stopcock. Use care when connecting the catheter to avoid damage.
6. Hold the syringe with the nozzle pointing downward and aspirate for 15 seconds - 20 seconds. Release the plunger.
7. Remove the syringe and evacuate all air from the barrel.
8. To prevent the possibility of air embolization, reconnect the syringe and aspirate for 15 seconds - 20 seconds until bubbles no longer appear during aspiration. Repeat process two more times. If air bubbles persist, discard the device. Release plunger and set aside until "Inflation Device Connection to Catheter" step 2.
9. Prepare the wire lumen of the catheter by attaching a syringe to the wire lumen port and flushing the lumen with approximately 5 ml sterile saline solution.
10. Immerse catheter in a sterile saline bath.

Inflation Device Connection to Catheter

1. To remove any air lodged in the distal luer fitting of the inflation device, purge approximately 1 ml of contrast medium.
2. Attach inflation device/syringe to stopcock; attach to inflation port.
3. With tip down, orient system vertically.
4. Open stopcock to balloon catheter; apply negative pressure for 15 seconds - 20 seconds.
5. Close stopcock to balloon catheter; purge inflation device/syringe of all air.
6. Repeat steps 3 - 5 until all air is expelled. If bubbles persist, do not use device.
7. If a syringe was used, attach a prepared inflation device to stopcock.
8. Open stopcock to balloon catheter.

Procedure

Use of the Coyote ES OTW PTA Balloon Dilatation Catheter

1. Insert introducer sheath and guide sheath/guide catheter using standard techniques. The choice of guiding sheaths or catheters depends upon anatomy and lesion location. Before insertion of the dilatation catheter, administer Heparin as appropriate.
2. Insert a guidewire through the hemostatic valve following the manufacturer's instructions or standard practice. Advance the guidewire carefully into the guiding sheath/catheter. When complete, withdraw the guidewire introducer, if used.
3. Attach a torque device to the wire, if desired. Under fluoroscopy, advance the guidewire to the desired vessel, then across the stenosis.
4. Thoroughly aspirate and flush the guiding sheath/catheter in preparation for introduction of the dilatation catheter.

- Backload the distal tip of the dilatation catheter onto the guidewire ensuring that the guidewire exits the wire port of the Coyote ES OTW PTA Balloon Dilatation Catheter manifold. When loading or exchanging the catheter, it is recommended to thoroughly wipe the guidewire clean for better catheter movement on the guidewire.

Note: To avoid kinking, advance the catheter slowly, in small increments, until the proximal end of the guidewire emerges from the catheter.

- Advance the dilatation catheter through the hemostatic valve slowly, while the balloon is fully deflated. If resistance is encountered, do not advance the dilatation catheter through the adapter. Care should be taken not to over-tighten the hemostatic adapter around the dilatation catheter shaft as lumen constriction may occur, affecting inflation/deflation of the balloon.
- Connect the side port of the guiding sheath/catheter hemostatic adapter to the proximal pressure recording/infusion line or manifold assembly, which permits proximal pressure recording or infusion through the guiding sheath/catheter.
- Position the balloon relative to the lesion to be dilated and inflate the balloon to the appropriate pressure (reference table 1). If difficulty is experienced during balloon inflation, do not continue; remove the catheter. It is strongly recommended that negative pressure is maintained on the balloon between inflations. Repeat inflation of balloon (not to exceed 10 total inflations) until the desired result is achieved. If the stenosis cannot be crossed with the desired dilatation catheter, use a smaller diameter catheter to predilate the lesion to facilitate passage of a more appropriately sized dilatation catheter.
- Apply negative pressure to fully deflate the balloon. Confirm that the balloon is fully deflated under fluoroscopy.
- Withdraw the balloon catheter until it is clear of the lesion. Maintain the guidewire across the stenosis. Perform angiography to confirm dilatation.
- While maintaining negative pressure, withdraw the deflated dilatation catheter and guidewire from the guiding sheath/catheter through the hemostatic valve. Tighten the hemostatic valve.

Table 1. Typical Coyote ES OTW PTA Balloon Dilatation Catheter Compliance

Balloon Diameter (mm)							
Pressure atm (kPa)	1.5 mm	2.0 mm	2.5 mm	3.0 mm	3.5 mm	4.0 mm	
3.0 (304)	1.39	1.86	2.28	2.76	3.19	3.66	
4.0 (405)	1.45	1.93	2.37	2.85	3.30	3.80	
5.0 (507)	1.48	1.99	2.44	2.93	3.39	3.88	
6.0 (608) Nominal	1.52*	2.03*	2.50*	3.00*	3.46*	3.96*	
7.0 (709)	1.55	2.07	2.55	3.06	3.52	4.04	
8.0 (811)	1.58	2.10	2.59	3.11	3.57	4.09	
9.0 (912)	1.60	2.13	2.62	3.15	3.61	4.14	
10.0 (1013)	1.62	2.15	2.65	3.18	3.64	4.18	
11.0 (1115)	1.64	2.18	2.67	3.21	3.68	4.22	
12.0 (1216)	1.65	2.19	2.69	3.23	3.72**	4.25**	
13.0 (1317)	1.66	2.21	2.72	3.26	---	---	
14.0 (1419)	1.68**	2.23**	2.74**	3.28**	---	---	
*Nominal Pressure = 6 atm							
**Rated Burst Pressure = 14 atm (1.5 mm) DO NOT EXCEED							
**Rated Burst Pressure and Stent Rated Burst Pressure = 14 atm (2.0 mm - 3.0 mm) and 12 atm (3.5 mm - 4.0 mm). DO NOT EXCEED							

Dilatation Catheter Exchange Procedure

To perform a catheter exchange:

- Apply negative pressure to fully deflate the balloon. Confirm that the balloon is fully deflated under fluoroscopy.
- Hold the distal end of the guidewire stationary to maintain position of wire and catheter in the targeted artery, and then begin withdrawing the deflated dilatation catheter out of the guiding sheath/catheter while monitoring the wire position under fluoroscopy.
- When entire device is withdrawn through the hemostatic valve, carefully remove the deflated balloon catheter from the guidewire while maintaining guidewire position across the lesion. Close the knurled knob on the hemostatic valve onto the wire to hold it in place, if applicable.
- Prepare the next balloon catheter to be used as previously described.
- Backload the new catheter onto the guidewire as previously described under the "Use of the Coyote ES OTW PTA Balloon Dilatation Catheter", step 5, and continue the procedure accordingly.

Note: Do not use the catheter if the shaft has been bent or kinked.

Disposal

To minimize the risk of infection or microbial hazards after use, dispose device and packaging as follows:

After use, device and packaging may contain biohazardous substances. Any device and packaging that came into contact with biohazardous substances should be treated and disposed of as biohazardous waste or be treated and disposed of in accordance with any applicable hospital, administrative, and/or local government regulations. Use of a biohazardous container with biological hazard symbol is recommended. Untreated biohazardous waste should not be disposed of in the municipal waste system.

Post-Procedure

Assess patient for hematoma and/or other signs of bleeding at the puncture site.

Any serious incident that occurs in relation to this device should be reported to the manufacturer and relevant local regulatory authority.

INFORMATION TO BRIEF THE PATIENT

The physician should consider the following points while counseling patients on the use of Coyote ES OTW PTA Balloon Dilatation Catheter in association with the interventional procedure:

Discuss the risks and benefits including review of potential adverse events listed in this document, both for the Coyote ES OTW PTA Balloon Dilatation Catheter and for other interventional treatments likely to be employed.

Discuss patient allergies in particular the risk for patients who may be allergic to contrast.

Discuss the post procedure instructions, including any follow-up requirements, lifestyle changes, medications, and home-care or rehabilitation guidelines.

WARRANTY

For device warranty information, visit (www.bostonscientific.com/warranty).

Coyote, Bioslide, and Xtra are trademarks of Boston Scientific Corporation or its affiliates.

All other trademarks are the property of their respective owners.

SYMBOL DEFINITIONS

Commonly used medical device symbols that appear on the labeling are defined at www.bostonscientific.com/SymbolsGlossary.

Additional symbols are defined at the end of this document.

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All rights reserved.
Boston Scientific Limited
Boston Scientific (Australia) Pty Ltd
© 2023 Boston Scientific Corporation or its affiliates.