EmergeTM

MONORAILTM

OVER-THE-WIRE

PTCA Dilatation Catheter

en	Instructions for Use	2
es	Instrucciones de uso	9
fr	Mode d'emploi	16
de	Gebrauchsanweisung	23
it	Istruzioni per l'uso	30
nl	Instructies voor gebruik	37
pt-EU	Instruções de Utilização	44

TABLE OF CONTENTS

REUSE WARNING STATEMENT 3

DEVICE DESCRIPTION 3

 Device Name 3

 Description 3

 Table 1. Emerge Balloon Coatings 3

 Contents 3

 Operating Principle 3

 Materials 4

 Non-Pyrogenic 4

 User Information 4

INTENDED USE / INDICATIONS FOR USE 4

 Clinical Benefit Statement 4

 Summary of Safety and Clinical Performance 4

CONTRAINDICATIONS 4

WARNINGS 4

PRECAUTIONS 4

ADVERSE EVENTS 5

HOW SUPPLIED 5

 Device Details 5

 Handling and Storage 5

OPERATIONAL INSTRUCTIONS 5

 Additional Items for Safe Use 5

 Preparation 5

 Procedure 6

 Disposal 7

 Post-Procedure 7

INFORMATION TO BRIEF THE PATIENT 7

REFERENCES 8

 Table 2. Typical Emerge PTCA Dilatation Catheter Balloon Compliance 8

WARRANTY 8

SYMBOL DEFINITIONS 8

 Table 3. Symbol Definitions 51

Emerge™

MONORAIL™

OVER-THE-WIRE

PTCA Dilatation Catheter

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

REUSE WARNING STATEMENT

For single patient use only. There has been no testing to validate this device beyond use in a single patient. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

DEVICE DESCRIPTION

Device Name

Boston Scientific Emerge Over-The-Wire PTCA Dilatation Catheter and Emerge Monorail PTCA Dilatation Catheter. The generic name of the device is Over-The-Wire Percutaneous Transluminal Coronary Angioplasty Dilatation Catheter/Rapid Exchange Percutaneous Transluminal Coronary Angioplasty Dilatation Catheter.

Description

The Emerge Over-The-Wire (OTW) and Emerge Monorail (MR) Percutaneous Transluminal Coronary Angioplasty (PTCA) Dilatation Catheters, from Boston Scientific, are Over-The-Wire and rapid exchange catheters, respectively, with a semi-compliant balloon near the distal tip. The distal section of both catheters (and the proximal section of the OTW catheter) is dual lumen and coaxial. The outer lumen is used for inflation of the balloon, and the inner lumen permits the use of ≤0.014 in (0.36 mm) guidewires to facilitate advancement of the catheter to and through the stenosis or stent to be dilated. The proximal section of the rapid exchange catheter is a single-lumen, stainless steel hypotube with a single luer port hub for inflation/deflation of the balloon. The OTW catheter has a dual luer port hub: one for inflation/deflation of the balloon, the other for guidewire lumen access. The balloon is designed to provide an inflatable segment of known diameter and length at recommended pressures. A balloon protector is placed over the balloon to maintain a low profile and a mandrel is placed into the inner lumen to protect the patency of the catheter. The catheter's tip is tapered to facilitate advancement of the catheter to and through the stenosis or stent. In addition to the standard design, the 1.20 mm and 1.50 mm catheters are offered in "Push" models that provide different performance characteristics. All shafts have ZGlide (hydrophilic) coating. For MR, the ZGlide is located from the guidewire port to just proximal of the proximal balloon waist. For OTW, the ZGlide is located from distal

of the proximal marks to just proximal of the proximal balloon waist. All balloons have Xtra (hydrophobic) coating and some balloons have ZGlide applied from the distal tip to just proximal of the balloon, per Table 1.

Table 1. Emerge Balloon Coatings

		Balloon Length				
		8 mm	12 mm	15 mm	20 mm	30 mm
Balloon Diameter	1.20 mm					
	1.20 mm Push					
	1.50 mm					
	1.50 mm Push					
	2.00 mm	Xtra Only				
	2.25 mm					
	2.50 mm					
	2.75 mm					
	3.00 mm					
	3.25 mm					
	3.50 mm					
	3.75 mm					
	4.00 mm					
		ZGlide and Xtra				

The effective length of the OTW is 143 cm and the rapid exchange catheter is 144 cm. Marks on the proximal portion of the catheter shaft (at 90cm and 100cm) indicate the exit of the balloon catheter tip out of the guide catheter.

Radiopaque marker bands, in conjunction with fluoroscopy, aid in the placement of the catheter's balloon segment. The 1.20 mm and 1.50 mm models have one radiopaque marker band, while all other models have two radiopaque marker bands. A CLIPIT Hypotube Clip is provided with the Emerge Monorail PTCA Catheter to aid in handling of the catheter.

Contents

Quantity Material

1	Emerge (Monorail or Over-The-Wire) PTCA Dilatation Catheter
1	CLIPIT Hypotube Clip (Monorail Catheter only)

Operating Principle

PTCA is a technique that involves enlargement of a narrowed coronary artery with an inflatable balloon attached to the distal tip of a catheter. During the procedure, while the balloon is deflated, the catheter is positioned into the narrowed coronary artery via a guiding catheter entered through a small opening in the femoral or radial artery. Once positioned inside the narrowed artery, the balloon is inflated, and the inner lumen of the artery is enlarged so that blood can flow more easily. The balloon is then deflated, and the catheter is removed. Similarly, the balloon catheters can also be used for the post-delivery expansion of balloon expandable stents.

Materials

Emerge Monorail PTCA Dilatation Catheter has a bare stainless steel shaft, a polymer distal section, a Pebax balloon, platinum-iridium markers and a hydrophilic (Zglide) and/or hydrophobic (Xtra) coating per Table 1.

Emerge Over-The-Wire PTCA Dilatation Catheter has a polymer proximal shaft, a polymer distal section, a Pebax balloon, platinum-iridium markers and a hydrophilic (Zglide) and/or hydrophobic (Xtra) coating per Table 1.

Contains cobalt: CAS No. 7440-48-4; EC No. 231-158-0. Defined as a 1B carcinogen and reproductive toxicant according to the European Commission in a concentration above 0.1 % weight by weight.

Note: This device is made with a metal alloy which contains cobalt. Current scientific evidence supports that metal alloys containing cobalt used in medical devices, do not cause an increased risk of cancer or adverse reproductive effects.

Non-Pyrogenic

The EmERGE PTCA Dilatation Catheter is non-pyrogenic.

User Information

Physicians experienced in interventional cardiology procedures.

INTENDED USE / INDICATIONS FOR USE

The EmERGE Over-The-Wire and EmERGE Monorail PTCA Dilatation Catheters (balloon models 1.20 mm - 4.00 mm) are intended for the balloon catheter dilation of the stenotic portion of a coronary artery or bypass graft stenosis for the purpose of improving myocardial perfusion.

Emerge Over-The-Wire and EmERGE Monorail PTCA Dilatation Catheters (balloon models 2.00 mm - 4.00 mm) are also intended for the post-delivery expansion of balloon expandable stents.

Note: Bench testing was conducted with EmERGE Over-The-Wire and EmERGE Monorail PTCA Dilatation Catheters and marketed Boston Scientific balloon expandable stents. Consideration should be taken when this device is used with other manufacturers' stents due to differences in stent design.

Clinical Benefit Statement

The EmERGE Over-The-Wire and EmERGE Monorail PTCA Dilatation Catheter is a minimally invasive medical device used for the purpose of dilating coronary artery stenoses to improve myocardial perfusion in patients with stenotic coronary lesions.

Summary of Safety and Clinical Performance

For customers in the European Union, use the device name found in the labeling to search for the device's Summary of Safety and Clinical Performance, which is available on the European database on medical devices (EUDAMED) website: (<https://ec.europa.eu/tools/eudamed>).

CONTRAINDICATIONS

- Unprotected left main coronary artery.
- Coronary artery spasm in the absence of a significant stenosis.

WARNINGS

- To reduce the potential for vessel damage, the inflated diameter of the balloon should approximate, or be less than, the diameter of the vessel just proximal and distal to the stenosis.
- PTCA in a patient who is not an acceptable candidate for coronary artery bypass graft surgery requires careful consideration, including possible hemodynamic support during PTCA.
- Use extreme caution and careful judgment in patients for whom anticoagulation is not indicated.
- Use extreme caution and careful judgment in patients who have severe reaction to contrast agents that cannot be adequately pre-medicated.
- When the catheter is exposed to the vascular system, it should only be manipulated while under high-quality fluoroscopic observation. Do not advance or retract the catheter unless the balloon is fully deflated and under vacuum. If resistance is met during manipulation, determine the cause of the resistance before proceeding.
- Do not exceed the rated balloon burst pressure. The rated burst pressure is based on the results of in vitro testing. At least 99.9 percent of the balloons (with 95 percent confidence) will not burst at or below their rated burst pressure. A pressure monitoring device is recommended during use to prevent over-pressurization.
- If difficulty is experienced during balloon inflation, do not continue; remove the catheter.
- PTCA should only be performed at hospitals where emergency coronary artery bypass graft surgery can be quickly performed in the event of a potentially injurious or life-threatening complication.
- Use only the recommended balloon inflation medium. Never use air or any gaseous medium to inflate the balloon.
- Before withdrawing the balloon catheter, visually confirm complete balloon deflation by fluoroscopy.
- Balloon catheter retrieval methods (use of additional wires, snares and/or forceps) may result in additional trauma to the vascular access site. Complications can include but are not limited to bleeding, hematoma, or pseudoaneurysm.
- Use the catheter prior to the "Use By" date specified on the package.
- Used devices may pose a biohazard risk and must be handled and disposed of properly.

Note: If balloon burst occurs, carefully remove the catheter and assess the lesion for potential vessel damage or embolization. If the device is kinked or damaged, remove it and replace with another device.

PRECAUTIONS

- The compatibility of the device has not been evaluated for the delivery of materials (e.g., drugs, alcohol, or stem cells) through the guidewire lumen, other than those required for normal use.
- Prior to angioplasty, the balloon catheter should be examined to verify functionality and to ensure that its size and shape are suitable for the specific procedure for which it is to be used.
- Before insertion of the balloon catheter, administer appropriate anticoagulant and coronary vasodilator therapy.
- Care should be taken to control the position of the guide catheter tip during manipulation of the balloon catheter.

- Caution should be taken not to overtighten a hemostatic adapter around the catheter shaft as lumen constriction may occur, affecting inflation/deflation of the balloon.
- When loading or exchanging the balloon catheter, it is recommended to thoroughly wipe the guidewire clean for better catheter movement on the guidewire.
- Do not expand the balloon if it is not properly positioned in the vessel.

ADVERSE EVENTS

Potential adverse events (in alphabetical order) that may be associated with the use of a PTCA Dilatation Catheter include, but are not limited to, the following:

- *abrupt closure, slow flow/no reflow*
- *arrhythmia, including ventricular fibrillation or heart block*
- *bleeding, hemorrhage or hematoma*
- *cardiac arrest/cardiogenic shock*
- *cardiac tamponade/pericardial effusion*
- *cerebrovascular accident (stroke or transient ischemic attack)*
- *death*
- *drug reactions including allergy*
- *embolism (air, tissue, device fragments, plaque)*
- *hemodynamic compromise including vasovagal reaction*
- *infection*
- *myocardial ischemia or infarction*
- *organ insufficiency/failure (heart, lung, kidney)*
- *pain (angina)*
- *prolonged procedure time including additional intervention or surgery*
- *radiation injury*
- *restenosis of the dilated vessel*
- *thrombosis*
- *vessel injury (spasm, dissection, perforation, rupture, arteriovenous fistula, pseudoaneurysm)*

HOW SUPPLIED

Device Details

Contents supplied sterile using an ethylene oxide (EO) process. Do not use if package is damaged or unintentionally opened before use. If damage is found, call your Boston Scientific representative. Do not use if labeling is incomplete or illegible.

Handling and Storage

This product has no special handling or storage requirements.

OPERATIONAL INSTRUCTIONS

Additional Items for Safe Use

One or more of each of the following materials are required for PTCA but not supplied with the Emerge Over-The-Wire (OTW) or Emerge Monorail (MR) PTCA Dilatation Catheter:

- Femoral or radial guide catheter access to the target vessel;
 - minimum I.D. of guide catheter = 0.066 in (1.68 mm) (Emerge OTW Catheter)
 - minimum I.D. of guide catheter = 0.056 in (1.42 mm) (Emerge MR Catheter)
- Contrast media
- Sterile saline or heparinized sterile saline

- Inflation device
- ≤0.014 in (0.36 mm) x 300 cm guidewire(s) (Emerge OTW Catheter)
- ≤0.014 in (0.36 mm) x 185 cm guidewire(s) (Emerge MR Catheter)
- 10 ml, 12 ml, or 20 ml (10 cm³ (cc), 12 cm³ (cc), or 20 cm³ (cc)) luer-lock syringe
- Hemostatic adapter
- Three-way stopcock

Prior to angioplasty, carefully examine all equipment to be used during the procedure, including the dilation catheter, to verify proper function. Verify that the catheter and sterile packaging have not been damaged. Verify that the catheter size is suitable for the specific procedure for which it is intended.

Note: Do not continue to use the catheter if damage occurs or sterility is compromised during use.

Preparation

1. Inflation Device Preparation

- Prepare the inflation device according to the manufacturer's instructions.
- Purge the system of air.

2. Dilation Catheter Selection

- The inflation diameter of the balloon must not exceed the diameter of the coronary artery proximal and distal to the stenosis. If the stenosis cannot be crossed with the desired dilation catheter, use a smaller diameter catheter to predilate the lesion to facilitate passage of a more appropriately sized dilation catheter.

3. Dilation Catheter Preparation

- Remove the catheter from the protective hoop. Use care when removing the catheter to avoid damage (e.g., shaft kink).
- Remove the balloon protector and mandrel by grasping the balloon dilation catheter just proximal to the balloon, and with the other hand, gently grasp the proximal section of the balloon protector and slide distally. For Emerge OTW Catheters, the mandrel will slide off with the balloon catheter protector. For Emerge MR Catheters, remove the mandrel distally after removing the balloon protector. If unusual resistance is felt during the product mandrel and/or balloon protector removal, do not use this product and replace with another.
- The Emerge MR Catheter may be coiled once and secured using the CLIPIT Clip provided in the catheter package. Only the proximal shaft should be inserted into the CLIPIT clip; the clip is not intended for the distal end of the catheter. Remove the CLIPIT Clip prior to inserting the catheter into the patient's body.

Note: Care should be taken not to kink the shaft of the catheter upon application or removal of the CLIPIT Clip.

- Prepare the balloon dilation catheter for purging. Fill a 10 ml, 12 ml, or 20 ml (10 cm³ (cc), 12 cm³ (cc), or 20 cm³ (cc)) luer-lock syringe with 3 ml (3 cm³ (cc)) of contrast medium. Use only appropriate balloon inflation medium (e.g., the equivalent of a 50:50 mixture of contrast medium and sterile normal saline). Never use air or any gaseous medium to inflate the balloon.
- Connect a three-way stopcock to the port fitting on the dilation catheter. Flush through the stopcock taking care to ensure that the balloon cannot be inflated. Connect syringe or inflation device to stopcock.

- f. Hold the syringe or inflation device with the nozzle pointing downward and aspirate for 5 seconds. Release the plunger or open stopcock to air.
- g. Remove the syringe or inflation device and evacuate all air from the barrel.
- h. Reconnect the syringe and aspirate until bubbles no longer appear during aspiration. If bubbles persist despite verification of a luer-lock seal, do not use the device.

4. Inflation Device Connection to Catheter

- a. To remove any air lodged in the distal luer fitting of the inflation device, purge approximately 1 ml (cm³ (cc)) of contrast medium while holding the inflation device pointing upwards.
- b. Disconnect the syringe used in preparation, taking care to completely fill the hub with fluid. A meniscus should appear in the balloon port when the syringe is removed. Verify that a meniscus of inflation medium is evident in both the dilation catheter balloon port and the inflation device connection. Securely couple the inflation device to the balloon port of the balloon dilation catheter.

Procedure

5. Guidewire Lumen Flush

- a. For Emerge MR Catheters, flush the guidewire lumen of the balloon catheter with sterile saline through the distal tip of the catheter.
- b. For Emerge OTW Catheters, flush the guidewire lumen of the catheter with sterile saline through the guidewire port of the catheter hub.
- c. Check for bends, kinks and other damage. Do not use if any defects are noted.

6. Catheter Advancement

- a. Prepare the vascular access site according to standard practice.
- b. Maintain neutral pressure on the inflation device attached to the catheter.
- c. Insert a guidewire through the hemostatic adapter following the manufacturer's instructions or standard practice. Advance the guidewire carefully into the guide catheter.
- d. Attach a torque device to the guidewire, if desired. Under fluoroscopic guidance, advance the guidewire to the desired vessel, then across the stenosis or stent.
- e. Backload the distal tip of the dilation catheter onto the guidewire ensuring that the guidewire exits the midsection opening in the Emerge MR Catheter or the wire port of the Emerge OTW Catheter manifold. When loading or exchanging the catheter, it is recommended to thoroughly wipe the guidewire clean for better catheter movement on the guidewire.

Note: To avoid kinking, advance the catheter slowly, in small increments, until the proximal end of the guidewire emerges from the catheter.

- f. Thoroughly aspirate and flush the guide catheter in preparation for introduction of the dilation catheter.
- g. Carefully advance the catheter through the hemostatic adapter while the balloon is fully deflated. If unusual resistance is felt, do not advance the catheter through the adapter. Caution should be taken not to over tighten the hemostatic adapter around the catheter shaft as lumen constriction may occur, affecting inflation/deflation of the balloon.
- h. Connect the side port of the guide catheter hemostatic adapter to the proximal pressure recording/infusion line or manifold assembly, which permits proximal pressure recording or infusion through the guide catheter.

- i. Advance the catheter over the guidewire under direct fluoroscopic visualization and position the balloon relative to the stenosis or stent to be dilated. Use the radiopaque marker bands as a reference point. The outside edges of the marker bands indicate the balloon shoulders on the 2.00 mm - 4.00 mm models; the 1.20 mm - 1.50 mm models have a single central marker band. Balloon inflation should not be undertaken if the balloon is not properly positioned within the stenosis or stent.

- j. Simultaneous Use of Two Balloon Catheters in a Guide Catheter (Kissing Balloon): Bench and preclinical testing has shown that one 4.00 mm x 30 mm (or smaller) Monorail balloon catheter and one 3.25 mm x 20 mm (or smaller) Monorail balloon catheter can be inserted simultaneously into a 6F (minimum 0.070 in ID) guide catheter and two 4.00 mm x 30 mm (or smaller) Over-the-Wire balloon catheters can be inserted into an 8F (minimum 0.088 in ID) guide catheter with acceptable compatibility. These tests did not account for all clinical situations and differing anatomy. Care should be used when attempting to use two balloon catheters simultaneously in a guide catheter. Balloon catheters with a diameter greater than those mentioned have not been tested for simultaneous use in a single guide catheter.

7. Catheter Inflation

- a. Inflate the balloon slowly to the appropriate pressure to perform dilation. Maintain negative pressure on the balloon between inflations. Do not exceed the rated balloon burst pressure. Refer to Table 2 or to the balloon compliance chart. If difficulty is experienced during balloon inflation, do not continue inflation; deflate and remove the catheter.
- b. After completion of PTCA or post-dilatation of a stent, deflate the balloon by pulling negative pressure on the inflation device until the balloon is fully deflated.
- c. Confirm angiographic results using standard angiographic techniques. Fluoroscopic visualization during balloon expansion should be used to properly judge the optimum expanded balloon diameter as compared to the proximal and distal coronary artery diameter(s). Repeat inflation of balloon until the desired result is achieved.
- d. If catheter exchange is necessary, proceed to step 9 – Catheter Exchange Procedure (Emerge Over-The-Wire Catheter) or step 10 – Catheter Exchange Procedure (Emerge Monorail PTCA Catheter). Otherwise, proceed to step 8 – Catheter Removal.

8. Catheter Removal

- a. Confirm with angiography that the lumen of the dilated artery has not abruptly occluded. Ensure balloon is fully deflated.
- b. While withdrawing the deflated catheter and guidewire from the guide catheter through the hemostatic adapter, tighten the hemostatic adapter.
- c. The Emerge MR Catheter may be coiled once and secured using the CLIPIT Clip provided in the catheter package. Only the hypotube should be inserted into the CLIPIT clip; the clip is not intended for the distal end of the catheter. Remove the CLIPIT Clip prior to the catheter being inserted in the patient's body.

Note: Care should be taken not to kink or bend the shaft upon application or removal of the coil clip.

9. Catheter Exchange Procedure (Emerge Over-The-Wire Catheter)

The Emerge OTW Catheters typically require two operators to exchange. To perform a catheter exchange, execute the following steps:

- a. Open the hemostatic adapter.
- b. The primary operator holds the hemostatic adapter in one hand, while grasping the catheter shaft in the opposite hand.
- c. The secondary operator is positioned near the foot of the patient and should maintain the guidewire position in the coronary artery by holding the guidewire stationary and confirming guidewire position at all times under fluoroscopy while the primary operator begins pulling the catheter out of the guide catheter.
- d. Withdraw the deflated catheter until the catheter tip exits the hemostatic adapter.
- e. Close the hemostatic adapter and remove the catheter from the guidewire while maintaining guidewire position across the stenosis or stent.
- f. Prepare the next therapeutic device to be used.

10. Catheter Exchange Procedure (Emerge Monorail Catheter)

The Emerge MR Catheters have been specifically designed for rapid, single operator catheter exchanges. To perform a catheter exchange, execute the following steps:

- a. Open the hemostatic adapter.
- b. Hold the guidewire and hemostatic adapter in one hand, while grasping the catheter shaft in the opposite hand.
- c. Maintain the guidewire position in the coronary artery by holding the guidewire stationary. Begin pulling the catheter out of the guide catheter while monitoring the guidewire position under fluoroscopy.
- d. Withdraw the deflated catheter until the opening in the guidewire lumen is reached (approximately 25 cm proximal to the balloon catheter tip).
- e. Carefully slide the flexible, distal portion of the catheter out of the hemostatic adapter, and close it onto the guidewire to hold it securely in place. Completely remove the catheter from the guidewire while maintaining guidewire position across the stenosis.
- f. Prepare the next therapeutic device to be used.

Disposal

To minimize the risk of infection or microbial hazards after use, dispose device and packaging as follows: After use, device and packaging may contain biohazardous substances. Any device and packaging that came into contact with biohazardous substances should be treated and disposed of as biohazardous waste or be treated and disposed of in accordance with any applicable hospital, administrative, and/or local government regulations. Use of a biohazardous container with biological hazard symbol is recommended. Untreated biohazardous waste should not be disposed of in the municipal waste system.

Post-Procedure

Any serious incident that occurs in relation to this device should be reported to the manufacturer and relevant local regulatory authority.

INFORMATION TO BRIEF THE PATIENT

The following sections/clauses in the IFU contain information that may be useful for briefing the patient prior to use of this device:

- Clinical Benefit Statement
- Adverse Events statements with potential relevance to the patient are in *italics*

Patients should be informed on post-procedure medications.

REFERENCES

Table 2. Typical Emerge PTCA Dilatation Catheter Balloon Compliance

Pressure atm (kPa)	Balloon Size										
	1.20 mm	1.50 mm	2.00 mm	2.25 mm	2.50 mm	2.75 mm	3.00 mm	3.25 mm	3.50 mm	3.75 mm	4.00 mm
3.0 (304)	1.09	1.39	1.86	2.06	2.28	2.53	2.76	3.01	3.19	3.45	3.66
4.0 (405)	1.13	1.45	1.93	2.14	2.37	2.61	2.85	3.10	3.30	3.58	3.80
5.0 (507)	1.16	1.48	1.99	2.20	2.44	2.68	2.93	3.18	3.39	3.67	3.88
6.0 (608) Nominal	1.20	1.52	2.03	2.26	2.50	2.75	3.00	3.26	3.46	3.75	3.96
7.0 (709)	1.23	1.55	2.07	2.31	2.55	2.81	3.06	3.33	3.52	3.81	4.04
8.0 (811)	1.26	1.58	2.10	2.34	2.59	2.85	3.11	3.37	3.57	3.86	4.09
9.0 (912)	1.28	1.60	2.13	2.38	2.62	2.88	3.15	3.41	3.61	3.90	4.14
10.0 (1013)	1.31	1.62	2.15	2.40	2.65	2.91	3.18	3.45	3.64	3.95	4.18
11.0 (1115)	1.33	1.64	2.18	2.42	2.67	2.94	3.21	3.48	3.68	3.98	4.22
12.0 (1216)	1.35	1.65	2.19	2.44	2.69	2.96	3.23	3.51	3.72**	4.02**	4.25**
13.0 (1317)	1.37	1.66	2.21	2.46	2.72	2.99	3.26	3.54	-	-	-
14.0 (1419)	1.39	1.68*	2.23**	2.48**	2.74**	3.02**	3.28**	3.57**	-	-	-
15.0 (1520)	1.41	-	-	-	-	-	-	-	-	-	-
16.0 (1621)	1.42	-	-	-	-	-	-	-	-	-	-
17 (1723)	1.43	-	-	-	-	-	-	-	-	-	-
18.0 (1824)	1.45*	-	-	-	-	-	-	-	-	-	-
19.0 (1925)	1.46	-	-	-	-	-	-	-	-	-	-
20.0 (2027)	1.48	-	-	-	-	-	-	-	-	-	-
* Rated Burst Pressure. Do not exceed.											
** Rated Burst Pressure and Stent Rated Burst Pressure. Do not exceed.											

The physician should consult current practice guidelines and medical literature on percutaneous coronary intervention, such as those published by ACC, AHA, SCAI, or ESC.

WARRANTY

For device warranty information, visit (www.bostonscientific.com/warranty).

EMERGE, Monorail and CLIPIT Clip are registered trademarks of Boston Scientific Corporation or its affiliates.

All other trademarks are the property of their respective owners.

EU Importer: Boston Scientific International B.V., Vestastraat 6, 6468 EX Kerkrade, The Netherlands.

SYMBOL DEFINITIONS

Commonly used medical device symbols that appear on the labeling are defined at www.bostonscientific.com/SymbolsGlossary.

Additional symbols are defined at the end of this document.

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All rights reserved.
Importador da UE: Boston Scientific International B.V., Vestastraat 6, 6468 EX Kerkrade, Países Baixos.
Boston Scientific Limited
Boston Scientific (Australia) Pty Ltd
© 2024 Boston Scientific Corporation or its affiliates.