RIGIFLEX™ II

Single Use Achalasia Balloon Dilator

REF M00554500, M00554510, M00554520

Instructions for Use	2
Instrucciones de uso	7
Mode d'emploi	12
Gebrauchsanweisung	17
Istruzioni per l'uso	22
Instructies voor gebruik	27
Instruções de Utilização	32

51376251-01 Rev. B

2022-04

TABLE OF CONTENTS

REUSE WARNING3

DEVICE DESCRIPTION3

 Contents.....3

 Operating Principle3

 User Information.....3

 Materials3

INTENDED USE / INDICATIONS FOR USE.....3

 Clinical Benefit Statement.....3

CONTRAINDICATIONS.....3

PRECAUTIONS.....4

ADVERSE EVENTS4

HOW SUPPLIED.....4

 Device Details.....4

 Handling and Storage4

OPERATIONAL INSTRUCTIONS.....4

 Additional Items For Safe Use.....4

 Inspection4

 Preparation and Procedure4

 Balloon Inflation5

 Catheter Withdrawal5

 Disposal5

 Post-Procedure5

PATIENT COUNSELING INFORMATION5

WARRANTY6

SYMBOL DEFINITIONS.....6

RIGIFLEX™ II

Single Use Achalasia Balloon Dilator

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

REUSE WARNING

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

DEVICE DESCRIPTION

Contents

One (1) Rigiflex II Achalasia Balloon Dilator.

One (1) 0.038 in (0.97 mm) floppy tip guidewire.

One (1) guidewire introducer.

Operating Principle

The Rigiflex II Achalasia Balloon Dilator is used for dilatation of the cardia Lower Esophageal Sphincter (LES) in patients with Achalasia. It is comprised of a catheter and a balloon. The balloon is composed of an inelastic material and will be filled with air to exert radial force on a narrowed segment in the alimentary tract resulting from achalasia. This balloon can be inflated with the Boston Scientific Achalasia Pneumatic Inflator.

User Information

Read the entire Instructions for Use before using the Rigiflex II Single Use Achalasia Balloon Dilator. The Rigiflex II Single Use Achalasia Balloon Dilator should only be used by or under the supervision of physicians thoroughly trained in endoscopic and/or fluoroscopic procedures. A thorough understanding of the technical principles, clinical applications, and risk associated with this procedure are necessary before using this device.

Materials

Contains cobalt: CAS No. 7440-48-4; EN No. 231-158-0. Defined as a 1B carcinogen and reproductive toxicant according to the European Commission in a concentration above 0.1% weight by weight.

Note: This device is made with stainless steel which contains cobalt. Current scientific evidence supports that stainless steels used in medical devices, do not cause an increased risk of cancer or adverse reproductive effects.

INTENDED USE / INDICATIONS FOR USE

The Rigiflex II Achalasia Balloon Dilator is indicated for dilatation of the cardia (LES) in patients with achalasia.

Clinical Benefit Statement

The clinical benefit of The Rigiflex II Achalasia Balloon Dilator is dilation of the lower esophageal sphincter (LES) to allow for food and liquid to pass through the esophagus to the stomach in patients with achalasia.

CONTRAINDICATIONS

- Acute or incompletely healed esophageal perforation.
- Bleeding disorders.

- Large thoracic aortic aneurysm.
- Pharyngeal or cervical deformity.
- Recent esophageal perforation or surgery.
- Recent laparotomy.
- Recent myocardial infarction.
- Severe pulmonary disease.

PRECAUTIONS

- Exercise care when handling the floppy tip of the guidewire to avoid damage.
- Insert the floppy tip of the guidewire first into the working channel of the endoscope.
- Do not exceed the maximum inflation pressure of the balloon.
- If the balloon is over inflated beyond the recommended pressure, balloon damage can occur.

ADVERSE EVENTS

Complications may include:

- Aspiration
- Bleeding
- Esophageal Intramural Hematoma
- Esophageal Tear
- Infection
- Inflammation
- Pain (e.g., throat, stomach, chest, etc)
- Perforation
- Pneumonia
- Reflux
- Respiratory Distress

HOW SUPPLIED

Device Details

The Rigiflex II Achalasia Balloon Dilator is supplied in a Tyvek pouch and is supplied sterile using an ethylene oxide (EO) process and is intended for single use only.

Do not use if package is damaged or unintentionally opened before use.

Do not use if labeling is incomplete or illegible.

Handling and Storage

This product has no special handling or storage requirements.

OPERATIONAL INSTRUCTIONS

Additional Items For Safe Use

Achalasia Pneumatic Inflator.

Inspection

1. Prior to opening package, inspect seal integrity and package surface integrity to ensure sterile barrier has not been compromised.
2. Open the package and remove the Rigiflex II Achalasia Balloon Dilator and guidewire hoop from the pouch.
3. Prior to use, carefully examine the catheter, balloon, and guidewire to verify that they were not damaged in shipment. Verify that the catheter tip is not frayed and that the catheter and guidewire are free of sharp ends and kinks.

Note: The guidewire is packaged with the floppy tip of the guidewire at the inner diameter of the hoop.

Preparation and Procedure

1. Attach a hospital supplied 60 ml syringe to the catheter.
2. Draw back on syringe and create a negative vacuum, drawing the balloon to its lowest profile.
3. Remove the guidewire introducer from the guidewire hoop.
4. Insert the guidewire introducer into the biopsy cap of the endoscope.

5. After full endoscopic exam, insert the floppy tip of the provided guidewire, indicated by the flag on the guidewire hoop, through the working channel of the scope and cross the LES.
6. Once the guidewire is in place, remove the endoscope, leaving the guidewire in place.
7. The balloon catheter is then backloaded over the guidewire and the balloon is centered in the LES.
8. The scope may be reinserted for endoscopic visualization.
9. The Rigiflex II Achalasia Balloon Dilator may also be used under fluoroscopic guidance. The catheter has 5 radiopaque markers: one at the distal end of the balloon, one at the proximal end of the balloon, two in the center of the balloon and one located under the proximal balloon sleeve. These positioning markers are easily seen under fluoroscopy as is the radiopaque tip. The soft atraumatic tip at the end of the catheter acts as a guide for the balloon catheter and gently dilates the esophagus and sphincter as the catheter is advanced.

Balloon Inflation

When the balloon is positioned in the esophagogastric junction, inflate, using the Achalasia Pneumatic Inflator. The inflator is attached to the balloon lumen on the catheter and is operated like the pump on a blood pressure cuff. Inflate the balloon per the recommended pressure until no balloon waisting is observed. Do not exceed the maximum inflation pressure of 20 psi (1.4 ATM, 138 kPa). Depending on the patient, you may wish to deflate and reinflate the balloon several times until the criteria for sufficient dilatation are met.

Warning: Check for proper position of the balloon catheter using endoscopic or fluoroscopic visualization. Balloon inflation in an improper location may lead to patient injury.

Warning: If the balloon is over inflated beyond the recommended pressure, patient injury such as aspiration, tissue damage, bleeding and/or perforation may occur.

Catheter Withdrawal

1. When dilation is complete, deflate the balloon with a 60 ml syringe to its lowest profile by applying suction to the balloon lumen.

Note: The greater the suction applied and held during withdrawal the smaller the deflated profile of the balloon. For maximum suction a 60 ml (60 cc) syringe is recommended.

2. Once the balloon is at its lowest profile, slowly remove the balloon and guidewire from the patient.

Disposal

To minimize the risk of infection or microbial hazards after use, dispose device and packaging as follows:

After use, device may contain biohazardous substances. The device and packaging should be treated and disposed of as biohazardous waste or have them treated and disposed of in accordance with any applicable hospital, administrative, and/or local government regulations. Use of a biohazardous container with biological hazard symbol is recommended. Untreated biohazardous waste should not be disposed of in the municipal waste system.

Post-Procedure

Any serious incident that occurs in relation to this device should be reported to the manufacturer and relevant local regulatory authority.

PATIENT COUNSELING INFORMATION

Brief the patient on any relevant post-procedure instructions, contraindications, warnings, precautions and/or adverse events found within this Instructions for Use (IFU), pertaining to the patient.

WARRANTY

For device warranty information, visit (www.bostonscientific.com/warranty).

Rigiflex is a trademark of Boston Scientific Corporation or its affiliates.

All other trademarks are the property of their respective owners.

SYMBOL DEFINITIONS

Commonly used medical device symbols that appear on the labeling are defined at www.bostonscientific.com/SymbolsGlossary.

Additional symbols are defined at the end of this document.

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or its affiliates. All rights reserved.
Boston Scientific (Australia) Pty Ltd
©2022 Boston Scientific Corporation