BLAZER™ Dx-20

Bidirectional Steerable Diagnostic Catheter

Instructions for Use	2
Instrucciones de uso	7
Mode d'emploi	12
Gebrauchsanweisung	17
Istruzioni per l'uso	22
Instructies voor gebruik	27
Instruções de Utilização	32

TABLE OF CONTENTS

REUSE WARNING 3

DEVICE DESCRIPTION 3

 Contents 3

 Operating Principle..... 3

 Non-pyrogenic..... 3

 User Information..... 3

 Figure 1. BLAZER Dx-20 Bidirectional Steerable Diagnostic Catheter 3

INTENDED USE 3

INDICATIONS FOR USE..... 3

 Clinical Benefit Statement 3

 Summary of Safety and Clinical Performance..... 3

CONTRAINDICATIONS..... 3

WARNINGS..... 3

PRECAUTIONS..... 4

ADVERSE EVENTS 4

HOW SUPPLIED..... 5

 Device Details 5

 Handling and Storage..... 5

OPERATIONAL INSTRUCTIONS..... 5

 Additional Required Items..... 5

 Preparation..... 5

 Procedure 5

 Disposal 5

 Post-Procedure 5

INFORMATION TO BRIEF THE PATIENT..... 5

 Figure 2. Diagnostic Catheter Connectivity Diagram A 5

 Figure 3. Diagnostic Catheter Connectivity Diagram B..... 6

WARRANTY..... 6

SYMBOL DEFINITIONS..... 6

BLAZER™ Dx-20

Bidirectional Steerable Diagnostic Catheter

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

REUSE WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

Carefully read all ancillary device instructions prior to use. Observe all contraindications, warnings and precautions noted in these instructions. Failure to do so may result in patient complications, failure of the device to perform as intended, or both.

DEVICE DESCRIPTION

The BLAZER Dx-20 Bidirectional Steerable Diagnostic Catheter (henceforth, referred to as the BLAZER Dx-20 Catheter) is a duodecapolar, bidirectional, steerable diagnostic catheter designed for use in intracardiac pacing and recording only. It has a distal electrode segment connected to an ergonomically designed handle via a 7F torquable shaft. The 2 mm tip electrode and nineteen ring electrodes are designed to carry electrical signals for the purpose of endocardial stimulation (pacing) and/or recording. The BLAZER Dx-20 Catheter is offered in multiple electrode pair configurations and different curve configurations to assist the physician in reaching the desired anatomical location. The BLAZER Dx-20 Catheter is compatible with most commercially available recording and mapping systems and connects to these systems via the BLAZER Dx Cable.

Contents

One (1) Diagnostic Catheter

Operating Principle

The electrode segment of the catheter is comprised of a distal tip electrode and nineteen ring electrodes. The tip and ring electrodes are designed to carry electrical signals for endocardial stimulation (pacing) and/or recording.

The degree of tip deflection of the BLAZER Dx-20 Catheter is controlled by the steering knob on the catheter handle (see Figure 1). If the steering knob is turned in a clockwise direction from its neutral position, the tip will curve proportionately in the corresponding direction. Turning the steering knob in the counterclockwise direction will cause the tip to deflect in the opposite direction. To hold the tip in the desired position, the catheter has an adjustable tension control knob. In the neutral '(-)' position, which is the minimum available tension, the catheter steers freely and will not hold a preset curve. In the '(+)' position, maximum tension is achieved. To increase steering resistance and/or hold a preset curve, the tension control knob can be rotated as necessary towards the '(+)' position.

The catheter is placed into point of interest positions in the heart and is guided to location by steering the distal tip area of the catheter. Tip and ring electrodes come into contact with the endocardium where electrical contact is made and pacing and recording of electrograms becomes possible. The BLAZER Dx-20 Catheter is designed to be used with ancillary devices such as: sheath, ECG Recording System, Pacing Stimulator and with the BLAZER Dx-20 cable accessory. The handle includes the electrical connector for the cable connection to the standard recording equipment.

Non-pyrogenic

This device meets pyrogen limit specifications for all patient-contacting parts.

User Information

The BLAZER Dx-20 Catheter is to be used only by physicians thoroughly trained and authorized to use intravascular catheters and in the specific approach to be used, and is to be used in a fully-equipped electrophysiology lab.


Figure 1. BLAZER Dx-20 Bidirectional Steerable Diagnostic Catheter

INTENDED USE

The BLAZER Dx-20 Catheter is intended for temporary use in electrophysiology studies for intracardiac stimulation (pacing) and/or recording of electrical potentials. The device is intended for use in adult (not pediatric) patients, with the exclusion of pregnant and/or nursing patients.

INDICATIONS FOR USE

The BLAZER Dx-20 Catheter is indicated for use to diagnose cardiac arrhythmia.

Clinical Benefit Statement

BLAZER Dx-20 Bidirectional Duodecapolar Steerable Diagnostic Catheter is intended for use in electrophysiology studies for intracardiac stimulation (pacing) and/or recording of electrical potentials. When operated by the intended user with commercially available recording and mapping systems, the Clinical Benefit of the BLAZER Dx-20 catheter is to aid in diagnosis of cardiac arrhythmias. Specifically, the BLAZER Dx-20 Catheter provides:

- Back Steering to facilitate positioning in the coronary sinus
- Torque response to aid in catheter handling and placement

Conversely, failure to accurately diagnose and treat a cardiac arrhythmia could result in symptoms such as shortness of breath, palpitations, dizziness, syncope, chest pain or sudden cardiac death.

Summary of Safety and Clinical Performance

For customers in the European Union, use the device name found in the labeling to search for the device's Summary of Safety and Clinical Performance, which is available on the European database on medical devices (EUDAMED) website: (<https://ec.europa.eu/tools/eudamed>).

CONTRAINDICATIONS

The BLAZER Dx-20 Catheter is contraindicated for use:

- in patients with a mechanical prosthetic heart valve through which the catheter must pass;
- in patients with active systemic infection;
- in patients with conditions where insertion into or manipulation in the cardiac chambers is unsafe as this may increase the risk of perforation or systemic embolic event, such as but not limited to, evidence of myxoma, Left Atrial (LA) thrombus, or intracardiac mural thrombus, or interatrial baffle or patch;
- from the femoral approach in patients who have vena cava embolic protection filter devices or known femoral thrombus;
- in patients who are hemodynamically unstable;
- in patients with a contraindication to an invasive electrophysiology procedure where insertion or manipulation of a catheter in the cardiac chambers is deemed unsafe, such as but not limited to, a recent previous cardiac surgery (e.g., ventriculotomy or atriotomy, Coronary Artery Bypass Graft [CABG], PTCA/PCI/coronary stent procedure/unstable angina) and/or in patients with congenital heart disease where the underlying abnormality increases the risk of the ablation (e.g. severe rotational anomalies of the heart or great vessels).

WARNINGS

- Cardiac catheterization procedures should be performed only by physicians thoroughly trained in invasive cardiology and in the specific approach to be used, in a fully-equipped electrophysiology lab.

- All equipment used in connection with the BSC Catheter must meet IEC 60601-1 requirements, including those for use with Type CF Defibrillation-Proof applied parts, as well as any other relevant regulatory and safety standards. The use of catheter with devices that do not comply with relevant standards may cause equipment damage or system malfunction, and/or harm to patient and/or user.
- Before using, inspect the BLAZER Dx-20 Catheter for any defects or physical damage, including electrical insulation on the cables and the catheter shaft that, if used, may cause patient and/or user injury. Do not use defective or damaged devices. Replace damaged equipment if necessary. No modification of this equipment is allowed.
- The use of catheters or cables with unprotected male pin connectors presents a risk of electrical hazard. Inadvertent attachment of pin connectors to power supply sockets of connectors could result in electrocution of the patient or operator.
- Do not use BLAZER Dx-20 Catheter as an internal defibrillation catheter. Doing so may result in perforation, arrhythmias, embolism, thrombus, and/or patient death.
- Cardiac catheterization procedures present the potential for significant radiation exposure, which can result in acute radiation injury as well as an increased risk for somatic and genetic effects, to both patients and laboratory staff due to the radiation beam intensity and duration of the fluoroscopic imaging. Cardiac catheterization should only be performed after adequate attention has been given to the potential radiation exposure associated with the procedure, and steps have been taken to minimize this exposure. Due to radiation exposure during cardiac catheterization, the safety and effectiveness of this device has not yet been established in pregnant and/or nursing women and pediatric patients.
- Careful catheter manipulation must be performed to avoid cardiac damage, perforation, or tamponade. Catheter advancement should be performed under fluoroscopic guidance. Do not use excessive force to advance or withdraw the catheter when resistance is encountered.
- Administer appropriate levels of peri-procedural anticoagulation therapy for selected patients undergoing right-sided procedures. Administer anticoagulation therapy during and post-procedure according to institution standards to minimize bleeding and thrombotic complications.
- Stimulation of cardiac tissues caused by pacing stimuli may lead to inadvertent induction of arrhythmias. These arrhythmias may require defibrillation that could also result in skin burns.
- Do not use if package is opened or damaged.
- Carefully read all equipment and ancillary device instructions required for the procedure prior to use. Observe all contraindications, warnings, and precautions noted in these instructions. Failure to do so may result in patient complications.
- This catheter is not indicated for use in Cardiac Ablation or Coronary Artery Mapping.
- Ensure that the cable/catheter connection remains dry throughout the procedure in order to prevent electric shock or other patient injuries as well as to prevent loss of device function.
- Fibrin may accumulate in or on the sheath/catheter assembly during the procedure. Aspirate when removing the dilator or catheter.
- In the presence of anticoagulation, there may be an increased risk of bleeding from all causes.
- Do not insert or withdraw the BLAZER Dx-20 Catheter without straightening the catheter tip (returning the steering lever to neutral position) in order to prevent entanglement/ entrapment within cardiac structure or other device that may result in myocardial trauma and/or may require additional medical/surgical intervention.
- Warnings for patients with implantable pacemakers and implantable cardioverter/ defibrillators (ICDs):
 - Have temporary external sources of pacing and defibrillation available.
 - Fluoroscopic or appropriate imaging guidance and care must be taken during catheter advancement, manipulation, and withdrawal to avoid lead dislodgement.
 - Monitor pre- and post-measurements for sensing and pacing thresholds and impedances to determine the integrity of the lead-patient function.
 - Temporarily Deactivate ICDs for the duration of the procedure as they could discharge and injure the patient or be damaged by the procedure.
 - Perform a complete analysis of the implanted device function after ablation.
 - Remember to reactivate the pulse generator after the procedure.
- At no time should a BLAZER Dx-20 Catheter be advanced or withdrawn when resistance is felt without determining the cause. Valve damage, vascular and/or cardiac perforation is a risk with any intracardiac catheter.
- There are no data to support the safety and effectiveness of this device in the pediatric population.
- Excessive bending or kinking of the catheter shaft may damage internal wires. Manual prebending of the distal curve can damage the steering mechanism and/or electrical wires and may cause patient injury.

- Use both fluoroscopy or other visualization technique, such as echocardiography, and electrograms to monitor the advancement of the catheter to the area of the endocardium under investigation to avoid conduction pathway injury, cardiac perforation or tamponade.
- Guiding catheters and/or long introducer sheaths present the potential for thromboembolic events. Pre-flush and maintain lumen patency with heparinized intravenous infusion.
- Do not wipe this catheter with organic solvents such as alcohol or immerse the handle cable connector in fluids. This may result in electrical or mechanical catheter failures. It may also result in an allergic reaction from the patient.
- Catheter entrapment within the heart or blood vessels is a possible complication of cardiac catheterization procedures. The potential for catheter entrapment may be increased when the catheter is over torqued and/or positioned in the chordae tendineae. The occurrence of this complication may necessitate surgical intervention and/or repair of injured tissue and/or valve damage.
- This catheter has not been validated in an MRI environment. Do not use BLAZER Dx-20 Catheter near MRI equipment. Doing so may result in movement or heating of the catheter, which may cause distortion/loss of the catheter image on electroanatomic mapping. This could lead to patient injuries such as cardiac trauma.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the BLAZER Dx-20 catheter, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
- Any serious incident that has occurred in relation to the device should be reported to BSC and the competent authority of the Member State in which the user and/or patient is established.

PRECAUTIONS

- Before using, check shelf life. Do not use catheter after expiration date.
- Before using, inspect for physical damage, including electrical insulation on the cables and the catheter shaft. Replace damaged equipment.

ADVERSE EVENTS

Potential adverse events associated with the use of the BLAZER Dx-20 Catheter include, but are not limited to:

- Pain or discomfort, for example:
 - Angina
 - Chest pain
 - Non-cardiovascular pain
- Cardiac arrest
- Death
- Hypertension
- Hypotension
- Infection/inflammation/exposure to biohazardous material
- Edema/heart failure/pleural effusion
- Procedural related side effects, for example:
 - Allergic reaction (including anaphylaxis)
 - Genitourinary complication
 - Side effects related to medication or anesthesia
 - Radiation injury
 - Renal failure/insufficiency
 - Gastrointestinal disorders
 - Vasovagal response
- Respiratory distress/insufficiency/failure/dyspnea
- Arrhythmia (new or exacerbated)
 - Conduction pathway injury (for example: heart block, injury to sinus or AV node, etc.)
- Vessel trauma, including:
 - Perforation
 - Dissection
 - Coronary artery injury
 - Vasospasm
 - Occlusion
 - Hemothorax
- Cardiac trauma, for example:
 - Cardiac perforation/cardiac tamponade/pericardial effusion
 - Valvular damage

- Injury related to tissue damage and/or adjacent structures, for example:
 - Catheter entrapment
 - Physical trauma
- Surgical and access complications, for example:
 - Hematoma/seroma
 - AV fistula
 - Bleeding
 - Pseudoaneurysm
 - Pneumothorax
 - Thrombosis
- Injury due to embolism/thromboembolism/air embolism/foreign body embolism:
 - Cerebrovascular Accident (CVA)/stroke
 - Transient Ischemia Attack (TIA)
 - Myocardial infarction
 - Pulmonary embolism
 - Asymptomatic cerebral embolism
 - Neurological impairment and its symptoms, for example:
 - Cognitive changes, visual disturbances, headache, motor impairment, sensory impairment, and speech impairment

The potential adverse events may be related to the catheter(s) and/or the interventional procedure. The severity and/or the frequency of these potential adverse events may vary and may result in prolonged procedure time and/or additional medical and/or surgical intervention, implantation of a permanent device such as a pacemaker, and in rare cases, may result in death.

HOW SUPPLIED

Device Details

One diagnostic catheter supplied sterile by an ethylene oxide process.

Do not use if package is damaged or unintentionally opened before use.

Do not use if labeling is incomplete or illegible.

Handling and Storage

Do not use if the BLAZER Dx-20 Catheter is exposed to environmental conditions outside of the following ranges:

Operating Environment

Ambient Temperature: 10 °C to 40 °C

Relative Humidity: 30% to 75%

Atmospheric Pressure: 70 kPa to 106 kPa

Transport Environment

Temperature: -29 °C to 60 °C

Relative Humidity: 30% to 85%

Atmospheric Pressure: Uncontrolled

Storage Environment

Ambient Temperature: 20 °C to 30 °C

Relative Humidity: Uncontrolled

Atmospheric Pressure: Uncontrolled

OPERATIONAL INSTRUCTIONS

Additional Required Items

Intracardiac electrophysiology procedures should be performed in a specialized clinical setting, equipped with a fluoroscopy unit, radiographic table, physiologic recorder, emergency equipment, and instrumentation for gaining vascular access. Ancillary materials required to perform EP studies include:

1. 8F (2.67 mm) Sheath
2. ECG Recording System (see Figures 2 and 3)
3. Pacing Stimulator (see Figures 2 and 3)

Preparation

Carefully inspect the package prior to use for any breach of the sterile barrier or damage to the contents. If the sterile barrier integrity is compromised or the contents damaged, do not use and contact your Boston Scientific representative.

Before using, inspect for physical damage, including electrical insulation on the cables and the catheter shaft. Replace damaged equipment.

Procedure

1. Patient should be connected to an ECG recording system prior to commencing the procedure in order to permit arrhythmia monitoring.
2. Remove components from package and place in a sterile work area.
3. Create a vascular access by either cut down or percutaneous technique. The catheter may be used from femoral, brachial, subclavian, or jugular access sites through an introducer sheath (8F (2.67 mm)).
4. Connect cable to catheter by inserting the cable connector into the handle. Ensure that the cable/catheter connection remains dry throughout the procedure.
5. Advance the catheter under fluoroscopic guidance (or other visualization technique) to the desired intracardiac position.
6. Once inside the vessel, the tip can be deflected as necessary to facilitate the advancement.
7. For intracardiac recordings, connect the appropriate quick connect or lead tip pins of the cable to the appropriate electronic recording equipment per the recording equipment Operator's Manual and perform the study. Pins should be inserted into the appropriate receptacles on the Catheter Input Module (pin-out box) of the recording equipment.
Note: Closer spaced electrodes may enhance localization of aberrant pathways.
8. For pacing, connect the distal catheter terminal to the negative terminal of an external pulse stimulator and any proximal ring electrode to the positive terminal. The standard pacing stimuli is set at twice the diastolic threshold with a duration of 1.5 to 2 msec. If the pacing threshold is greater than 1 to 2 mA, repositioning of the catheter should be considered.
9. Prior to removing the catheter, ensure that the distal end of the catheter is straightened completely.
10. Withdraw the catheter from the vessel.
11. Remove the introducer sheath, and then follow standard practice for management of insertion site.

For questions regarding the use of this catheter, consult your local representative.

Disposal

To minimize risk of infection or microbial hazards after use, dispose of device and packaging as follows:

After use, the catheter may contain biohazardous substances. The catheter and packaging should be treated and disposed of as biohazardous waste in accordance with any applicable hospital, administrative, and/or local government regulations. Use of a biohazardous container with biological hazard symbol is recommended. Untreated biohazardous waste should not be disposed of in the municipal waste system.

Post-Procedure

In the event that a serious incident occurred in relation to the device, including all patient deaths for procedures where the BSC product was used, the event should be reported to BSC and the competent authority of the Member State or other relevant local authority in which the user and/or patient is established. Return any catheter related to a complaint, patient harm, injury, or death to Boston Scientific using a BSC Returned Product Kit.

Returning products for analysis and providing product performance observations helps drive reliability higher on an ongoing basis.

Be sure to follow the instructions with regard to packaging and shipping any biohazard devices taking care not to expose the handle or connector to fluid that may compromise the catheter and limit analysis.

INFORMATION TO BRIEF THE PATIENT

The physician should consider the following points while counseling patients on the use of the BLAZER Dx-20 Catheter in association with the Electrophysiological cardiac interventional procedure:

- Discuss the risks and benefits including review of potential adverse events listed in this document.
- Discuss post procedure instructions, including any lifestyle changes, medications, when to call the Healthcare Provider (HCP) and any post procedure follow-up that might be needed.


Figure 2. Diagnostic Catheter Connectivity Diagram A


Figure 3. Diagnostic Catheter Connectivity Diagram B

WARRANTY

For device warranty information, visit (www.bostonscientific.com/warranty).

EU Importer: Boston Scientific International B.V., Vestastraat 6, 6468 EX Kerkrade, The Netherlands

BLAZER is a trademark of Boston Scientific Corporation or its affiliates.

All other trademarks are the property of their respective owners.

SYMBOL DEFINITIONS

Commonly used medical device symbols that appear on the labeling are defined at www.bostonscientific.com/SymbolsGlossary.

Additional symbols are defined at the end of this document.

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All rights reserved.
Boston Scientific (Australia) Pty Ltd
© 2024 Boston Scientific Corporation or its affiliates.