Comet™ II

Pressure Guidewire

en	Instructions for Use	2
es	Instrucciones de uso	7
fr	Mode d'emploi	12
de	Gebrauchsanweisung	18
it	Istruzioni per l'uso	24
nl	Instructies voor gebruik	30
pt-EU	Instruções de Utilização	36

TABLE OF CONTENTS

REUSE WARNING3

DEVICE DESCRIPTION3

 Contents3

 Operating Principle.....3

 Materials.....3

 Non-Pyrogenic.....3

 User Information.....3

 Intended Patient Population3

INTENDED USE3

INDICATIONS FOR USE.....3

CONTRAINDICATIONS.....3

CLINICAL BENEFIT STATEMENT.....3

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE.....3

 Figure 1. Comet II Pressure Guidewire.....4

WARNINGS.....4

PRECAUTIONS.....4

ADVERSE EVENTS4

HOW SUPPLIED.....5

 Device Details5

 Handling and Storage.....5

OPERATIONAL INSTRUCTIONS5

 Preparation.....5

 Zeroing the Pressure Guidewire.....5

 Diagnostic Procedure.....5

 Intervention6

 How to Finish the Procedure.....6

 Disposal6

 Post-Procedure6

PATIENT COUNSELING INFORMATION.....6

 Technical Specifications6

WARRANTY.....6

SYMBOL DEFINITIONS.....6

 Table 2. Symbol Definitions42

Comet™ II

Pressure Guidewire

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

REUSE WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death.

Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

Carefully read all instructions prior to use. Observe all contraindications, warnings, precautions and adverse events noted in these directions. Failure to do so may result in patient complications. Boston Scientific relies on the physician to determine, assess, and communicate to each patient all foreseeable risks of the procedure.

DEVICE DESCRIPTION

The Comet II Pressure Guidewire is a single use, hydrophilic-coated, steerable wire with a custom optical cable. A pressure sensor is mounted approximately 3 cm from the distal end of the radiopaque and shapeable straight tip. There is a hydrophobic coating on the proximal end, and hydrophilic coating on the distal end and tip to facilitate lubricity and guidewire performance. For product specifications, including wire diameter, length, and radiopaque tip length, please refer to the product label.

Radial and femoral markers are located on the proximal segment of the pressure guidewire to aid in estimating the pressure guidewire's position relative to the distal guide catheter tip. The proximal markers are compatible with radial and femoral guide catheters that are at least 90 cm or 100 cm long, respectively.

Contents

Quantity	Material
(1)	Pressure Guidewire
(1)	Torque Device
(1)	Optical Cable

Operating Principle

The Comet II device utilizes a pressure sensor mounted within the distal hypotube to measure pressure gradients across stenoses to assess vascular physiology. Pressure readings are transmitted from the sensor to the proximal connector using an optical fiber. The device contains a radiopaque tip with radial and femoral proximal markers, along with lubricious hydrophobic and hydrophilic coatings to facilitate use as a coronary guidewire for interventional treatments.

Materials

Contains cobalt: CAS No. 7440-48-4; EC No. 231-158-0. Defined as a 1B carcinogen and reproductive toxicant according to the European Commission in a concentration above 0.1 % weight by weight.

Note: This device is made with a metal alloy which contains cobalt. Current scientific evidence supports that metal alloys containing cobalt used in medical devices, do not cause an increased risk of cancer or adverse reproductive effects.

Non-Pyrogenic

The Comet II Pressure Guidewire is non-pyrogenic.

User Information

Intended users of the pressure guidewire are primarily interventional cardiologists experienced in performing endovascular procedures in the sterile field. Catheterization lab technicians and/or nurses experienced in assisting endovascular procedures are also intended users for preparation and system utilization.

Intended Patient Population

The intended patient population are adult males and females undergoing intravascular interventional procedure(s) and/or measurement of intravascular physiological signals.

INTENDED USE

The Comet II Pressure Guidewire measures blood pressure gradient across coronary lesions during endovascular procedures. The Comet II Pressure Guidewire may also be used as a coronary guidewire for interventional treatments.

INDICATIONS FOR USE

The Comet II Pressure Guidewire is indicated to direct a catheter through a blood vessel and to measure physiological parameters in the coronary blood vessels.

CONTRAINDICATIONS

The Comet II Pressure Guidewire is contraindicated for use in the cerebral vasculature.

CLINICAL BENEFIT STATEMENT

To obtain physiological information about the blood pressure gradient to allow for physician informed treatment, optimizing outcomes for the patient during a percutaneous coronary interventions (PCI) procedure.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

For customers in the European Union, use the device name found in the labeling to search for the device's Summary of Safety and Clinical Performance, which is available on the European database on medical devices (Eudamed) website: <https://ec.europa.eu/tools/eudamed>.


Figure 1. Comet II Pressure Guidewire

WARNINGS

- The Comet II Pressure Guidewire should be used only by physicians experienced in angiography and percutaneous interventions including transluminal coronary angioplasty (PTCA). Vessel trauma (dissection, perforation, rupture or injury) may result from the improper use of this device.
- Resulting pressure guidewire fractures/separations might require additional percutaneous intervention or surgery.
- Use caution in patients for whom anticoagulation is not indicated.
- Severe reaction may occur in response to contrast agents that cannot be adequately premedicated.

PRECAUTIONS

1. Remove the pressure guidewire carefully from the carrier tube to reduce the possibility of damage to the distal tip. Refer to the Diagnostic Procedure section.
2. The pressure guidewire should be wiped clean only with heparinized saline prior to use to prevent thrombus accumulation.
3. Maintain control of the distal tip at all times during an intervention to avoid vessel dissections and perforations. The pressure guidewire should be manipulated in the body only under fluoroscopy. The pressure guidewire should never be advanced, removed, and/or torqued against resistance without first determining the reason for resistance. If ability to torque the pressure guidewire is diminished during procedure and the source of torque loss is not known, remove the guidewire without further torquing.
4. Exercise care in handling the pressure guidewire during a procedure to reduce the possibility of accidental breakage, bending, kinking, tip or wire separation or other damage to the pressure guidewire.
 - Do not use a pressure guidewire that has been damaged as it may decrease performance characteristics and/or lead to vessel damage.
 - Do not attempt to straighten a pressure guidewire that has been kinked or bent.
 - Do not leave the pressure guidewire in a prolapsed condition, as damage to the pressure guidewire may occur.
 - Wire damage may occur when the pressure guidewire is manipulated in a sharp bend, such as that caused by incomplete coaxial alignment of the delivery catheter and ostium, creating a sharp bend of the pressure guidewire between the delivery catheter and the vessel wall. Sharp bends are more likely to occur when using a less supportive delivery catheter, such as a diagnostic catheter.
5. When crossing a stent, exercise care to avoid entanglement between the pressure guidewire and the stent.

6. Avoid abrasion of the pressure guidewire coating.
 - To avoid damage to the hydrophilic coating, do not withdraw or manipulate the pressure guidewire in a metal cannula or sharp object.
 - Due to the variations of certain catheter tip inner diameters, abrasion of the hydrophilic coating may occur during manipulation. If any resistance is felt during introduction of the catheter, use a different catheter.
 - Excessive tightening of the torque device onto the pressure guidewire may result in abrasion of the coating on the pressure guidewire.

7. Exercise care in handling the optical cable during a procedure to reduce the possibility of accidental breakage, bending, kinking, or other damage.
8. Use only the optical cable provided to connect the pressure guidewire to FFR Link. Use of a different optical cable will produce inaccurate pressure readings.
9. Keep the proximal end of the pressure guidewire clean to ensure proper connection to the optical cable.
10. The accuracy of the diagnostic information is affected by, but not limited to:
 - Failure to achieve maximum coronary and myocardial hyperemia if using FFR (Fractional Flow Reserve) modality.
 - Interventional devices, such as balloon catheters, which are positioned so as to affect the blood flow or guidewires that stretch the vessel.
 - Pressure wire positioning relative to the lesion.
 - Microvascular resistance.
11. Carefully check and match therapeutic device compatibility to the pressure guidewire prior to use. Refer to instructions supplied with any interventional devices to be used in conjunction with pressure guidewire for their intended uses, contraindications and potential complications.
12. Do not use the pressure guidewire in conjunction with atherectomy catheters. It may damage the pressure guidewire.

ADVERSE EVENTS

Potential adverse events which may result from the use of the device include but are not limited to:

- Abrupt closure
- Angina or unstable angina
- Allergic reaction
- Arrhythmias
- Cardiac tamponade/pericardial effusion
- Contrast induced renal insufficiency or renal failure

- Death
- Embolism
- Exposure to biohazardous material
- Infection
- Myocardial infarction or ischemia
- Radiation injury
- Restenosis (reocclusion)
- Spasm
- Stroke/cerebral vascular accident (CVA)/transient ischemic attack (TIA)
- Vascular thrombus
- Vessel trauma (dissection, perforation, rupture or injury)

Some of the potential adverse events may require additional intervention including surgery.

HOW SUPPLIED

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

Device Details

Do not use if package is opened or damaged.

Do not use if labeling is incomplete or illegible.

Do not use if the use by date has been exceeded.

Handling and Storage

Transport and Storage Environment

Environment Temperature: -10 °C to 50 °C (+14 °F to 122 °F)

Relative Humidity: 30% RH to 85% RH

Atmospheric Pressure: 70 kPa to 106 kPa

OPERATIONAL INSTRUCTIONS

The following are required for use with the Comet II Pressure Guidewire:

- iLab, iLab Polaris System, AVVIGO branded Multi-Modality Guidance Systems, or an ANSI/AAMI-BP22 compatible hemodynamic receiving system
- FFR Link

Additional items for safe use:

- Introducer/Guide catheter of appropriate size and length
- Syringe for saline flushing of device
- Heparinized saline
- Hemostatic valve
- Guidewire insertion tool

Preparation

1. Ensure the receiving system and FFR Link are connected, powered on, and functioning properly.

Note: Refer to the receiving system User Guide throughout this procedure for instructions on using the system.

2. Remove the sterile pouch from the carton. Verify the device is suitable for use by checking the expiration date and ensuring the packaging and sterile barrier are intact. Do not use if the sterile package is damaged or the product is expired.

3. Use standard sterile technique to open the pouch and remove the tray lid. Leave the pressure guidewire in the carrier tube.
4. Flush the carrier tube only with heparinized saline by injecting into the hub end of the carrier tube to hydrate the hydrophilic segment of the pressure guidewire.

Zeroing the Pressure Guidewire

5. Place the carrier tube on a flat surface.
6. Insert the optical cable connector into the Pd IN port on the FFR Link. Upon recognizing the pressure guidewire, the FFR Link will perform an automatic zeroing.

Note: If the pressure guidewire is not recognized, attempt to disconnect and reconnect the optical cable connector to the FFR Link, ensuring it is fully inserted into the FFR link. In addition, ensure that the pressure guidewire is securely connected to the optical cable. If the system still will not recognize the pressure guidewire, use a new pressure guidewire.

Diagnostic Procedure

7. Carefully remove the pressure guidewire from the carrier tube and inspect it for damage.

Note: Do not use a damaged device. If the pressure guidewire cannot be removed easily, repeat the injection of heparinized saline into the hub end of the carrier tube and attempt to remove the pressure guidewire again.

8. If desired, the tip may be carefully shaped according to standard tip shaping practices. Do not use a shaping instrument with a sharp edge

Note: Take care while handling the distal tip during removal and shaping. Excessive handling may lead to sensor damage.

9. Open the hemostatic valve and the flush line of the coronary manifold. Insert a guidewire insertion tool through the valve and into the guide catheter.
10. Wipe the wire with heparinized saline prior to insertion into the insertion tool.
11. Carefully insert the distal tip of the pressure guidewire through the insertion tool and into the guide catheter.
12. Remove the insertion tool from the hemostatic valve and continue to advance the guidewire. Tighten the hemostatic valve knurled knob so that the valve seals around the guidewire, but does not inhibit intentional wire movement.

Note: The radial/femoral markers may be used to facilitate advancement.

13. Position the pressure guidewire sensor, which is 3 cm from the tip of the pressure guidewire, just distal to the tip of the guide catheter.
14. Flush the guide catheter with heparinized saline until all contrast is removed.
15. Ensure the aortic pressure transducer is at the height of the heart.
16. Allow the waveforms and trendlines to stabilize.
17. Use the system to equalize the pressure guidewire and aortic pressure transducer to 1.00 +/- 0.01. During equalization, do not disconnect the optical cable from the pressure guidewire or FFR Link.

Note: If the equalization is not successful, repeat the equalization procedure. Refer to receiving system Instructions for Use/User Manual for troubleshooting instructions.

18. Advance the pressure guidewire to the desired location for measurement.

Note: If a different tip shape or pressure guidewire is required, carefully withdraw and remove the pressure guidewire while observing pressure guidewire movement under fluoroscopy.

19. If using FFR modality, use standard technique to induce maximum hyperemic condition.
20. Begin recording pressure using the receiving system. If using FFR modality, record pressure until a steady state of maximum hyperemic condition is reached, or until the hyperemic effect starts to decrease.
21. Use the receiving system to stop recording, and review the data.
22. After every recording, withdraw the pressure guidewire to the equalization position, so that the sensor is just distal to the guide catheter tip, and recheck the pressure equalization. If necessary, repeat steps 17-21.
23. If additional measurements are to be taken from other locations, reposition the pressure guidewire and repeat steps 17-22.
24. If intervention is required, refer to the Intervention instructions below.

Intervention

The Comet II Pressure Guidewire or a different guidewire may be used for delivery of the interventional device. If using the pressure guidewire for the intervention, perform the following steps:

25. Disconnect the pressure guidewire from the optical cable handle.
26. Remove the torque device.
27. Clean the proximal end of the pressure guidewire with a saline-soaked gauze.
28. Carefully advance the interventional device over the pressure guidewire according to manufacturer's instructions.
29. Perform the intervention using standard practice according to manufacturer's instructions.
30. If additional measurements are needed following intervention, carefully clean the proximal end of the pressure guidewire and reconnect the pressure guidewire to the optical cable handle.
31. Verify proper connection by confirming a pressure signal has been achieved.
32. Return to step 17 above.

How to Finish the Procedure

33. When the procedure is complete, remove and discard the pressure guidewire and optical cable.
34. If desired, use the receiving system to review the physiologic data.

Disposal

To minimize the risk of infection or microbial hazards after use, dispose device and packaging as follows:

After use, device and packaging may contain biohazardous substances. Any device and packaging that came into contact with biohazardous substances should be treated and disposed of as biohazardous waste or be treated and disposed of in accordance with any applicable hospital, administrative, and/or local government regulations. Use of

a biohazardous container with biological hazard symbol is recommended. Untreated biohazardous waste should not be disposed of in the municipal waste system.

Post-Procedure

Any serious incident that occurs in relation to this device should be reported to the manufacturer and relevant local regulatory authority.

PATIENT COUNSELING INFORMATION

Ensure the patient is informed of the following important points:

- Contraindications, Warnings, Precautions and Adverse Events that are directly relevant to the patient.

Technical Specifications

Parameter	Value
Pressure Accuracy	+/-3% or +/-3 mmHg, whichever is greater
Zero Thermal Effect	0.3 mmHg/°C
Zero Drift	< 3 mmHg/H
Frequency Response	> 25 Hz
Operating Range	-45 mmHg to 300 mmHg

WARRANTY

For device warranty information, visit (www.bostonscientific.com/warranty).

The following are trademarks of Boston Scientific Corporation or its affiliates: Comet, iLab and AVVIGO. All other trademarks are the property of their respective owners.

SYMBOL DEFINITIONS

Commonly used medical device symbols that appear on the labeling are defined at www.bostonscientific.com/SymbolsGlossary

Additional symbols are defined at the end of this document.

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All rights reserved.
Boston Scientific Limited
Boston Scientific (Australia) Pty Ltd
© 2025 Boston Scientific Corporation or its affiliates.