RX Cytology Brush

BILIARY

Wireguided Cytology Brush

REF M00545000

en	Instructions for Use	2
es	Instrucciones de uso	6
fr	Mode d'emploi	10
de	Gebrauchsanweisung	14
it	Istruzioni per l'uso	19
nl	Instructies voor gebruik	23
pt-EU	Instruções de Utilização	27

51596217-01

2023-04

RX Cytology Brush

BILIARY

Wireguided Cytology Brush

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

REUSE WARNING

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

DEVICE DESCRIPTION

Contents

The package includes:

One (1) RX Biliary Wireguided Cytology Brush

Operating Principle

The RX Biliary Wireguided Cytology Brush (Figure 1) is a 8F (2.7 mm) dual-lumen, wire-guided catheter with a 2.1 mm cytology brush. The guidewire and brush both exit the catheter at the distal tip. The brush-lumen may also be used to inject contrast media. The RX Biliary Wireguided Cytology Brush is capable of accepting a 0.035 in (0.89 mm) Boston Scientific guidewire. The open channel allows for the quick exchange of a guidewire, which is completely isolated from injection agents. The device may be placed with or without the aid of a guidewire. The distal radiopaque marker aids in the catheter visualization during fluoroscopy/cholangiography.


Figure 1

Visual Markers	Distance From Tip
Brown band	5 cm (50 mm)
Blue band	10 cm (100 mm)
Green band	15 cm (150 mm)

Catheter Shaft

The guidewire channel of the catheter accommodates a 0.035 in (0.89 mm) Boston Scientific guidewire through the length of the RX Biliary Wireguided Cytology Brush. Contrast medium injection may be achieved through the proximal lumen by attaching a contrast-filled Luer-Lok Syringe to the proximal luer. The RX Biliary Wireguided Cytology Brush may be placed with or without a guidewire.

Catheter Tip

The distal catheter contains a radiopaque (RO) marker. This RO marker can be used to fluoroscopically monitor and/or confirm position.

Materials

Contains Cobalt: CAS No. 7440-48-4; EN No. 231-158-0.
Defined as a 1B carcinogen and reproductive toxicant according to the European Commission in a concentration above 0.1% weight by weight.

Note: Current scientific evidence supports that metal alloys containing cobalt used in medical devices, do not cause an increased risk of cancer or adverse reproductive effects.

Compatibility

The RX Biliary Wireguided Cytology Brush is designed for use with an endoscope working channel that has a diameter of at least 2.8 mm. Due to various tolerances among manufacturers, it is recommended that the RX Biliary Wireguided Cytology Brush be used with a Boston Scientific 0.035 in (0.89 mm) guidewire. Any use of this device, other than those indicated in these instructions, is not recommended.

User Information

The RX Biliary Wireguided Cytology Brush should only be used by or under the supervision of physicians thoroughly trained in endoscopic biliary procedures. A thorough understanding of the technical principles, clinical applications, and risks associated with ERCP and ES is necessary before using this device.

INTENDED USE/INDICATIONS FOR USE

The RX Biliary Wireguided Cytology Brush is used endoscopically to collect cells and/or to aid in the delivery of contrast during a cholangiogram in the biliary system.

Clinical Benefit Statement

The clinical benefit of the RX Biliary Wireguided Cytology Brush is to allow for endoscopic collection of tissue to aid in the diagnosis of pathology in the biliary tree or to aid in the delivery of contrast during a cholangiogram.

CONTRAINDICATIONS

Contraindications for this device are those specific to endoscopic retrograde cholangiopancreatography (ERCP) and endoscopic sphincterotomy (ES).

WARNINGS

- While brushing, monitor the catheter position with contrast media injection and/or fluoroscopy in order to avoid patient injury.
- Any use of this device, other than those indicated in these instructions, is not recommended and could result in patient injury.

ADVERSE EVENTS

Potential complications include, but may not be limited to:

- Allergic reaction to contrast medium
- Cholangitis
- Hematoma
- Hemorrhage
- Pancreatitis
- Perforation
- Septicemia/Infection

HOW SUPPLIED

The RX Biliary Wireguided Cytology Brush is supplied sterile using an ethylene oxide (EO) process. The packaging and device should be inspected prior to use.

Device Details

Do not use if package is damaged or unintentionally opened before use.

Do not use if labeling is incomplete or illegible.

Handling and Storage

This product has no special handling or storage requirements. Check device label for expiration date.

OPERATIONAL INSTRUCTIONS

Preparation

Precaution: Before using, inspect the package for any breach to confirm the product’s sterility. If the package has been opened, do not use the catheter.

1. Remove the RX Biliary Wireguided Cytology Brush from its package and inspect for any damage.
2. If contrast medium will be injected, flush catheter with medium to remove all bubbles before passing catheter through scope.
3. Remove protector prior to placement.
4. Ensure a RX Locking Device has been attached to the scope.

Procedure

1. Confirm guidewire placement, both endoscopically and fluoroscopically, above the stricture, stone or injury.
2. After a Boston Scientific 0.035 in (0.89 mm) guidewire has been positioned and locked into place, back load the RX Biliary Wireguided Cytology Brush over the proximal end of the guidewire. Slide the RX Biliary Wireguided Cytology Brush onto the guidewire and through the first 20-30 cm (200-300 mm) until the guidewire exits the distal U-channel port.
3. Grasp the guidewire and slide the catheter over the entire length of the guidewire until it reaches the guidewire locking device. Unlock the guidewire from the guidewire locking device. Insert the 20-30 cm (200-300 mm) section using the standard exchange method. Lock guidewire into place and continue advancing the RX Biliary Wireguided Cytology Brush through the scope until cannulation is achieved.
4. Following cannulation, inject contrast medium through the injection port if contrast medium is being used, and fluoroscopically confirm position in the common bile duct.
5. Position the RX Biliary Wireguided Cytology Brush at the desired brushing site using fluoroscopic guidance.
6. Brush the site several times. After the brushing is complete retract the brush into the catheter by inserting thumb into thumb ring and pulling the thumb ring back relative to the handle to protect sample.
7. The catheter can then be removed.
8. Prepare specimen for examination per institutional guidelines.

Device Removal

The RX Biliary Wireguided Cytology Brush may be removed with or without leaving the guidewire in place. If the guidewire is to remain in place while the RX Biliary Wireguided Cytology Brush is withdrawn, utilize the following steps:

1. Place guidewire into locking channel of the Guidewire Locking Device.
2. Withdraw the RX Biliary Wireguided Cytology Brush from the scope until the distal U-channel marker is visible at the Guidewire Locking and resistance is felt.

3. Perform standard exchange procedure over the final 20-30 cm (200-300 mm) of the RX Biliary Wireguided Cytology Brush. Once the distal blue tip is exposed from the Guidewire Locking Device, re-lock the guidewire into the Guidewire Locking Device and pull the RX Biliary Wireguided Cytology Brush off the wire.

Disposal

To minimize the risk of infection or microbial hazards after use, dispose device and packaging as follows:

After use, device may contain biohazardous substances. The device and packaging should be treated and disposed of as biohazardous waste or have them treated and disposed of in accordance with any applicable hospital, administrative, and/or local government regulations. Use of a biohazardous container with biological hazard symbol is recommended. Untreated biohazardous waste should not be disposed of in the municipal waste system.

This device contains a sharps hazard. Take precautions to ensure that sharps are handled properly. Dispose of all sharps directly into a sharps disposal container labeled with a biological hazard symbol. Sharps waste should be safely disposed of using available sharps waste channels in accordance with hospital, administrative, and/or local government policy.

Post-Procedure

Any serious incident that occurs in relation to this device should be reported to the manufacturer and relevant local regulatory authority.

PATIENT COUNSELING INFORMATION

Brief the patient on any relevant post-procedure instructions, contraindications, warnings, precautions and/or adverse events found within this Instructions for Use (IFU), pertaining to the patient.

WARRANTY

For device warranty information, visit (www.bostonscientific.com/warranty).

All other trademarks are the property of their respective owners.

SYMBOL DEFINITIONS

Commonly used medical device symbols that appear on the labeling are defined at www.bostonscientific.com/SymbolsGlossary.

Additional symbols are defined at the end of this document.

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or its affiliates. All rights reserved.
Boston Scientific (Australia) Pty Ltd
Boston Scientific Limited
Boston Scientific Corporation
© 2023 Boston Scientific Corporation