WallFlex™ Colonic

Stent System with Anchor Lock
Delivery System

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2025-04

TABLE OF CONTENTS

REUSE WARNING	3
DEVICE DESCRIPTION	3
Contents	3
Operating Principle	3
Figure A	4
Materials	4
User Information	4
Instructions for Use	4
INDICATIONS FOR USE / INTENDED USE	4
Clinical Benefit Statement	4
Summary of Safety and Clinical Performance	4
CONTRAINDICATIONS	5
WARNINGS	5
PRECAUTIONS	5
MR CONDITIONAL	5
ADVERSE EVENTS	6
HOW SUPPLIED	6
Device Details	6
Handling and Storage	6
EQUIPMENT REQUIRED	6
OPERATIONAL INSTRUCTIONS	6
Procedure	6
Figure B	7
Figure C	8
Figure D	8
Disposal	9
Post-Procedure	9
PATIENT COUNSELING INFORMATION	9
Implantable Device Patient Information	9
Implant Card Instructions	9
WARRANTY	9
SYMBOL DEFINITIONS	9
Table 1. Symbol Definitions	9

WallFlex™ Colonic

Stent System with Anchor Lock Delivery System

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

REUSE WARNING

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious diseases) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

DEVICE DESCRIPTION

Contents

The package includes:

One (1) WallFlex Colonic Stent System with Anchor Lock Delivery System.

Operating Principle

The WallFlex Colonic Stent System with Anchor Lock Delivery System consists of two components: the implantable metallic stent and the Anchor Lock Delivery System (See Figure A). The stent is composed of nitinol wire braided in a tubular mesh configuration. This design configuration results in a stent that is both flexible and self-expanding. The delivery system consists in part of coaxial tubes. The exterior tube serves to constrain the stent until retracted during deployment. The interior tube contains a central lumen that accommodates a 0.035 in (0.89 mm) guidewire. The device may be inserted through the working channel of an endoscope (minimum channel diameter 3.7 mm).

There are 3 radiopaque marker bands; The Exterior Tube Marker Band, the Deployment Limit Marker Band and the Post Deployment Marker Band. The Exterior Tube Marker Band is located adjacent to the leading end of the stent and the Deployment Limit Marker Band is located adjacent to the trailing end of the stent. The Deployment Limit Marker Band identifies the stent deployment limit, the point beyond which the stent cannot be reconstrained. The Post Deployment Marker Band helps to facilitate accurate placement of the stent as it represents the approximate final position of the trailing end of the stent when fully expanded.


Figure A

Materials

Implantable Material	% Weight
Nitinol	100

Warning: The stent contains nickel, which may cause an allergic reaction in individuals with nickel sensitivity.

User Information

The WallFlex Colonic Stent System with Anchor Lock Delivery System should only be used by or under the supervision of physicians thoroughly trained in the placement of colonic stents. A thorough understanding of the techniques, principles, clinical applications and risks associated with this procedure is necessary before using the device. Read the entire instructions for use thoroughly before using the WallFlex Colonic Stent System with Anchor Lock Delivery System.

Instructions for Use

(See Figure A for Device Diagram)

The exterior tube is pulled back by immobilizing the hub handle in one hand, grasping the exterior tube handle with the other hand, and gently sliding the exterior tube handle along the stainless-steel tube towards the hub handle. Retraction of the exterior tube releases the stent.

The deployment process can be reversed if repositioning is desired. The stent can be reconstrained by the exterior tube if the stent Deployment Limit Marker Band has not been passed during deployment. To reposition, first fully reconstrain the stent by gently pushing the exterior tube handle forward while holding the hub handle completely immobile. Once reconstrained, the stent can be repositioned either distally or proximally and the deployment process restarted. Reversing the deployment process can be repeated twice, allowing a total of three deployment attempts.

INDICATIONS FOR USE / INTENDED USE

The device is indicated for the palliative treatment of Colonic strictures produced by malignant neoplasm and to relieve large bowel obstruction prior to colectomy in patients with malignant strictures.

Clinical Benefit Statement

To treat colonic strictures produced by malignant neoplasm and relieve large bowel obstruction prior to colectomy in patients with malignant strictures.

Summary of Safety and Clinical Performance

For customers in the European Union, use the device name found in the labeling to search for the device's Summary of Safety and Clinical Performance, which is available on the European database on medical devices (EUDAMED) website: (<https://ec.europa.eu/tools/eudamed>).

CONTRAINDICATIONS

Contraindications associated with the use of the WallFlex Colonic Stent System with Anchor Lock Delivery System include:

- Enteral ischemia
- Suspected or impending perforation
- Intra-abdominal abscess/perforation
- Strictures that do not allow passage of a guidewire
- Patients for whom endoscopic techniques are contraindicated
- Any use other than those specifically outlined under indications for use

WARNINGS

- As perforation is a known risk, the stent should be used with caution and only after careful consideration in patients who are:
 - undergoing radiation therapy and/or chemotherapy
 - in advanced stages of cancer
- The device should be used with caution and only after careful consideration in patients with elevated bleeding times, coagulopathies, or in patients with radiation colitis or proctitis as this may result in bleeding.
- Stents cannot be repositioned after complete deployment as this may result in a patient injury such as perforation.
- The safety and effectiveness of this device for use in benign strictures have not been established.
- Chemoradiation therapy or radiotherapy alone may lead to tumor shrinkage and subsequent stent migration.
- To minimize pain and tenesmus, the proximal stent end should be placed 2 cm above the anal canal or 6 cm from the anus.

PRECAUTIONS

- The device is intended for single use only. Do not attempt to reload deployed stents onto the delivery system.
- Use of fluoroscopy is recommended. Not using fluoroscopy can result in misplacement of the stent.

MR CONDITIONAL

 MRI Safety Information A person with the <i>WallFlex Stent(s)</i> may be safely scanned under the following conditions. Failure to follow these conditions may result in injury.	
Device Name	<i>WallFlex Stent</i>
Static magnetic Field Strength (B ₀)	1.5T or 3.0T
Maximum Spatial Field Gradient	30 T/m (3,000 gauss/cm)
RF Excitation	Circularly Polarized (CP)
RF Transmit Coil Type	Cylindrical Whole-body Coil Cylindrical Head Coil
Operating Mode	Normal Operating Mode
Maximum Whole-Body SAR	2 W/kg (Normal Operating Mode)
Maximum Head SAR	3.2 W/kg (Normal Operating Mode)
Scan Duration	Under scan conditions defined above, <i>WallFlex Stent</i> can be used for 60 minutes of continuous RF (a sequence or back-to-back series/scan without breaks)

 MRI Safety Information A person with the <i>WallFlex Stent(s)</i> may be safely scanned under the following conditions. Failure to follow these conditions may result in injury.	
MR Image Artifact	Image artifact caused by device may extend approximately 20 mm from the perimeter and 5 mm from the end of the stent with a spin and gradient echo pulse sequence

ADVERSE EVENTS

Complications associated with the use of the WallFlex Colonic Stent System with Anchor Lock Delivery System may include:

- Allergic Reaction
- Death (other than that due to normal disease progression)
- Diarrhea
- Erosion (e.g. ulceration)
- Gastrointestinal Symptoms
- Hemorrhage
- Infection
- Inflammation
- Migration
- Obstruction
- Pain (e.g. discomfort, foreign body sensation)
- Perforation
- Peritonitis
- Stent misplacement or inadequate expansion
- Tenesmus (e.g. symptoms of tenesmus or urgency/incontinence)
- Tissue Damage
- Tissue Overgrowth/Ingrowth

HOW SUPPLIED

Device is supplied sterile using an ethylene oxide (EO) process. The packaging and device should be inspected prior to use.

Device Details

Do not use if package is damaged or unintentionally opened before use.

Do not use if labeling is incomplete or illegible.

Handling and Storage

This product has no special handling and storage requirements. Check device label for expiration date.

EQUIPMENT REQUIRED

- Endoscope (with a channel size of 3.7 mm or greater)
- Fluoroscope
- 0.035 in (0.89 mm) guidewire

OPERATIONAL INSTRUCTIONS

Procedure

The WallFlex Colonic Stent System with Anchor Lock Delivery System can be placed with the aid of fluoroscopy, or under direct visualization alongside an endoscope or the combination of both.

1. Prepare and sedate the patient as per endoscopy protocol for the site.

A. Fluoroscopy Procedure

Pass a 0.035 in (0.89 mm) guidewire to the level of the stricture. The guidewire is maneuvered until the wire traverses the obstructed area.

B. Endoscopic Procedure

Pass an endoscope to the site of the stricture. Under direct visualization pass a 0.035 in (0.89 mm) guidewire through the working channel of the endoscope and maneuver the guidewire across the site of the stricture.

2. Identify the stricture location, and the proximal and distal stricture margins, using fluoroscopic and/or endoscopic visualization. Fluoroscopic markers (e.g. external metal markers taped to the skin or internal injection of radiopaque dye) may be used to indicate the stricture margins.
3. Estimate stricture length by one of several methods (e.g. withdrawal of catheter with radiopaque marker at its tip or withdrawal of marked guidewire).
4. Allowing for possible further tumor development and post implant stent shortening (due to continued expansion), determine the length and number of stents necessary to cross the stricture. (Note: The fully expanded stent length should be at least 4 cm longer than the length of the stricture.) Should multiple stents be required, the leading stent should be placed first followed by the trailing stent, allowing for adequate overlapping.

If the obstruction happens to fall near a bend, it is recommended that the stent(s) total length should be such that the two free ends fall beyond the bend and open freely into the lumen.

5. Carefully remove the delivery system from its protective packaging.

Precaution: Visually inspect the entire system for damage. Visually check that the exterior tube covers the leading end of the stent. Ensure that the delivery system does not have any kinked sections.

6. Pass the WallFlex Colonic Stent System with Anchor Lock Delivery System over the guidewire, through the endoscope working channel and to the site of the stricture.
7. Advance the stent across the stricture until the Post Deployment Marker Band is at the outermost point of the proximal end of the stricture. (See Figure B).

Note: Accounting for maximum shortening, the Post Deployment Marker Band identifies the approximate furthest position of the proximal end of the deployed stent. Because shortening will vary based on anatomy, the Post Deployment Marker Band should only be used to approximate the final proximal end of the stent. The Exterior Tube Marker Band should be at least 4 cm beyond the distal end of the stricture (See Figure B). If the Exterior Tube Marker Band is not at least 4 cm beyond the distal end of the stricture, a longer stent may be required or a second overlapping stent may be used to cross the stricture adequately.


Figure B

8. To begin stent deployment, immobilize the hub handle in one hand and grasp the exterior tube handle with the other hand. Gently slide the exterior tube handle back along the stainless-steel tube towards the hub handle. At the same time and the same rate, the endoscopist should gently withdraw the exterior tube to prevent inadvertent stent advancement.

Note: The stent can be reconstrained at any time up to when the Exterior Tube Marker Band reaches the Deployment Limit Marker band (See Figure C).


Figure C

Warning: Do not push forward or pull backward on the hub handle with the stent partially deployed. The hub handle must be securely immobilized. Inadvertent movement of the hub handle may cause misalignment of the stent and possible intestinal wall damage.

Warning: The stent should deploy easily. Do not deploy the stent if unusual force is required. Using such force may cause a patient injury such as perforation, bleeding, or tissue damage. To reconstrain and remove the device, see the Caution in Step 10 below.

Precaution: Do not reconstrain around tortuous anatomy as it may cause damage to the device.

9. Assess stent position and reposition if desired. To reposition, first fully reconstrain the stent by gently pushing the exterior tube handle while holding the hub handle completely immobile. When fully reconstrained, the delivery system can be moved either proximally or distally and the deployment process restarted. Repositioning can be repeated twice, allowing a total of three deployment attempts.

Precaution: A stent cannot be repositioned after the Deployment Limit Marker Band has been passed.

10. To complete stent deployment, immobilize the hub handle with one hand, grasp the exterior tube handle with the other hand, and gently slide the exterior tube handle along the stainless-steel tube towards the hub handle until the stent is fully deployed (See Figure D). As before, the endoscopist should gently withdraw the exterior tube at the same time and at the same rate to prevent inadvertent stent advancement.


Figure D

Precaution: If the stent cannot be fully deployed or fully reconstrained, the entire delivery system can be pulled into the endoscope, using the working channel of the endoscope to reconstrain the stent. This may cause damage to the working channel of the endoscope. In addition, the stent could inadvertently be forced off of the delivery system.

11. After the stent is correctly positioned and fully deployed, the delivery system should be closed, by pushing the exterior tube handle forward, and removed. The implanted stent length should allow for adequate overlapping into the non-obstructed anatomy to compensate for further tumor progression and stent shortening. In the event the stent does not adequately

cover the stricture, a second stent should be implanted providing adequate overlapping of the initially placed stent.

Precaution: Do not place stents end to end (without overlap) as this may cause kinking. Do not use in combination with stents from other manufacturers.

Precaution: If the delivery system is not fully closed prior to withdrawal, there is a possibility that the tip of the delivery system will get caught in the stent.

12. Perform routine post implant procedures.

Warning: Do not dilate the stent after placement, this may result in perforation.

Disposal

To minimize the risk of infection or microbial hazards after use, dispose device and packaging as follows:

After use, device may contain biohazardous substances. The device and packaging should be treated and disposed of as biohazardous waste or have them treated and disposed of in accordance with any applicable hospital, administrative, and/or local government regulations. Use of a biohazardous container with biological hazard symbol is recommended. Untreated biohazardous waste should not be disposed of in the municipal waste system.

Post-Procedure

Any serious incident that occurs in relation to this device should be reported to the manufacturer and relevant local regulatory authority.

PATIENT COUNSELING INFORMATION

Brief the patient on any relevant post-procedure instructions, contraindications, warnings, precautions and/or adverse events found within this Instructions for Use (IFU), pertaining to the patient.

Implantable Device Patient Information

Advise the patient that additional information may be available to them on the Boston Scientific website (www.bostonscientific.com/patientlabeling).

Implant Card Instructions

- Apply the peel-off label from the product onto the supplied patient implant card.
- Fill in the implant date, patient name, healthcare institution and/or physician information.
- Give the completed card to the patient.

WARRANTY

For device warranty information, visit (www.bostonscientific.com/warranty).

The following are trademarks of Boston Scientific Corporation or its affiliates: WallFlex.

All other trademarks are the property of their respective owners.

SYMBOL DEFINITIONS

Commonly used medical device symbols that appear on the labeling are defined at www.bostonscientific.com/SymbolsGlossary.

Table 1. Symbol Definitions

	Contents
---	----------

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or its affiliates. All rights reserved.
Boston Scientific (Australia) Pty Ltd
Boston Scientific Limited
© 2025 Boston Scientific Corporation