# Ranger™

**MONORAIL™**

## Paclitaxel-Coated PTA Balloon Catheter

|       |                          |    |
|-------|--------------------------|----|
| en    | Instructions for Use     | 2  |
| es    | Instrucciones de uso     | 8  |
| fr    | Mode d'emploi            | 14 |
| de    | Gebrauchsanweisung       | 20 |
| it    | Istruzioni per l'uso     | 26 |
| nl    | Instructies voor gebruik | 32 |
| pt-EU | Instruções de Utilização | 38 |

## TABLE OF CONTENTS

|                                                                |   |
|----------------------------------------------------------------|---|
| <b>WARNING</b>                                                 | 3 |
| <b>REUSE WARNING</b>                                           | 3 |
| <b>DEVICE DESCRIPTION</b>                                      | 3 |
| Materials                                                      | 3 |
| Table 1: Paclitaxel Dose per Balloon                           | 3 |
| Contents                                                       | 3 |
| Operating Principle                                            | 3 |
| Non-pyrogenic                                                  | 3 |
| User Information                                               | 3 |
| <b>INTENDED USE / INDICATIONS FOR USE</b>                      | 3 |
| <b>CLINICAL BENEFIT STATEMENT</b>                              | 3 |
| <b>SUMMARY OF SAFETY AND CLINICAL PERFORMANCE</b>              | 3 |
| <b>CONTRAINDICATIONS</b>                                       | 3 |
| <b>WARNINGS</b>                                                | 3 |
| Pre- and Post-Procedure Antiplatelet Therapy                   | 4 |
| <b>PRECAUTIONS</b>                                             | 4 |
| Pregnancy / Lactation                                          | 4 |
| Drug Information                                               | 4 |
| Drug Interaction                                               | 4 |
| Carcinogenicity, Genotoxicity, and Reproductive Toxicology     | 4 |
| <b>ADVERSE EVENTS</b>                                          | 4 |
| <b>CLINICAL STUDIES</b>                                        | 4 |
| <b>HOW SUPPLIED</b>                                            | 5 |
| <b>Device Details</b>                                          | 5 |
| Handling and Storage                                           | 5 |
| <b>OPERATIONAL INSTRUCTIONS</b>                                | 5 |
| Additional Items for Safe Use                                  | 5 |
| Inspection Prior to Use                                        | 5 |
| <b>PREPARATION</b>                                             | 5 |
| Inflation Device Preparation                                   | 5 |
| <b>Balloon Catheter Selection</b>                              | 5 |
| <b>Balloon Catheter Preparation</b>                            | 5 |
| <b>PROCEDURE</b>                                               | 5 |
| Insertion Procedure                                            | 5 |
| Figure 1: Inserting the Balloon Catheter with the Loading Tool | 6 |
| Table 2: Typical Ranger Balloon Catheter Compliance            | 7 |
| Disposal                                                       | 7 |
| Post-Procedure                                                 | 7 |
| <b>PATIENT COUNSELING INFORMATION</b>                          | 7 |
| <b>WARRANTY</b>                                                | 7 |
| <b>SYMBOL DEFINITIONS</b>                                      | 7 |
| Table 3. Symbol Definitions                                    | 7 |

# Ranger™

## MONORAIL™

### Paclitaxel-Coated PTA Balloon Catheter

#### Rx ONLY

**Caution:** Federal Law (USA) restricts this device to sale by or on the order of a physician.

#### WARNING

A signal for increased risk of late mortality has been identified following the use of paclitaxel-coated balloons and paclitaxel-eluting stents for femoropopliteal arterial disease beginning approximately 2-3 years post-treatment compared with the use of non-drug coated devices. There is uncertainty regarding the magnitude and mechanism for the increased late mortality risk, including the impact of repeat paclitaxel-coated device exposure. Physicians should discuss this late mortality signal and the benefits and risks of available treatment options with their patients.

#### REUSE WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and / or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and / or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and / or local government policy. See section "Disposal" for additional details regarding disposal of this medical device.

#### DEVICE DESCRIPTION

The Ranger Paclitaxel-Coated PTA Balloon Catheter is a Monorail (MR) Percutaneous Transluminal Angioplasty (PTA) Balloon Catheter with a semicompliant balloon coated with a formulation of paclitaxel (drug) and an excipient. The Ranger Drug-Coated Balloon (DCB) is designed to inhibit restenosis by delivering drug to diseased arterial tissue.

#### Materials

##### Device Component Description

The Ranger DCB has a coaxial shaft design with a Nylon / Polyamide outer shaft. The outer lumen is used for inflation of the balloon, and the wire lumen permits the use of guidewires 0.014 in (0.36 mm) or 0.018 in (0.46 mm) to facilitate advancement of the catheter. The Polyether / Polyamide balloon is designed to provide an inflatable segment of known diameter and length at recommended pressures. The catheter includes a tapered tip to facilitate advancement of the catheter to and through the stenosis.

The Ranger DCB has two radiopaque marker bands (one proximal and one distal) which, in conjunction with fluoroscopy, aid in the placement of the balloon. The proximal shaft marks define the length to the distal end of the catheter.

The working length of the Ranger DCB is 200 cm. The 200 cm working length catheter has one mark at 90 cm and two marks at 100 cm.

The proximal portion of the Ranger DCB includes one female luer-lock port connected to the inflation lumen. In addition, the Ranger DCB is equipped with a loading tool to help protect the drug coating as it enters the hemostatic valve.

##### Drug Component Description

The Ranger DCB has a drug coating formulation consisting of a 4:1 ratio of paclitaxel (the active pharmaceutical ingredient) and an excipient (the inactive ingredient).

##### Paclitaxel

The active pharmaceutical ingredient in the balloon coating is paclitaxel (PTx).

The dose density of paclitaxel is 2.0 µg per mm² of the balloon surface area. The principal mechanism by which PTx inhibits neointimal growth is through the stabilization of microtubules by preventing depolymerization during the final G2/M phase of cell division. PTx is a semi-synthetic paclitaxel synthesized from precursor compounds, isolated from a spectrum of Taxus species and hybrids.

Paclitaxel is a diterpenoid with a characteristic taxane skeleton of 20 carbon atoms, a molecular weight of 853.91 g/mol and a molecular formula of C<sub>47</sub>H<sub>51</sub>NO<sub>14</sub>. It is highly lipophilic, insoluble in water, but freely soluble in methanol, ethanol, chloroform, ethyl acetate and dimethyl sulfoxide.

#### Acetyl Tributyl Citrate

The coating utilizes the inactive ingredient Acetyl Tributyl Citrate (ATBC) as an excipient to facilitate the release and transfer of paclitaxel into the arterial wall.

ATBC is a carboxylic acid ester with a molecular weight of 402.48 g/mol.

Table 1: Paclitaxel Dose per Balloon

|                       |   | Balloon Length (mm) |         |         |         |         |         |         |
|-----------------------|---|---------------------|---------|---------|---------|---------|---------|---------|
|                       |   | 40                  | 60      | 80      | 100     | 120     | 150     | 200     |
| Balloon Diameter (mm) | 4 | 1043 µg             | 1564 µg | 2086 µg | 2607 µg | 3100 µg | 3878 µg | 5165 µg |
|                       | 5 | 1301 µg             | 1951 µg | 2601 µg | 3251 µg | 3978 µg | 4942 µg | 6592 µg |
|                       | 6 | 1467 µg             | 2200 µg | 3108 µg | 3889 µg | 4809 µg | 5975 µg | 7966 µg |
|                       | 7 | 1711 µg             | 2714 µg | 3626 µg | 4538 µg | 5564 µg | 6958 µg | 9266 µg |

#### Contents

| Qty | Material                                  |
|-----|-------------------------------------------|
| 1   | Ranger Drug-Coated Balloon (DCB) Catheter |

#### Operating Principle

The Ranger Drug Coated PTA Balloon Catheter consists of a balloon catheter coated with paclitaxel as the active pharmaceutical ingredient and ATBC as the excipient. The Ranger DCB consists of a shaft with a balloon near the distal tip.

The balloon is inflated to expand the luminal diameter of a diseased vessel and deliver the drug coating to reduce the rate of neointimal proliferation of the vessel wall.

The paclitaxel coating is distributed evenly across the working length of the balloon.

The ATBC excipient facilitates the release and transfer of the active pharmaceutical ingredient into the arterial vessel wall.

#### Non-pyrogenic

This device meets pyrogen limit specifications.

#### User Information

The balloon catheter should be used only by physicians experienced in the performance of percutaneous transluminal angioplasty.

#### INTENDED USE / INDICATIONS FOR USE

The Ranger DCB is indicated for percutaneous transluminal angioplasty (PTA) in adults with de novo or restenotic lesions located in native superficial femoral and proximal popliteal arteries (SFA / PPA), and infrapopliteal arteries.

#### CLINICAL BENEFIT STATEMENT

The Ranger DCB is designed to treat peripheral artery disease and improve lumen diameter through percutaneous transluminal angioplasty (PTA) in the target lesion. The clinical benefit can be measured by overall clinical outcomes as demonstrated by primary patency rates, freedom from amputation and freedom from TLR, among others, with respect to other existing therapies.

#### SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

For customers in the European Union, use the device name found in the labeling to search for the device's Summary of Safety and Clinical Performance, which is available on the European database on medical devices (EUDAMED) website:

(<https://ec.europa.eu/tools/eudamed>).

#### CONTRAINDICATIONS

Use of the Ranger DCB is contraindicated in patients with known hypersensitivity to paclitaxel (or structurally-related compounds).

#### WARNINGS

- This product should not be used in patients with uncorrected bleeding disorders or patients who cannot receive anticoagulation or antiplatelet aggregation therapy.
- Never use air or any gaseous medium to inflate the balloon (a 50:50 mixture of contrast medium and sterile saline is recommended).
- When the balloon catheter is exposed to the vascular system, it should be manipulated while under high-quality fluoroscopic observation.

- Do not advance or retract the catheter unless the balloon is fully deflated under vacuum.
- If unusual resistance is felt during manipulation, determine the cause of the resistance before proceeding. If the source of resistance cannot be determined, it is recommended to extract the entire system with the guide or introducer sheath.
- If difficulty is experienced during balloon inflation, do not continue. Deflate the balloon and remove the catheter.
- To reduce the potential for vessel damage, the inflated diameter of the balloon should approximate the diameter of the vessel segment to be treated. The inflated length of the balloon (shoulder to shoulder) may exceed the length of the lesion / stenosis by approximately 10 mm on either side within the targeted artery.
- Do not exceed the balloon rated burst pressure. Use of an inflation device is recommended to prevent over-pressurization.
- If difficulty is experienced removing the balloon catheter, other retrieval methods (use of additional wires, snares and / or forceps) may result in additional trauma to the vessels.
- The safety of using multiple Ranger DCBs with a total drug dosage exceeding 9266 µg of paclitaxel in a patient has not been studied.
- Using a drug-eluting stent in conjunction with Ranger DCB at the same treatment site has not been studied.

#### Pre- and Post-Procedure Antiplatelet Therapy

The device carries an associated risk of thrombosis. It is strongly advised that the treating physician follow the Inter-Society Consensus (TASC II) Guidelines recommendations (or other applicable country guidelines) for antiplatelet therapy pre- and post-procedure to reduce the risk of thrombosis.

#### PRECAUTIONS

- Any use for procedures other than those indicated in these instructions is not recommended.
- Use the balloon catheter prior to the "Use By" date specified on the package.
- The balloon catheter should be used with caution for procedures involving calcified lesions due to the abrasive nature of these lesions.
- The balloon catheter is not intended for injection of contrast medium.
- Significant stenosis (> 50 %) proximal to the target lesion must be pretreated to prevent delamination of the balloon coating while crossing the lesion.
- Carefully inspect the balloon catheter prior to use to verify that it has not been damaged during shipment and that its size, shape and condition are suitable for the procedure for which it is to be used.
- If unusual resistance is felt during positioning of the loading tool, do not use the balloon catheter. Replace with a new Ranger DCB.
- The distal tip of the balloon catheter should be visible when backloading the distal tip onto the guidewire.
- Do not touch, wipe, bend, or squeeze the balloon. Do not allow it to contact any liquids including organic solvents such as alcohol or detergents prior to insertion. Damage to the balloon coating or premature release of the drug may occur.
- To minimize the possible introduction of air into the system, it is imperative that prior to proceeding, careful attention is paid to the maintenance of tight catheter connections and thorough aspiration and flushing of the system.
- If using a Tuohy-Borst type adapter, take care not to over-tighten the hemostatic valve around the catheter shaft as lumen constriction may occur, affecting inflation / deflation of the balloon or damaging the drug coating.
- Never advance the balloon catheter without the guidewire extending from the distal tip.
- If treating a long lesion (longer than the maximum balloon length available), each individual segment should be treated only once with a drug-coated balloon. Treat each segment with a new balloon, and try to minimize overlapping of treated segments.

#### Pregnancy / Lactation

This product has not been tested in pregnant or breastfeeding women or in men intending to father children; effects on the developing fetus have not been studied and the risks and reproductive effects remain unknown.

It is not recommended that the Ranger DCB be used in women attempting to conceive, or who are pregnant.

Prior to use, careful consideration should be given to the continuation of breastfeeding, taking into account the importance of the procedure to the mother. It is not known whether paclitaxel is distributed in human milk. In lactating rats, milk concentrations appeared to be higher than maternal plasma levels and declined in parallel with the maternal levels. Mothers should be advised of the potential for serious adverse reactions to paclitaxel in nursing infants.

#### Drug Information

The mechanism of action by which paclitaxel reduces or reverses neointima formation and proliferation, leading to restenosis, as demonstrated in clinical studies has not been established. It is known that paclitaxel promotes the assembly of microtubules from tubulin dimers and stabilizes microtubules by preventing depolymerization.

This stability results in the inhibition of the normal dynamic reorganization of the microtubule network that is essential for vital interphase and mitotic cellular functions.

#### Drug Interaction

Possible interactions of paclitaxel with concomitantly administered medications have not been formally investigated. Drug interactions of systemic chemotherapeutic levels of paclitaxel with possible concomitant medications are outlined in the labeling for finished pharmaceuticals containing paclitaxel, such as TAXOL.

#### Carcinogenicity, Genotoxicity, and Reproductive Toxicology

Published literature has reported paclitaxel to be aneugenic at concentrations similar to tissue concentrations recorded after deployment of a Ranger DCB in preclinical studies. The aneugenic effect is due to its pharmacodynamics action, interfering with microtubule disassembly, which also is the basis for the pharmacodynamic action preventing restenosis in vascular tissue following treatment with a Ranger DCB.

The relevance of both this observation and the aneugenic mechanism of genotoxicity for human carcinogenicity is currently not known.

#### ADVERSE EVENTS

Potential adverse events which may be associated with the use of Ranger Drug-Coated Balloon include, but are not limited to, the following:

- Allergic reaction (device, contrast medium, medications)
- Arteriovenous fistula
- Death
- Embolism (air, plaque, thrombus, device, tissue, or other)
- Hematoma
- Hemorrhage / Bleeding
- Hypotension / Hypertension
- Infection / Sepsis
- Restenosis of treated vessel
- Thrombus / Thrombosis
- Vessel occlusion
- Vessel trauma (e.g., dissection, perforation, pseudoaneurysm, rupture)
- Vasospasm

Potential adverse events not captured above that may be unique to the paclitaxel drug coating:

- Allergic / immunologic reaction to drug (paclitaxel or structurally-related compounds) or coating or its individual components
- Alopecia
- Anemia
- Gastrointestinal symptoms
- Hematologic dyscrasia (including leukopenia, neutropenia, thrombocytopenia)
- Hepatic enzyme changes
- Histologic changes in vessel wall, including inflammation, cellular damage or necrosis
- Myalgia / Arthralgia
- Peripheral neuropathy

There may be other potential adverse events that are unforeseen at this time. Apart from hypersensitivity reactions (allergic / immunologic reactions), the likelihood of paclitaxel related adverse events is low due to the low exposure.

#### CLINICAL STUDIES

##### Late Mortality Signal for Paclitaxel-Coated Devices

A meta-analysis of randomized controlled trials published in December 2018 by Katsanos et. al. identified an increased risk of late mortality at 2 years and beyond for paclitaxel-coated balloons and paclitaxel-eluting stents used to treat femoropopliteal arterial disease. In response to these data, FDA performed a patient level meta-analysis of long-term follow-up data from the pivotal premarket randomized trials of paclitaxel-coated devices used to treat femoropopliteal disease using available clinical data through May 2019. The meta-analysis also showed a late mortality signal in study subjects treated with paclitaxel-coated devices compared to patients treated with uncoated devices. Specifically, in the 3 randomized trials with a total of 1090 patients and available 5-year data, the crude mortality rate was 19.8 % (range 15.9 % - 23.4 %) in patients treated with paclitaxel-coated devices compared to 12.7 % (range 11.2 % - 14.0 %) in subjects treated with uncoated devices. The relative risk for increased mortality at 5 years was 1.57 (95 % confidence interval 1.16 - 2.13), which corresponds to a 57 % relative increase in mortality in patients treated with paclitaxel-coated devices. As presented at the June 2019 FDA Advisory Committee Meeting, an independent meta-analysis of similar patient-level data provided by VIVA Physicians, a vascular medicine organization, reported similar findings with a hazard ratio of 1.38 (95 % confidence interval 1.06 - 1.80). Additional analyses have been conducted and are underway that are specifically designed to assess the relationship of mortality to paclitaxel-coated devices.

The presence and magnitude of the late mortality risk should be interpreted with caution because of multiple limitations in the available data, including wide confidence intervals due to a small sample size, pooling of studies of different paclitaxel-coated devices that were not intended to be combined, substantial amounts of missing study data, no clear evidence of a paclitaxel dose effect on mortality, and no identified pathophysiologic mechanism for the late deaths.

Paclitaxel-coated balloons and stents improve blood flow to the legs and decrease the likelihood of repeat procedures to reopen blocked blood vessels compared to uncoated devices. The benefits of paclitaxel-coated devices (e.g., reduced reinterventions) should be considered in individual patients along with potential risks (e.g., late mortality).

The RANGER SFA first-in-man study enrolled 105 subjects. Given the limited sample size with the 2:1 randomization in the RANGER SFA study, no conclusion can be drawn about this potential mortality risk at 2 years and 3 years. Kaplan Meier mortality estimates at 2 years and 3 years are 9.0 % (95 % CI: 2.1%, 15.9%) and 13.6 % (95 % CI: 5.4 %, 21.8 %) for the Ranger treatment device and 10 % (95 % CI: 0.0 %, 20.8 %) and 10 % (95 % CI: 0.0 %, 20.8 %), respectively, for the standard PTA control device, which fall within the expected mortality rates for this patient population. All causes of death (ACD) at one year: 2.8 % (2/71) for Ranger treatment device and 2.9 % (1/34) for Standard PTA device with the difference between Ranger treatment and Standard PTA devices -0.1 % [95 % CI: -7.0 %, 6.7 %]. ACD at two years: 8.7 % (6/69) for Ranger treatment device and 9.7 % (3/31) for Standard PTA device with the difference between Ranger treatment and Standard PTA -1.0 % [95 % CI: -13.3 %, 11.4 %]. ACD at three years: 13.8 % (9/65) for Ranger treatment device and 10.7 % (3/28) for Standard PTA device with the difference between Ranger treatment device and Standard PTA 3.1 % [95 % CI: -11.1 %, 17.3 %].

The RANGER II SFA trial enrolled 376 randomized control subjects. Given the limited sample size with the 3:1 randomization in the RANGER II SFA study, no conclusion can be drawn about this potential mortality risk at one year. Kaplan Meier mortality estimates at one and two years are 1.9 % (95 % CI: 0.3 %, 3.5 %) and 5.9 % (95 % CI: 3.2 %, 8.6 %)\* for the Ranger treatment device and 2.1 % (95 % CI: 0.0 %, 5.0 %) and 3.1 % (95 % CI: 0.0 %, 6.4 %)\* for the standard PTA control device, which fall within the expected mortality rates for this patient population. All causes of death (ACD) at one year 2.2 % (6/274) for Ranger treatment device and 2.0 % (2/98) for Standard PTA control device with the difference between Ranger treatment and Standard PTA devices 0.1 % [95 % CI: -3.1 %, 3.4 %]\*. ACD at two years 5.9 % (16/269) for Ranger treatment device and 3.1 % (3/98) for Standard PTA control device with the difference between Ranger treatment and Standard PTA devices 2.9 % [95 % CI: -1.5 %, 7.3 %].\*

\* Survival status data were collected as part of the Ranger II SFA 24 month mortality dataset which served as the basis for KM survival analysis. Survival status data were not collected for the 12 Month Ranger II SFA study or the Ranger SFA FIH study.

## HOW SUPPLIED

The Ranger DCB has been sterilized with ethylene oxide.

## Device Details

One (1) Ranger DCB with loading tool is loaded into a coiled carrier tube which is placed inside a sterile barrier pouch and packed within a shelf carton.

Do not use if package is damaged or unintentionally opened before use.

Do not use if labeling is incomplete or illegible.

## Handling and Storage

This product has no special handling or storage requirements

## OPERATIONAL INSTRUCTIONS

### Additional Items for Safe Use

One or more of each of the following materials are required for PTA with the Ranger Balloon Catheter:

- Guidewire(s) of appropriate size for advancement of guide sheath
- Appropriate introducer sheath and dilator set
- Guide sheaths in the appropriate size and configuration to select the targeted artery
- Vial of contrast medium
- Vial of sterile saline solution
- PTA balloon catheter for pre-dilatation
- Guidewire(s) that are of appropriate length relative to the procedure requirements
- 10 ml, 12 ml, or 20 ml luer-lock inflation device with manometer
- Hemostatic valve
- Three-way stopcock

**Note:** Select guidewire, introducer, and guide sheaths as appropriate per the label.

## Inspection Prior to Use

Prior to angioplasty, carefully examine all equipment to be used during the procedure, including the balloon catheter, to verify proper function. Verify that the balloon catheter and sterile packaging have not been damaged. Verify that the catheter size is suitable for the specific procedure for which it is intended. Do not use if sterile package is damaged.

**Caution:** Do not touch, wipe, bend, or squeeze the balloon, or allow it to come into contact with any liquids prior to insertion. Damage to the balloon coating or premature release of the drug may occur.

**Note:** Do not use the balloon catheter if damage occurs or sterility is compromised.

## PREPARATION

### Inflation Device Preparation

1. Prepare the inflation device according to the manufacturer's instructions.
2. Purge the system of air.

### Balloon Catheter Selection

The inflation diameter of the balloon catheter must not exceed the diameter of the targeted artery proximal and distal to the stenosis. Similarly, the inflated length of the balloon (shoulder to shoulder) may exceed the length of the lesion / stenosis by approximately 10 mm on either side within the targeted artery. Careful consideration should also be given to select an appropriate catheter shaft length that considers the relationship between the arterial access point and the location of the targeted lesion or stenosis.

**Caution:** Observe the following precautions to prevent local overdosing:

- Do not use a second drug-coated balloon at the same treatment site.
- If treating a long lesion (longer than the maximum balloon length available), each individual segment should be treated only once with a drug-coated balloon. Treat each segment with a new balloon, and try to minimize overlapping of treated segments.
- Avoid implanting a drug-eluting stent at the same treatment site since an overdosage or interaction between the active agents may occur.

### Balloon Catheter Preparation

1. Remove the balloon catheter from the protective coil. Use care when removing the balloon catheter to avoid damage (e.g., shaft kink).
2. Prepare the balloon catheter for purging. Fill a 10 ml, 12 ml, or 20 ml luer-lock inflation device with appropriate balloon catheter inflation medium.

**Warning:** Use only the recommended balloon inflation medium (50:50 mixture of contrast medium and sterile saline solution). Never use air or any gaseous medium to inflate the balloon.

3. Connect a three-way stopcock to the inflation port fitting on the balloon catheter. Close the stopcock to the balloon catheter and flush through the stopcock. Use care when connecting the balloon catheter to avoid damage (e.g., shaft kink).
4. Connect the inflation device to the stopcock. Assure luer connections are properly aligned to avoid stripping the luer thread and causing subsequent leakage.
5. Open the stopcock to the balloon catheter. Hold the inflation device with the nozzle pointing downward and aspirate for 15 seconds to 20 seconds. Release the plunger or open the stopcock to air.
6. Remove the inflation device and evacuate all air from the barrel.
7. To prevent the possibility of air embolization, reconnect the inflation device and aspirate for 15 seconds to 20 seconds until bubbles no longer appear. If bubbles persist, check the luer connections. Discard and replace the balloon catheter if the problem cannot be resolved.
8. Disconnect the inflation device used in preparation. To remove any air lodged in the distal luer fitting of the inflation device, purge approximately 1 ml of contrast medium while pointing the inflation device upwards.
9. Verify that a meniscus of contrast medium is evident in both the inflation port and the inflation device connection to ensure a fluid-to-fluid connection. Adding a drop of inflation medium to the port may be necessary. Securely couple the inflation device to the inflation port.
10. Open the stopcock to the balloon catheter and maintain neutral pressure.
11. Move the loading tool proximally until only the distal tip of the balloon catheter is visible. If unusual resistance is encountered, first move the loading tool slightly distal (no more than 1 cm).
12. With the tip of the catheter pointed down, flush the guidewire lumen of the balloon catheter with sterile saline. Check for bends, kinks and other damage. Do not use if any defects are noted.

## PROCEDURE

### Insertion Procedure

**Note:** It is recommended to always use a smaller diameter balloon dilatation catheter to pre-dilate the stenosis to facilitate passage of the drug-coated balloon catheter.

**Caution:** Significant stenosis (> 50 %) proximal to the target lesion must be pretreated to prevent delamination of the balloon coating while crossing the lesion.

#### 1. Catheter Advancement

- Prepare the vascular access site according to standard practice.
- Insert a guidewire through the hemostatic valve following the manufacturer's instructions or standard practice. Advance the guidewire carefully into the guide or introducer sheath. When complete, withdraw the guidewire introducer, if used.
- Attach a torque device to the guidewire, if desired. Under fluoroscopy, advance the guidewire to the desired vessel, then across the stenosis.
- Thoroughly aspirate and flush the guide or introducer sheath in preparation for introduction of the balloon catheter.
- Backload the distal tip of the balloon catheter onto the guidewire ensuring that the guidewire exits the wire port of the catheter.

**Caution:** The distal tip of the balloon catheter should be visible when backloading the distal tip onto the guidewire.

**Caution:** When loading or exchanging the balloon catheter, it is recommended to thoroughly wipe the guidewire clean for better catheter movement on the guidewire.


**Caution:** Do not touch, wipe, bend, or squeeze the balloon, or allow it to come into contact with any liquids prior to insertion. Damage to the balloon coating or premature release of the drug may occur.

**Note:** To avoid kinking, advance the catheter slowly, in small increments.

- Holding the proximal end of the loading tool with the balloon catheter, carefully insert both of them into the hemostatic valve of the guide or introducer sheath (refer to Figure 1). Ensure that the loading tool and balloon catheter are advanced together.

**Caution:** If unusual resistance is felt during loading tool insertion into sheath valve, do not insert loading tool into sheath valve.

**Caution:** When using a sheath with an inner diameter  $\geq$  8F (2.67 mm), do not insert the loading tool more than 10 mm to 20 mm into the sheath.

**Figure 1: Inserting the Balloon Catheter with the Loading Tool**

- While holding the loading tool in place, slowly advance the catheter through the introducer sheath until the balloon is past the hemostasis valve. If unusual resistance is felt, do not advance the balloon catheter through the valve.

**Caution:** If using a Tuohy-Borst type adapter, take care not to over-tighten the hemostatic valve around the catheter shaft as lumen constriction may occur, affecting inflation / deflation of the balloon or damaging the drug coating.

- Remove the loading tool from the hemostatic valve of the guiding sheath and slide proximally along the balloon catheter shaft.

**Note:** If the loading tool meets the guidewire exchange port prior to the balloon fully advancing into the sheath, peel or partially peel the loading tool to access the additional length required to fully insert the balloon through the hemostasis valve.

- Connect the side port of the guiding sheath hemostatic valve to the proximal pressure recording / infusion line or manifold assembly, which permits proximal pressure recording or infusion through the guiding sheath.
- Advance the balloon catheter over the guidewire under direct fluoroscopic visualization, and position the balloon relative to the lesion to be treated. Use the radiopaque marker bands as a reference point. The outside edges of the marker bands indicate the balloon shoulders. Balloon inflation should not be undertaken if the balloon is not properly positioned within the stenosis.

**Note:** To access the entire working length of the Ranger DCB, the loading tool should be removed. Slide the loading tool to the guidewire exchange port, find the slit in the loading tool, and then simultaneously peel the two sides away.

#### 2. Balloon Inflation

- Slowly inflate the balloon to the appropriate pressure (refer to Table 2). Fluoroscopic visualization during balloon expansion should be used to properly judge the optimum expanded balloon diameter as compared to the proximal and distal vessel diameter(s). The recommended balloon inflation time is three (3) minutes.

**Warning:** Do not exceed the balloon rated burst pressure. Use of an inflation device is recommended to prevent over-pressurization.

**Warning:** If difficulty is experienced during balloon inflation, do not continue. Deflate the balloon and remove the catheter.

- Apply negative pressure to fully deflate the balloon. Confirm that the balloon is fully deflated under fluoroscopy.
- Withdraw the balloon catheter until it is clear of the lesion. Maintain the guidewire across the stenosis.

**Note:** Therapeutic drug levels will be administered only during the first inflation.

#### 3. Catheter Removal

- Confirm with angiography that the lumen of the dilated vessel has not abruptly occluded. Also ensure the balloon is fully deflated.
- While maintaining negative pressure, withdraw the deflated balloon catheter and guidewire from the guide or introducer sheath through the hemostatic valve.

**Warning:** Do not advance or retract the catheter unless the balloon is fully deflated under vacuum. If resistance is felt during manipulation, determine the cause of the resistance before proceeding. If the source of resistance cannot be determined, it is recommended to extract the entire system with the guide or introducer sheath.

**Table 2: Typical Ranger Balloon Catheter Compliance**

| Pressure |      | Mean Balloon Diameter (Balloon Lengths) |                |                 |               |                |                 |               |                |                 |
|----------|------|-----------------------------------------|----------------|-----------------|---------------|----------------|-----------------|---------------|----------------|-----------------|
|          |      | 4.0 mm                                  | 5.0 mm         | 5.0 mm          | 6.0 mm        | 6.0 mm         | 6.0 mm          | 7.0 mm        | 7.0 mm         | 7.0 mm          |
| atm      | kPa  | all lengths                             | 30 mm - 100 mm | 120 mm - 200 mm | 30 mm - 60 mm | 80 mm - 100 mm | 120 mm - 200 mm | 30 mm - 40 mm | 60 mm - 100 mm | 120 mm - 200 mm |
| 6        | 608  | 3.96                                    | 5.02           | 4.89            | 6.01          | 6.22           | 6.06            | 6.92          | 7.13           | 7.01            |
| 7        | 709  | 4.04                                    | 5.11           | 5.00            | 6.09          | 6.32           | 6.16            | 7.02          | 7.25           | 7.12            |
| 8        | 811  | 4.11                                    | 5.19           | 5.07            | 6.17          | 6.40           | 6.23            | 7.11          | 7.35           | 7.20            |
| 9        | 912  | 4.16                                    | 5.25           | 5.15            | 6.24          | 6.47           | 6.30            | 7.18          | 7.43           | 7.27            |
| 10       | 1013 | 4.22                                    | 5.31           | 5.21            | 6.29          | 6.53           | 6.36            | 7.25          | 7.49           | 7.34            |
| 11       | 1115 | 4.26                                    | 5.36           | 5.27            | 6.35          | 6.57           | 6.42            | 7.31          | 7.56           | 7.40            |
| 12       | 1216 | 4.30                                    | 5.41           | 5.32            | 6.39          | 6.62           | 6.46            | 7.36          | 7.61           | 7.46            |
| 13       | 1317 | 4.33                                    | 5.45           | 5.36            | 6.43          | 6.67           | 6.51            | 7.41          | 7.66           | 7.52            |
| 14       | 1419 | 4.37                                    | 5.49           | 5.41            | 6.48          | 6.72           | 6.56            | 7.46          | 7.73           | 7.58            |

#### Disposal

To minimize the risk of infection or microbial hazards after use, dispose of device and packaging as follows:

After use, device may contain biohazardous substances. The device and packaging should be treated and disposed of as biohazardous waste or have them treated and disposed of in accordance with any applicable hospital, administrative, and / or local government regulations. Use of a biohazardous container with biological hazard symbol is recommended. Untreated biohazardous waste should not be disposed of in the municipal waste system. Biohazardous substances should be removed via incineration prior to disposal. Alternatively, biohazardous waste may be disposed of utilizing a certified facility for biohazardous waste for proper treatment in accordance with hospital, administrative, and / or local government policy.

#### Post-Procedure

Assess patient for hematoma and / or other signs of bleeding at the puncture site. Determine appropriate antiplatelet therapy based on the Inter-Society Consensus (TASC II) Guidelines recommendations (or other applicable country guidelines) for antiplatelet therapy post-procedure to reduce the risk of thrombosis.

Long-term follow-up of ongoing studies to further support the safety and effectiveness of the Ranger Monorail Paclitaxel-Coated PTA Balloon Catheter are being conducted.

Any serious incident that occurs in relation to this device should be reported to the manufacturer and relevant local regulatory authority.

#### PATIENT COUNSELING INFORMATION

The physician should consider the following points while counseling patients on the use of Ranger DCB in association with the interventional procedure:

- Discuss the risks and benefits including review of potential adverse events listed in this document, both for Ranger DCB and for other interventional treatments likely to be employed.
- Discuss patient allergies in particular the risk for patients who may be allergic to paclitaxel or contrast.
- Discuss the risks and benefits of antiplatelet therapy including risk of thromboembolism should the patient be allergic or discontinue use.
- Discuss post procedure instructions, including any follow-up appointments, lifestyle changes, medications, and home-care or rehabilitation guidelines.

#### WARRANTY

For device warranty information, visit ([www.bostonscientific.com/warranty](http://www.bostonscientific.com/warranty)).

#### SYMBOL DEFINITIONS

Commonly used medical device symbols that appear on the labeling are defined at [www.bostonscientific.com/SymbolsGlossary](http://www.bostonscientific.com/SymbolsGlossary)

The following are trademarks of Boston Scientific Corporation or its affiliates: MONORAIL, Ranger.

TAXOL is a trademark of Bristol-Myers Squibb Company.

All other trademarks are the property of their respective owners.

**Table 3. Symbol Definitions**

|                                                                                     |                               |
|-------------------------------------------------------------------------------------|-------------------------------|
|  | Contents                      |
|  | Recommended Guidewire         |
|  | Recommended Introducer Sheath |
|  | Maximum Crossing Profile      |

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All rights reserved.
Boston Scientific Limited
Boston Scientific (Australia) Pty Ltd
©2025 Boston Scientific Corporation or its affiliates.