VersaCross Connect™ Access Solution for FARADrive™

Access Solution

en	Instructions for Use	2
es	Instrucciones de uso	10
fr	Mode d'emploi	18
de	Gebrauchsanweisung	26
it	Istruzioni per l'uso	34
nl	Instructies voor gebruik	42
pt-EU	Instruções de Utilização	50

TABLE OF CONTENTS

REUSE WARNING 3

DEVICE DESCRIPTION 3

 Contents 3

 Table 1. Product Specifications 3

 Materials 4

 Non-pyrogenic 4

 User Information 4

INTENDED USE 4

INDICATIONS FOR USE 4

 Clinical Benefit Statement 4

 Summary of Safety and Clinical Performance 4

CONTRAINDICATIONS 4

WARNINGS 4

PRECAUTIONS 5

ADVERSE EVENTS 6

HOW SUPPLIED 6

 Device Details 6

 Handling and Storage 6

OPERATIONAL INSTRUCTIONS 6

 Prior to Use 6

 Additional Required Items 6

 Installation 6

 Figure 1. Connections 7

 Preparation 7

 Procedure 7

 Disposal 8

 Post-Procedure 8

PATIENT COUNSELING INFORMATION 8

TROUBLESHOOTING 9

 Table 2. Troubleshooting 9

WARRANTY 9

SYMBOL DEFINITIONS 9

 Table 3. Symbol Definitions 59

VersaCross Connect™ Access Solution for FARADrive™

Access Solution

Notice (EU): Any serious incident as defined in the EU Medical Device Regulation (2017/745) that has occurred in relation to the VersaCross Connect Access Solution for FARADrive should be reported to Boston Scientific Corporation and the competent authority of the EU Member State in which the physician and/or patient is established.

REUSE WARNING

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness, or death.

Reuse, reprocessing, or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

Carefully read all ancillary device instructions prior to use. Observe all contraindications, warnings and precautions noted in these directions. Failure to do so may result in patient complications.

DEVICE DESCRIPTION

The VersaCross Connect Access Solution for FARADrive consists of 3 components: a 13 F dilator (4.4 mm maximum OD) VersaCross Connect Transseptal Dilator, a VersaCross RF Wire and a single-use Connector Cable (Connector Cable). The VersaCross RF Wire must be used with an approved RFP-100A Radiofrequency (RF) Puncture Generator (Generator) and the Connector Cable.

The VersaCross Connect Transseptal Dilator is designed for safe and easy catheterization and angiography of specific heart chambers and locations. The dilator provides torque control and is flexible. The dilator features a tapered tip and a shaft that can be reshaped manually. The echogenic shaft and tip and radiopaque tip maximize visualization of the dilator during manipulation.

The VersaCross Connect Transseptal Dilator is for use with a 13 F ID FARADrive Steerable Sheath which is 74 cm in length, specifically, models: 21M401, 21M402, M004PFCE21M402.

Carefully read the applicable FARADrive Steerable Sheath instructions prior to use. Observe all warnings and precautions noted in these instructions. Failure to do so may result in patient complications.

The VersaCross RF Wire delivers radiofrequency (RF) power in a monopolar mode between its distal electrode and a commercially available external Disposable Indifferent (Dispersive) Patch (DIP) Electrode, which is in compliance with current IEC 60601-2-2 requirements. The Connector Cable connects the Generator to the VersaCross RF Wire. This Connector Cable enables RF power to be delivered from the Generator


to the VersaCross RF Wire. Detailed information concerning the RFP-100A RF Puncture Generator is contained in a separate manual that accompanies the Generator.

The dimensions of the VersaCross RF Wire and the Connector Cable can be found on the product labels or in Table 1. The insulation on the body of the VersaCross RF Wire facilitates smooth advancement of the device and provides electrical insulation. The floppy distal portion of the VersaCross RF Wire has a curve and the active tip is rounded to be atraumatic to cardiac tissue unless RF energy is applied. A radiopaque and echogenic marker coil is positioned on the distal section for visualization during manipulation. The main body of the VersaCross RF Wire provides a stiff rail for advancing ancillary devices into the left atrium following the creation of an atrial septal defect. The VersaCross RF Wire features visible markers along its length to assist with aligning the wire tip in a compatible transseptal sheath and/or dilator assembly (e.g., the VersaCross Transseptal Sheath kit). The proximal end of the VersaCross RF Wire is bare metal to connect only with the provided Connector Cable and not with electrocautery or electrosurgery devices. The other end of the Connector Cable connects to the Generator.

Contents

- One (1) VersaCross Connect Transseptal Dilator
- One (1) VersaCross RF Wire
- One (1) Connector Cable (Single-Use)

Table 1. Product Specifications

Device	Specification
Transseptal Dilator	Outer Diameter: 13 F (4.4 mm) Inner Diameter: 0.036 in (0.9 mm) Usable Length: 93 cm
RF Wire	Length: 180 cm or 230 cm Outer Diameter: 0.035 in (0.89 mm) Distal Type: - J-Tip (Radius: 4.5 mm, Diameter: 9 mm) - Pigtail (Diameter 24 mm) Rated Voltage: 300 V _{rms}
Connector Cable	Length: 10 ft (3 m) Device Connector: Push Button Generator Connector: 4-Pin (3-Pin) Rated Voltage: 300 V _{rms}  IP Rating: IP50 (Dust-protected, not protected against water)

Materials

Contains cobalt; CAS No. 7440-48-4; EC No. 231-158-0. Defined as a 1B carcinogen and mutagen according to the European Commission in a concentration about 0.1% weight by weight.

Note: This device is made with a metal alloy which contains cobalt. Current scientific evidence supports that metal alloys containing cobalt used in medical devices do not cause an increased risk of cancer or adverse reproductive effects.

Non-pyrogenic

This device meets pyrogen limit specifications for all patient-contacting parts. The RF wire and dilator supplied are non-pyrogenic unless package is opened or damaged.

User Information

Only physicians with a thorough understanding of angiography, aseptic technique and percutaneous interventional procedures should use this device. It is recommended that physicians avail themselves of pre-clinical training, a review of pertinent literature and other appropriate education before attempting new interventional procedures.

INTENDED USE

The VersaCross Connect Access Solution for FARADrive is indicated for creation of an atrial septal defect in the heart and for the percutaneous introduction of various types of cardiovascular catheters and guidewires to all heart chambers, including the left atrium via transseptal perforation / puncture. The VersaCross Connect Access Solution for FARADrive is intended for use in percutaneous transseptal procedures in patients requiring delivery of cardiovascular catheters to the heart, including the left side via transseptal puncture. The VersaCross Connect Access Solution for FARADrive is intended for use in adult (not pediatric) patients, with the exclusion of pregnant and/or nursing patients. VersaCross Connect Access Solution for FARADrive is a single-use device. Only physicians trained in the techniques of the approach to be used should perform the interventional procedures. The RF Wire included is intended to deliver radiofrequency power in a monopolar mode between its distal electrode and a return electrode. The RF Wire can be used without exchanges, as a guidewire, as a transseptal device or as an exchange rail for delivering therapy sheaths. The RF Wire included is not recommended for use with any conditions that do not require the creation of an atrial septal defect.

INDICATIONS FOR USE

The VersaCross Connect Access Solution for FARADrive is indicated for creation of an atrial septal defect in the heart and for the percutaneous introduction of various types of cardiovascular catheters and guidewires to all heart chambers, including the left atrium via transseptal perforation / puncture.

Clinical Benefit Statement

The clinical benefit of using the VersaCross Connect Access Solution for FARADrive is the successful creation of an atrial septal defect in the heart and the percutaneous introduction of various types of cardiovascular catheters and guidewires to all heart chambers, including the left atrium via transseptal perforation / puncture.

Summary of Safety and Clinical Performance

For customers in the European Union, this information is available on the European database on medical devices (Eudamed) website:

<https://ec.europa.eu/tools/eudamed>

Use the manufacturer and device name found in the labeling to search for the device's

Summary of Safety and Clinical Performance.

CONTRAINDICATIONS

The included VersaCross RF Wire is not recommended for use with any conditions that do not require the creation of an atrial septal defect. The Connector Cable is not recommended for use with any other RF Generator or any other device.

The VersaCross RF Wire is not intended for use with pediatric patients.

WARNINGS

- Only physicians with a thorough understanding of angiography, aseptic technique and percutaneous interventional procedures should use this device. It is recommended that physicians avail themselves of pre-clinical training, a review of pertinent literature and other appropriate education before attempting new interventional procedures.
- Laboratory staff and patients can undergo significant x-ray exposure during interventional procedures due to the continuous usage of fluoroscopic imaging. This exposure can result in acute radiation injury as well as increased risk for somatic and genetic effects. Therefore, adequate measures must be taken to minimize this exposure. The use of echocardiography is recommended.
- The VersaCross Connect Access Solution for FARADrive is intended for single patient use only. Do not attempt to sterilize and reuse any component of the VersaCross Connect Access Solution for FARADrive. Reuse can cause the patient injury and/or the transmission of infectious disease(s) from one patient to another. Failure to follow this instruction may result in patient complications.
- The VersaCross Connect Access Solution for FARADrive is supplied sterile using an ethylene oxide process. Do not use if the package is damaged. If damage is found, call your Boston Scientific representative.
- Manual shaping of the dilator distal curve shall be done with smooth motions along the curve. Do not use excessive force and/or pressure when reshaping.
- The VersaCross RF Wire must be used with the Connector Cable provided. Attempts to use it with other connector cables can result in electrocution of the patient and/or operator.
- Do not use the VersaCross RF Wire with electrocautery or electrosurgery generators, connector cables or accessories as attempted use can result in patient and/or operator injury.
- The Connector Cable must only be used with the RFP-100A RF Puncture Generator and the included VersaCross RF Wire. Attempts to use it with other RF Generators and devices can result in electrocution of the patient and/or operator.
- The VersaCross RF Wire must be used with 0.035 in (0.89 mm) compatible transseptal sheath and/or dilator devices. Use of incompatible accessory devices may damage the integrity of the VersaCross RF Wire or accessory devices and may cause patient injury.
- The VersaCross RF Wire has only been validated for transseptal puncture use through VersaCross dilators which have been demonstrated to provide the required support for optimal function.
- The active tip and distal curve of the VersaCross RF Wire are fragile. Be careful not to damage the tip or the distal curve while handling the VersaCross RF Wire. If the tip or the distal curve becomes damaged at any time during its use, discard the VersaCross RF Wire immediately. Do not attempt to straighten the active tip if bent. Damage to device can lead to patient injury.
- The VersaCross RF Wire is not intended for use with pediatric patients. Do not

attempt to treat pediatric patients with the VersaCross RF Wire.

- Do not attempt to insert or retract the VersaCross RF Wire through a metal cannula or a percutaneous needle, which may damage the device and may cause patient injury.
- Care should be taken to ensure that all air is removed from the dilator before infusing through the proximal hub.
- Care should be taken when inserting or removing the dilator from access sheaths or introducer sheaths.
- Care should be taken when inserting or removing compatible guidewires from the dilator lumen.
- Direct percutaneous insertion of the dilator without a sheath is not recommended.
- Do not attempt direct percutaneous insertion of the dilator without a guidewire as this may cause vessel injury.
- Careful manipulation of the dilator must be performed to avoid cardiac damage or tamponade. Dilator advancement should be performed under imaging guidance. Fluoroscopic or echocardiographic imaging is recommended. If resistance is encountered, do not use excessive force to advance or withdraw the device.

PRECAUTIONS

- Do not attempt to use the VersaCross Connect Access Solution for FARADrive before thoroughly reading the accompanying Instructions for Use.
- The sterile barrier system, and all devices should be visually inspected prior to use. Do not use if the sterile barrier integrity or devices have been compromised or damaged.
- Careful manipulation must be performed to avoid cardiac damage, tamponade or vessel trauma. Dilator and wire advancement should be performed under imaging guidance, such as fluoroscopy or echocardiography. If resistance is encountered, do not use excessive force to advance or withdraw the VersaCross RF Wire or VersaCross Connect Transseptal Dilator. Excessive force may lead to bending or kinking of the VersaCross RF Wire limiting advancement and retraction of sheath and/or dilator device.
- Only physicians thoroughly trained in the techniques of the approach to be used should perform interventional procedures.
- RF puncture procedures should be performed only by physicians thoroughly trained in the techniques of RF-powered puncture in a fully equipped catheterization laboratory.
- Do not use the VersaCross Connect Access Solution for FARADrive after its "Use By" date.
- The VersaCross Connect Transseptal Dilator is compatible with introducer sheaths 13 F or larger.
- The VersaCross Connect Transseptal Dilator is for use with specified models of 13 F ID FARADrive Steerable Sheath which are 74 cm in length.
- The VersaCross Connect Transseptal Dilator is compatible with 0.035 in (0.89 mm) transseptal devices and guidewires or smaller.
- The VersaCross Connect Access Solution for FARADrive is not compatible with transseptal needles such as the "NRG Transseptal Needle".
- Visually inspect the VersaCross RF Wire and Connector Cable prior to use to ensure there are no cracks or damage to the insulating material. Do not use the wire or the cable if there is any damage.

- The VersaCross RF Wire and Connector Cable are intended for use with only those devices listed in Additional Required Items.
- Read and follow the manufacturer's Instructions For Use for the DIP electrode. Always use DIP electrodes that meet or exceed IEC 60601-2-2 requirements.
- Placement of the DIP electrode on the thigh could be associated with higher impedance.
- In order to prevent the risk of ignition, ensure that flammable materials are not present in the room during RF power application.
- Take precautions to limit the effects that the electromagnetic interference (EMI) produced by the Generator may have on the performance of other equipment. Check the compatibility and safety of combinations of other physiological monitoring and electrical apparatus to be used on the patient in addition to the Generator.
- Adequate filtering must be used to allow continuous monitoring of the surface electrocardiogram (ECG) during RF power applications.
- Do not attempt to insert and use the proximal end of the VersaCross RF Wire as the active tip.
- Do not bend the VersaCross RF Wire or the Connector Cable. Excessive bending or kinking of the wire shaft, distal curve of the wire and/or the Connector Cable may damage the integrity of the device components and may cause patient injury. Care must be taken when handling the VersaCross RF Wire and Connector Cable.
- VersaCross RF Wire and ancillary sheath and/or dilator assembly advancement should be done under imaging guidance. The use of visible markers on the wire body are only an approximate guide for positioning the wire tip with the distal end of the dilator.
- Do not attempt to deliver RF energy until the active tip of the VersaCross RF Wire is confirmed to be in good contact with the target tissue.
- Avoid RF energy delivery of the VersaCross RF Wire with incompatible dilator or cannula devices, which may lead to patient burns, ineffective puncture or failure to puncture.
- It is recommended not to exceed five (5) RF power applications per VersaCross RF Wire.
- Never disconnect the Connector Cable from the Generator while RF power is being delivered.
- Never disconnect the Connector Cable from the Generator by pulling on the cable. Failure to disconnect the cable properly may result in damage to the cable.
- Do not twist the Connector Cable while inserting or removing it from the Isolated Patient Connector on the Generator. Twisting the cable may result in damage to the pin connectors.
- The Generator is capable of delivering significant electrical power. Patient or operator injury can result from improper handling of the VersaCross RF Wire and/or DIP electrode, particularly when operating the device.
- During power delivery, the patient should not be allowed to come in contact with ground metal surfaces.
- Apparent low power output or failure of the equipment to function properly at normal settings may indicate faulty application of the DIP electrode, failure to an electrical lead, or poor tissue contact at the active tip. Check for obvious equipment defects or misapplication. Attempt to better position the active tip of the VersaCross RF Wire against the atrial septum. Only increase the power if low power output persists.

- If using electroanatomical mapping guidance, it is recommended to use it along with alternative imaging modality in the event there is loss of visibility of the device.

ADVERSE EVENTS

Adverse events that may occur while using the VersaCross Connect Access Solution for FARADrive include, but are not limited to, the following:

- | | |
|------------------------------|-------------------------------------|
| • Allergic Reaction | • Hypotension |
| • Arrhythmias | • Infection |
| • Burn | • Myocardial Infarction |
| • Cardiac Arrest | • Nerve Injury |
| • Cardiac Tamponade | • Pain |
| • Cardiac Trauma | • Pericardial Effusion |
| • Cerebral Vascular Accident | • Pneumothorax |
| • Death | • Radiation Exposure |
| • Electric Shock | • Respiratory Insufficiency/Failure |
| • Embolism | • Thrombosis |
| • Fistula | • Transient Ischemic Attack |
| • Heart Failure | • Vasovagal Response |
| • Hematoma | • Vessel Trauma |
| • Hemorrhage | |

There may be other potential adverse events that are unforeseen at this time.

HOW SUPPLIED

Device Details

The VersaCross Connect Access Solution for FARADrive is supplied sterile using an Ethylene Oxide (EO) process, and for single use only.

Do not use if package is damaged or unintentionally opened before use.

Do not use if labeling is incomplete or illegible.

Do not use the device if past the "Use By" date.

Handling and Storage

Keep away from sunlight.

OPERATIONAL INSTRUCTIONS

Prior to Use

Prior to use of the VersaCross Connect Access Solution for FARADrive, the individual components should be carefully examined for damage or defects, as should all equipment used in the procedure. Do not use defective equipment. Do not reuse the device.

Additional Required Items

RF transseptal procedures should be performed in a specialized clinical setting equipped with appropriate imaging equipment and compatible examination table, echocardiography imaging, physiologic recorder, emergency equipment and instrumentation for gaining vascular access. Ancillary materials required to perform this procedure include:

- RFP-100A RF Puncture Generator
- FARADrive Steerable Sheath


- DIP electrode, meeting or exceeding IEC 60601-2-2 requirements for electrosurgical electrodes (not included)
- 0.035 in (0.89 mm) maximum guidewire outer diameter or smaller (optional)
- DuoMode Cable for use with electroanatomic mapping systems (optional)

The following devices are compatible with the VersaCross Connect Access Solution for FARADrive:

- Male luer lock syringes
- Introducer sheaths 13 F or larger
- 13 F ID FARADrive Steerable Sheath, which is 74 cm in length
- 0.035 in (0.89 mm) or smaller Transseptal Puncture Devices and Guidewires (such as VersaCross RF wire)
- DIP electrode, meeting or exceeding IEC 60601-2-2 requirements for electrosurgical electrodes
- DuoMode Cable for use with electroanatomic mapping systems
- RFP-100A RF Puncture Generator

Installation

- The VersaCross Connect Access Solution for FARADrive is supplied sterile using an ethylene oxide process. Use aseptic technique when opening the packaging and handling the products in the sterile field.
- Connect the generator connector end of the Connector Cable to the isolated patient connector port on the Generator as per the RFP-100A RF Puncture Generator Instructions for Use. Gently line up the connector pins with the socket and push in until the connector fits firmly into the socket. Any attempt to connect the cable otherwise will damage the pins on the connector.
- Do not use excessive force in connecting the Connector Cable to the Generator. Use of excessive force may result in damage to the connector pins.

[1] Generator [2] Single-Use Connector Cable [3] DIP Grounding Pad (not included) [4] VersaCross RF Wire (Pigtail model shown) [5] Footswitch (optional)

Figure 1. Connections

Preparation

- Carefully read all instructions prior to use. Failure to do so may result in complications.
- Carefully read the applicable FARADRIE Steerable Sheath instructions prior to use. Failure to do so may result in complications.
- Carefully read the applicable compatible guidewire instructions prior to use. Observe all warnings and precautions noted in these instructions. Failure to do so may result in patient complications.
- Use of VersaCross Connect is expected to reduce the number of exchanges in the

procedure resulting in a more efficient transseptal puncture. This should be taken into account when estimating the timing for heparin administration to ensure appropriate ACT levels after transseptal puncture.

- Thoroughly flush the dilator with heparinized saline solution prior to use.

Procedure

- Gain access to the right femoral vein using standard methods.
- Note that a compatible access introducer sheath may be used at the venous cutaneous puncture site if desired. Refer to the compatible introducer sheath's Instructions for Use for details and directions.
- Introduce a compatible guidewire through the vasculature access point and advance to required depth (the VersaCross RF Wire may be used).
- Enlarge the cutaneous puncture site as necessary.
- Dilator can be inserted fully into the compatible access sheath.
- The distal curvature of the dilator or sheath and dilator assembly may be adjusted manually if desired prior to insertion into the body. Manual shaping of the distal curve shall be done with smooth motions along the curve. Do not use excessive force and/or pressure when reshaping.
- Thread the dilator and compatible access sheath over the VersaCross RF Wire (or guidewire if used) using a slight twisting motion into the superior vena cava (SVC) under imaging guidance such as fluoroscopy or echocardiography. If resistance is encountered, do not use excessive force to advance or withdraw the dilator over the guidewire. Determine the cause of resistance before proceeding.
- Use standard technique to position the sheath, dilator and guidewire assembly into the desired heart chamber.
- If the VersaCross RF Wire was not used to advance the sheath and dilator assembly to the SVC, remove the guidewire and exchange for the VersaCross RF Wire with the provided tip straightener.
- Advance the VersaCross RF Wire through the sheath and dilator assembly until the wire tip is just within the dilator tip. The visible markers on the wire body can be used to assist with the positioning of the wire tip with the distal end of the dilator.
- Firmly grasp the catheter connector end of the Connector Cable in one hand. Using your thumb, depress the red button on the top of the connector. Slowly insert the proximal end of the VersaCross RF Wire into the opening of the catheter connector. Once the exposed portion of the proximal end of the device is no longer visible, release the red button on the connector. Gently tug on the device to ensure that you have a secure connection.
- Position the tip of the transseptal assembly (RF wire, sheath, dilator) in the right atrium against the fossa ovalis under appropriate imaging guidance including but not limited to fluoroscopic, echocardiographic and/or electroanatomic mapping guidance using standard technique.

NOTE: If using electroanatomical mapping guidance, it is recommended to confirm tip placement and septal tenting with echocardiographic imaging or another imaging modality.

- Apply pressure to the dilator to tent the septum at the fossa ovalis.
- Advance the VersaCross RF Wire so that the active tip is engaging the septum at the fossa ovalis but still within the dilator.
- Once appropriate positioning has been achieved, deliver RF power via the Generator to the active tip. This results in puncture of the targeted cardiac tissue.

Please refer to the Generator Instructions for Use for the correct operation of the generator.

- Apply firm pressure to the VersaCross RF Wire during the application of RF energy to successfully advance the VersaCross RF Wire through the tissue.

NOTE: Use the lowest appropriate RF settings to achieve the desired puncture.

For RFP-100A: An initial RF setting between one (1) second on “PULSE” mode to two (2) seconds on “CONSTANT” mode has been shown to be sufficient for successful puncture.

- RF power delivery can be terminated by pressing the RF ON/OFF button on the Generator if the timer has not expired.
- Entry into the left atrium can be confirmed by monitoring the VersaCross RF Wire under appropriate imaging guidance. Echocardiographic guidance is also recommended.
- If septal puncture is not successful after five (5) RF power applications, it is advised that the user utilize an alternate method for the procedure.
- Once the puncture is successfully completed, the VersaCross RF Wire should be mechanically advanced without any RF power. Positioning in the left atrium is sufficient when the full distal curve and floppy section have crossed the septum and are observed in the left atrium. Echocardiographic guidance is also recommended.
- The dilator can then be advanced over the VersaCross RF Wire to enlarge the puncture.
- To disconnect the VersaCross RF Wire from the Connector Cable, depress the red button on the catheter connector and gently remove the proximal end of the VersaCross RF Wire from the Connector Cable.
- To disconnect the Connector Cable from the Generator, grasp the connector firmly and gently pull it straight out of the socket.
- To use the VersaCross RF Wire to facilitate advancement of ancillary devices, slowly retract the transseptal sheath and/or dilator assembly over the wire, holding the wire in place to ensure left atrium access is maintained.
- Advance the required ancillary device over the VersaCross RF Wire until the distal tip reaches the desired positioning in the left atrium.
- Retract the VersaCross RF Wire slowly through the transseptal sheath and dilator assembly.
- Ensure the dilator is clear of air. To aspirate blood, use the dilator hub.
- Monitor the location of the radiopaque tip frequently under imaging guidance, such as fluoroscopy or echocardiography.
- Deliver a continuous heparinized solution infusion or aspirate periodically. This may help reduce the risk of thromboembolic complications due to thrombus formation, as there may be a possibility of thrombus development at the distal dilator tip or inside the dilator lumen. Also aspirate when removing the transseptal device or dilator.
- After removal of the dilator, use standard technique to achieve hemostasis.

Disposal

To minimize risk of infection or microbial hazards after use, dispose of device and packaging as follows:

After use, the device(s) may contain biohazardous substances. The device(s) and packaging should be treated and disposed of as biohazardous waste in accordance with any applicable hospital, administrative, and/or local government regulations. Use of a biohazardous container with biological hazard symbol is recommended. Untreated biohazardous waste should not be disposed of in the municipal waste system.

Post-Procedure

Any serious incident that occurs in relation to this device should be reported to Boston Scientific and the relevant local regulatory authority.

PATIENT COUNSELING INFORMATION

The physician should consider the following points while counseling patients on the use of the VersaCross Connect Access Solution for FARADRIVE in association with the electrophysiological cardiac interventional procedure:

- Discuss the risks and benefits including review of potential adverse events listed in this document.
- Discuss post procedure instructions, including any lifestyle changes, medications, when to call the Healthcare Provider (HCP) and any post procedure follow-up that might be needed.
- Any serious incident that occurs in relation to this device should be reported to Boston Scientific and the relevant local regulatory authority.

TROUBLESHOOTING

The following table is provided to assist the user in diagnosing potential problems.

Table 2. Troubleshooting

PROBLEM	COMMENTS	TROUBLESHOOTING
Connector Cable does not fit into the Isolated Patient Connector on the front panel of the generator	The connectors are designed to connect in a specific way for safety reasons. If the connector "keys" are out of line, the connectors won't fit together.	Check that the connector keys are lined up in the proper orientation.
Generator Error Messages	In order to successfully perforate tissue using RF energy, all devices must be properly connected and in good working order.	Ensure that all connections are made i.e.: • VersaCross RF Wire to Connector Cable • Connector Cable to Generator • Generator to power outlet • Generator to grounding pad (not included) Visually inspect the VersaCross RF Wire and Connector Cable for damage. Immediately discard any damaged devices. If problem persists, discontinue use. For error messages encountered while attempting RF puncture, refer to the Instructions for Use that accompanies the Generator.
Wire breaks or kinks	Breaks and kinks in the VersaCross RF Wire are a potential cause of patient injury.	Discard immediately.

WARRANTY

For device warranty information, visit www.bostonscientific.com/warranty.

EU Importer: Boston Scientific International B.V., Vestastraat 6, 6468 EX, Kerkrade, The Netherlands
The following are trademarks of Boston Scientific Corporation or its affiliates: VersaCross Connect, VersaCross, FARADrive, RFP-100A, NRG, and DuoMode. All other trademarks are the property of their respective owners.

SYMBOL DEFINITIONS

Commonly used medical device symbols that appear on the labeling are defined at www.bostonscientific.com/SymbolsGlossary.

Additional symbols are defined at the end of this document.

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All rights reserved.
Boston Scientific (Australia) Pty Ltd
© 2025 Boston Scientific Corporation or its affiliates.