Speedband Superview Super 7™

Multiple Band Ligator

en	Instructions for Use	2
es	Instrucciones de uso	11
fr	Mode d'emploi	21
de	Gebrauchsanweisung	31
it	Istruzioni per l'uso	41
nl	Instructies voor gebruik	51
pt-EU	Instruções de Utilização	61

52073322-01

2025-09

TABLE OF CONTENTS

REUSE WARNING	3
DEVICE DESCRIPTION	3
Figure 1	3
Figure 2	3
Contents	4
Operating Principle	4
User Information	4
INTENDED USE/INDICATIONS FOR USE	4
Clinical Benefit Statement	4
CONTRAINDICATIONS	4
WARNINGS	4
PRECAUTIONS	4
ADVERSE EVENTS	4
HOW SUPPLIED	5
Device Details	5
Handling and Storage	5
OPERATIONAL INSTRUCTIONS	5
Preparation	5
Initial setup	5
Figure 3A. (Olympus™)	5
Figure 3B. (Pentax™)	5
Figure 3C. (Fujinon™)	5
Figure 3D	6
Figure 4A	6
Figure 4B	6
Figure 4C	7
Figure 5	7
Figure 6A	8

Figure 6B	8
Figure 7A	8
Figure 7B	8
Method for ligating an esophageal varix	9
Figure 8A	9
Figure 8B	9
Additional Clinical Information for Esophageal Variceal Ligation	9
Method for Ligating an Anorectal Hemorrhoid	9
Figure 9A	10
Figure 9B	10
Removing Speedband Superview Super 7	
Multiple Band Ligator from Endoscope	10
Disposal	10
Post-Procedure	10
PATIENT COUNSELING INFORMATION	10
WARRANTY	10
SYMBOL DEFINITIONS	10
Table 1. Symbol Definitions	10

Speedband Superview Super 7™

Multiple Band Ligator

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

REUSE WARNING


For single use only. Do not reuse, reprocess or sterilize. Reuse, reprocessing or sterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or sterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

DEVICE DESCRIPTION

The Speedband Superview Super 7 Multiple Band Ligator is composed of two main components.


1. Ligating Unit

The Ligating Unit is a cylinder which fits at the distal end of the endoscope. Elastic bands are stretched around the cylinder (Figure 1). The pulling loop on the Ligating Unit is affixed to the trip wire loop of the Handle Unit during initial setup.

[1] Elastic Bands [2] Pulling Loop

Figure 1.

[1] Scope Fastener [2] Trip Wire [3] Trip Wire Loop

Figure 2.

2. Handle Unit with Trip Wire and Scope Fastener

The Handle Unit contains a plastic spool, which turns only in the clockwise direction (Figure 2). The trip wire is attached to the plastic spool. When the plastic spool is rotated a half turn (180°), the handle will make a distinct "click" sound, and one band will be fired automatically. The Handle Unit also incorporates a scope fastener to secure the handle onto the endoscope. An irrigation

valve is located on the side of the Handle Unit. This valve is closed unless a syringe, or the irrigation tube provided, is fitted to its luer-tapered fitting.

The Speedband Superview Super 7 Multiple Band Ligator requires an endoscope with a scope O.D. of 8.6 mm – 11.5 mm.

Contents

- (1) Super 7 Ligating unit.
- (1) Handle Unit with trip wire and scope fastener.
- (1) Extension tube.

Operating Principle

Endoscopic band ligation treats esophageal varices and anorectal hemorrhoids. This device is capable of releasing 7 bands; a single band is released with each 180 ° clockwise rotation of the handle, and the handle makes a "click" sound with each 180 ° clockwise rotation. The bands are designed for strength and retention on the varix throughout the procedure.

User Information

The Speedband Superview Super 7 Multiple Band Ligator should only be used by or under the supervision of physicians thoroughly trained in endoscopic ligation of esophageal varices and anorectal hemorrhoids. A thorough understanding of the technical principles, clinical applications, and risks associated with this procedure is necessary before using this device.

INTENDED USE/INDICATIONS FOR USE

The Speedband Superview Super 7 Multiple Band Ligator is used for endoscopic ligation of esophageal varices and anorectal hemorrhoids.

Clinical Benefit Statement

Speedband Superview Super 7 aids in stopping or preventing hemorrhage of esophageal varices or anorectal hemorrhoids.

CONTRAINDICATIONS

The Speedband Superview Super 7 Multiple Band Ligator is contraindicated for patients with bleeding disorders, unless the bleeding disorder is first identified and treated appropriately, and is contraindicated for any other condition which would otherwise contraindicate gastrointestinal endoscopy. The Speedband Superview Super 7 Multiple Band Ligator is not intended for ligation of esophageal varices below the gastroesophageal junction.

WARNINGS

- Ensure the trip wire is secured; failing to do so could lead to deployment failures which could lead to a delay in treatment.
- Do not remove the shrink wrap until specified, premature removal of the shrink wrap could lead to ligator thread becoming dislodged and result in deployment issues which could lead to a delay in treatment.

PRECAUTIONS

- The Speedband Superview Super 7 Multiple Band Ligator should only be used by or under the supervision of physicians trained in variceal and hemorrhoidal band ligation. A thorough understanding of the technical principles, clinical applications and risks associated with variceal and hemorrhoidal band ligation is necessary before using this device.
- The materials used in this device have been found to be an irritant in animal studies.
- Any use of this device, other than those indicated in these instructions, is not recommended.

ADVERSE EVENTS

Complications associated with any use of the Speedband Superview Super 7 Multiple Band Ligator may include:

- Allergic Reaction
- Discomfort
- Dysphagia
- Embolism
- Gastrointestinal Symptoms
- Hemorrhage
- Infection
- Inflammation
- Obstruction
- Pain
- Perforation
- Stricture Formation
- Tissue Damage

- Treatment-related bleeding secondary to ulceration at banding sites
- Ulceration
- Vasovagal Symptoms

Additional complications associated with Anorectal Hemorrhoidal ligation may include:

- Incontinence
- Prolapse of Organ
- Urinary Retention

HOW SUPPLIED

Contents supplied NON-STERILE. Do not use if damaged and/or opened.

Device Details

Do not use if package is damaged or unintentionally opened before use.

Do not use if labeling is incomplete or illegible.

Handling and Storage

Store at 25 °C (77 °F); excursions permitted to 15 °C - 30 °C (59 °F - 86 °F). Do not use if the device has been exposed to environmental conditions outside of those ranges, as device performance may be impacted. Do not expose to organic solvents, ionizing radiation or ultraviolet light. Rotate inventory so that products are used prior to the expiration date on package label.

OPERATIONAL INSTRUCTIONS

The following subsections address the preparation, initial setup and the usability steps for ligating an esophageal varix:

Preparation

1. Read all instructions thoroughly.
2. Verify that the packaging is undamaged.

Precaution: A thorough understanding of the technical principles, clinical applications, and risks associated with endoscopic variceal and anorectal hemorrhoidal band ligation is necessary **before using this device.**

Precaution: Do not use any component if the package is damaged or opened.

Precaution: Do not use if components are broken or damaged.

Initial setup

1. Prior to initial setup, check for the following:
 - a. that the endoscope is straight;
 - b. that the spool of the Handle Unit is at its "starting position" (i.e. the arrow on the spool should align with the arrow on the bracket);
 - c. that the trip wire is connected to the Handle Unit and through the hole in the spool; and
 - d. that the scope fastener is attached to the endoscope handle.
2. Feed the trip wire through the biopsy valve and then through the working channel of the endoscope until the wire loop exits the distal end of the endoscope.
3. Firmly install the biopsy valve onto the Stem of the Handle Unit (Figure 3A, 3B and 3C).

Figure 3A. (Olympus™)

Figure 3B. (Pentax™)

Figure 3C. (Fujinon™)

Install the Handle Unit (with the biopsy valve on the Stem) into the biopsy opening of the endoscope. Make sure the scope fastener is next to the endoscope handle. Wrap the fastener

around the endoscope and secure the connection between the endoscope and the Handle Unit by securing the velcro scope fastener (Figure 3D).

Attach the pulling loop of the Ligating Unit to the wire loop connector on the distal end of the trip wire by threading the pulling loop through the trip wire loop, then threading the Ligating Unit through the pulling loop (Figure 4A and 4B).

4. Take up slack by pulling the proximal end of the trip wire on the Handle Unit.
The pulling loop of the Ligating Unit will advance into the working channel of the endoscope. Verify that the trip wire does not fall into the gap between the handle unit spool and the bracket (Figure 4C).
5. Stop pulling on the trip wire when the soft proximal end of the Ligating Unit cylinder contacts the distal tip of the endoscope.


Figure 3D.


Figure 4A.



Figure 4B.

[1] Spool [2] Gap [3] Bracket [4] Proximal end of trip wire

Figure 4C.

6. Line up the black stripe on the adapter with the working channel of the endoscope (Figure 5).

Warning: If the Black Stripe is not aligned with the working channel, the bands may fire improperly or not at all, and the field of vision may be impaired leading to a delay in treatment.

Secure the Ligating Unit to the endoscope by twisting or rocking the soft proximal adapter of the Ligating Unit while pushing it onto the distal tip of the endoscope. The Ligating Unit should be pushed until the tip of the endoscope is up against the rigid part of the Ligating Unit cylinder.

7. Pull gently on the trip wire to take up the rest of the slack. Stop applying tension when resistance occurs.

8. Secure the trip wire in the slot on the Handle Unit (Figure 6A).
9. Orient excess trip wire such that it cannot tangle up on the handle unit and endoscope (Figure 6B).
10. Carefully remove the shrink wrap by pulling the marked tab such that the ligator bands and threads are not disturbed (Figure 7A and 7B).

Precaution: Unless the shrink wrap is removed, the bands will not fire.

11. For irrigation, attach a syringe or the irrigation tube provided to the luer fitting on the irrigation valve. A sterile syringe is not required.

The Speedband Superview Super 7 Multiple Band Ligator is now ready for use.



Figure 5.

[1] Slot [2] Trip Wire

Figure 6A.


Figure 7A.


Figure 6B.


Figure 7B.

Method for ligating an esophageal varix

(See below for instructions on ligating anorectal hemorrhoids.)

Note: This device may impede the selective aspiration of blood and secretions.

1. Introduce the endoscope, with the Ligating Unit attached, into the alimentary tract.
2. Visualize the target varix.
3. Press the distal tip of the endoscope with the Ligating Unit against the varix.
4. Once complete contact is made around the varix, activate suction of the endoscope, drawing the lesion into the housing.
5. When the lesion is completely in the housing, a "red out" will be observed.
6. Rotate the handle clockwise 180° until you hear a "click", and feel the indent. This action deploys one band and ligates the lesion (Figure 8A and 8B).
7. Release suction.
8. To ligate another lesion, visualize the new target and repeat steps 3, 4, 5, 6 and 7.
9. If irrigation is desired, attach a syringe, directly or via the irrigation tube, to the luer fitting on the irrigation valve and infuse the desired volume. The valve is closed unless a syringe or the irrigation tube is fully inserted into the valve.

Figure 8A.

Figure 8B.

Additional Clinical Information for Esophageal Variceal Ligation

Caution: Ligating bands may be difficult to apply to small (Grade 1) varices.

Warning: Esophageal variceal ligation may not achieve complete fibrosis of inner wall of esophagus which may be necessary to prevent rebleeding or recurrence of varices.

Method for Ligating an Anorectal Hemorrhoid

(See above for instructions on ligating esophageal varices.)

Note: This device may impede selective aspiration of blood and secretions.

1. Introduce the endoscope, with the Ligating Unit attached, into the rectum.
2. Visualize the target hemorrhoid. Select the largest hemorrhoid first.

Warning : Hemorrhoidal ligation should only be performed on internal hemorrhoids above the dentate line. Treatment below (caudal to) this level may result in severe pain.

3. Press the distal tip of the endoscope with the Ligating Unit against the hemorrhoid.
4. Once complete contact is made around the hemorrhoid, activate suction of the endoscope, drawing the lesion into the housing.
5. When the lesion is completely in the housing, a "red out" will be observed.
6. Rotate the handle clockwise until you hear a "click", and feel the indent. This action deploys one band and ligates the lesion (Figure 9A and 9B).
7. Release suction.
8. To ligate another lesion, visualize a new target and repeat steps 3, 4, 5, 6 and 7.
9. If irrigation is desired, attach a syringe, directly or via the irrigation tube, to the luer fitting on the irrigation valve and infuse the desired volume. The valve is closed unless a syringe is fully inserted into the valve.


Figure 9A.


Figure 9B.

Removing Speedband Superview Super 7 Multiple Band Ligator from Endoscope

1. Remove the endoscope from the patient.
2. Disconnect the scope fastener. Remove the Handle Unit and the trip wire from the endoscope.
3. Remove the Ligating Unit from the endoscope.
4. If desired, another Speedband Superview Super 7 Multiple Band Ligator can be used on the patient for additional banding.

Disposal

To minimize the risk of infection or microbial hazards after use, dispose device and packaging as follows: After use, device may contain biohazardous substances. The device and packaging should be treated and disposed of as biohazardous waste or have them treated and disposed of in accordance with any applicable hospital, administrative, and/or local government regulations. Use of a biohazardous container with biological hazard symbol is recommended. Untreated biohazardous waste should not be disposed of in the municipal waste system.

Post-Procedure

Any serious incident that occurs in relation to this device should be reported to the manufacturer, and to the relevant local regulatory authority.

PATIENT COUNSELING INFORMATION

Brief the patient on any relevant post-procedure instructions, contraindications, warnings, precautions and/or adverse events found within this Instructions for Use (IFU), pertaining to the patient.

WARRANTY

For device warranty information, visit (www.bostonscientific.com/warranty).

The following is a trademark of Boston Scientific Corporation or its affiliates: Speedband Superview Super 7

All other trademarks are the property of their respective owners.

SYMBOL DEFINITIONS

Commonly used medical device symbols that appear on the labeling are defined at www.bostonscientific.com/SymbolsGlossary.

Table 1. Symbol Definitions

	Contents
	Maximum Scope Length
	Outer Diameter
	Store at 25 °C (77 °F); excursions permitted to 15-30 °C (59-86 °F).
	Includes Irrigation Tube

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© 2025 Boston Scientific Corporation or its affiliates. All rights reserved.
A seguinte é uma marca comercial da Boston Scientific Corporation ou respectivas afiliadas: Speedband Superview Super 7
Boston Scientific Limited