Ultraflex™
Tracheobronchial
Stent System

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2024-07

TABLE OF CONTENTS

REUSE WARNING	3
DEVICE DESCRIPTION	3
Figure 1. Ultraflex Tracheobronchial Stent System	3
Contents	4
Operating Principle	4
Materials	4
Uncovered Stents	4
Covered Stents	4
User Information	5
INTENDED USE/INDICATIONS FOR USE	5
Clinical Benefit Statement	5
Summary of Safety and Clinical Performance	5
CONTRAINDICATIONS	5
WARNINGS	5
ADVERSE EVENTS	6
HOW SUPPLIED	6
Device Details	6
Equipment Required	6
Handling and Storage	7
OPERATIONAL INSTRUCTIONS	7
Initial Preparation of Delivery System	7
Procedure	7
Ultraflex Tracheobronchial Stent System (Uncovered) – Distal Release and Proximal Release	7
Figure 2. Delivery Catheter Positioning	8
Figure 3A. Stent Deployment – Distal Release System	8
Figure 3B. Stent Deployment-Proximal Release System	9
Figure 4. Delivery System Removal	9
Ultraflex Tracheobronchial Stent System (Covered) – Distal Release	9
Figure 5. Delivery Catheter Positioning	10
Figure 6. Stent Deployment – Distal Release System	11
Figure 7. Delivery System Removal	11
Disposal	11
Post-Procedure	12
MRI SAFETY INFORMATION	12
PATIENT COUNSELING INFORMATION	12
Implantable Device Patient Information	13
Implant Card Instructions	13
WARRANTY	13
SYMBOL DEFINITIONS	13
Table 1. Symbol Definitions	13

Ultraflex™ Tracheobronchial Stent System

Rx ONLY


Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

REUSE WARNING

For single use only. Do not reuse, reprocess or re-sterilize. Reuse, reprocessing or re-sterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or re-sterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

DEVICE DESCRIPTION

The Ultraflex Tracheobronchial Stent Systems are permanently implanted expandable metal stents designed to serve as an intraluminal support to keep open the inner wall of the tracheobronchial tree. They consist of a flexible delivery catheter preloaded with an expandable metallic stent (Figure 1).

[1] Ultraflex Tracheobronchial Stent System (Uncovered) [1.1] Uncovered Stent [1.2] RO Markers Identify Expanded Stent Margins [1.3] Crocheted Suture Cord [1.4] Finger Ring
[2] Ultraflex Tracheobronchial Stent System (Covered) [2.1] Covered Segment [2.2] Uncovered Segment [2.3] Uncovered Segment [2.4] Outer RO Markers Identify Expanded Stent Margins [2.5] Crocheted Suture Cord [2.6] Inner RO Markers Identify Covered Segment Margins after Expansion [2.7] Finger Ring

Figure 1. Ultraflex Tracheobronchial Stent System

The stent is an open-ended cylindrical mesh constructed from a single strand of nitinol wire. The wire is configured into a series of circumferential interwoven loops, with the number of loops being dependent on the diameter of the stent. The stent is elongated and compressed onto a plastic delivery catheter. The stent is held onto the delivery catheter by a crocheted nylon suture wrapped around the stent. The delivery catheter has a flush taper tip at the distal end, and a round hub handle at the proximal end.

Both an uncovered and covered stent are available. The covered stent has a single layer of translucent silicone that covers the midsection of the stent.

Contents

One (1) Ultraflex Tracheobronchial Stent System

Operating Principle

Covered stents are available with a distal release system only. Uncovered stents are available with either a proximal or distal release system. The distal release system begins stent deployment from the lower (distal) end of the delivery catheter. The proximal release system begins stent deployment from the upper (proximal) end of the delivery catheter.

The radiopaque (RO) markers on the delivery system facilitate fluoroscopic placement.

The uncovered stent delivery system has two (2) RO markers. These RO markers indicate the estimated final position of the ends of the deployed stent (Figure 1).

The covered stent has four RO markers. The outer two (2) RO markers indicate the estimated final position of the ends of the deployed stent. The inner two (2) RO markers indicate the estimated final position of the margins of the deployed stent cover (Figure 1).

The delivery system accepts a 0.035 in (0.89 mm) or 0.038 in (0.97 mm) guidewire. The stent delivery system is passed over the guidewire into the tracheobronchial lumen. The stent is positioned appropriately using the RO markers for guidance under fluoroscopy and by bronchoscopic visualization of the stent.

The delivery system is a suture release design with a finger ring at the proximal end that enables removal of the nylon suture. When ready to be implanted, the stent is deployed by holding the handle hub in the palm of one hand and grasping the finger ring with the other hand. By retracting the finger ring the suture crochet knots are unraveled in a circular manner along the length of the stent, gradually deploying the stent. When using the proximal release system, stent deployment begins from the upper (proximal) end of the delivery system and continues to release toward the tip of the delivery system (distally) as the entire suture unravels. When using the distal release system, stent deployment begins from the lower (distal) end of delivery system and continues to release towards the operator (proximally) as the entire suture unravels.

After the stent is completely released, and the nylon suture has been completely removed, the delivery system catheter can be removed.

Materials

Materials and substances for which the patient can be exposed to, by the implantable portion of the medical device, are the following:

Uncovered Stents

Implantable Material	% Weight (0.108 g - 1.020 g)
Nitinol	65.7-97.1
Polyester	2.9-34.3

Covered Stents

Implantable Material	% Weight (0.420 g - 3.948 g)
Nitinol	31.0-32.4
Silicone	30.9-45.2

Implantable Material	% Weight (0.420 g - 3.948 g)
Silicone Adhesive	16.7-35.6
Polyester	1.1-7.1

Warning: This device contains nitinol, an alloy of nickel and titanium. Persons with allergic reactions to these metals may suffer an allergic reaction to this implant. Prior to implantation, patients should be counseled on the materials contained in the device, as well as potential for allergy/hypersensitivity to these materials.

User Information

Read the entire Instructions for Use in detail before using the Ultraflex Tracheobronchial Stent System. The Ultraflex Tracheobronchial Stent System should only be used by or under the supervision of physicians thoroughly trained in bronchoscopic stent placement technique. A thorough understanding of the technical principles, clinical applications, and risk associated with this procedure is necessary before using this device.

INTENDED USE/INDICATIONS FOR USE

The Ultraflex Tracheobronchial Stent System is indicated for use in the treatment of tracheobronchial strictures produced by malignant neoplasms.

Warning: The safety and effectiveness of this device for use in the vascular system has not been established and can result in serious harm and/or death.

Clinical Benefit Statement

The clinical benefit of the Ultraflex Tracheobronchial Stent System is the treatment of tracheobronchial strictures produced by malignant neoplasms.

Summary of Safety and Clinical Performance

For customers in the European Union, use the device name found in the labeling to search for the device’s Summary of Safety and Clinical Performance, which is available on the European database on medical devices (EUDAMED) website: (<https://ec.europa.eu/tools/eudamed>).

CONTRAINDICATIONS

The Ultraflex Tracheobronchial Stent System is contraindicated for:

- Concurrent fistula of the tracheobronchial tree, unless a covered stent is being used.
- Tracheobronchial obstruction with a luminal diameter which cannot be dilated to and maintained at least 4 mm, or preventing passage of either a rigid or flexible bronchoscope.
- Patients for whom bronchoscopic procedures are contraindicated.
- Any use other than those specifically outlined under Indications for Use.

WARNINGS

The Ultraflex Tracheobronchial Stent System should be used with caution and only after careful consideration in patients with: Compromised immune system; elevated bleeding times or coagulopathies; prior pneumonectomy; concurrent acute inflammation in the lumen, as this may potentiate granuloma formation and fibrosis; a tumor stricture adjacent to a major vessel, as this may potentiate fistula formation.

- Placement of an improperly sized stent can lead to a higher incidence of granulation tissue formation.

- Laser ablation to treat excessive granulation tissue post stent placement may result in loss of stent integrity and/or cause airway fire.
- Placement of any type of stent in the subglottic trachea is associated with a higher rate of granuloma formation and should be avoided.

ADVERSE EVENTS

The potential complications associated with tracheobronchial stent placement may include:

- Allergic reaction
- Aphonia
- Aspiration
- Bleeding/Hemorrhage
- Burn
- Bronchospasm
- Death
- Dysphagia
- Dyspnea
- Fistula formation
- Halitosis
- Hemoptysis
- Inadequate stent expansion
- Infection and sepsis
- Inflammation
- Obstructive atelectasis (even with a well-positioned stent)
- Pain/Discomfort
- Perforation
- Pneumoperitoneum
- Pneumothorax
- Respiratory distress
- Respiratory failure
- Restenosis
- Sedation or anesthesia complications
- Stent fracture
- Stent migration
- Stent misplacement
- Stent occlusion due to mucous, tissue ingrowth/overgrowth
- Tissue damage
- Tissue ingrowth OR overgrowth
- Tracheitis
- Ulceration and fistula formation

HOW SUPPLIED

The Ultraflex Tracheobronchial stent system is supplied sterile using an ethylene oxide (EO) process.

Device Details

Do not use if package is damaged or unintentionally opened before use.

Do not use if labeling is incomplete or illegible.

Equipment Required

- Bronchoscope, flexible or rigid.
- 0.038 in (0.97 mm) guidewire or a 0.035 in (0.89 mm) Jagwire Pulmonary Guidewire.
- Ultraflex Tracheobronchial Stent System of the appropriate stent length and diameter.
- Syringe for irrigation.
- Grasping forceps and/or retrieval snare

- Fluoroscopic capability for pre-stent placement and stent placement confirmation
- Dilatation balloon (as necessary).

Handling and Storage

This product has no special handling or storage requirements. Check device label for expiration date.

OPERATIONAL INSTRUCTIONS

Initial Preparation of Delivery System

- Carefully remove the delivery system from the protective packaging.
- Visually inspect the device for damage or defects

Warning: Visually inspect the system for any signs of damage. DO NOT USE if the system has any visible signs of damage. Failure to observe this warning may result in patient injury.

Procedure

Precaution: Read the entire Procedural Instructions for the specific stent system being employed:

Ultraflex Tracheobronchial Stent System (Uncovered) – Distal Release and Proximal Release

1. Locate Stricture and Dilate

Pass a bronchoscope into the airway and through the tracheobronchial stricture. If necessary, dilate the stricture using a balloon dilator to the minimum diameter of the preloaded stent delivery catheter diameter, or until a bronchoscope can be passed.

Warning: Physicians should use judgement based on experience in dilating airway strictures. Perforation, bleeding and/or migration could result from dilation of a tumor.

Warning: Placement of the Ultraflex Tracheobronchial Stent System is contraindicated in patients with strictures that cannot be dilated to at least 4 mm or cannot pass a bronchoscope in order to avoid patient injuries such as perforation, hemorrhage and/or obstruction.

2. Measure the Stricture Length and Luminal Diameter

This may be done by visual inspection on bronchoscopy. To determine the stricture length, measure the distance while pulling back on the bronchoscope. To determine the lumen diameter, estimate the diameter of the normal appearing tracheobronchial lumen proximal to the stricture. An open biopsy forceps may be used for a reference guide. Alternatively, the stricture length and luminal diameter may be measured by reviewing a recent CT Scan of the strictured tracheobronchial lumen.

3. Examine Stricture

Bronchoscopically visualize below the stricture, noting the distance to any bronchial branches. Examine the stricture location fluoroscopically. The stricture should be dilated to approximately 75 % of the normal lumen diameter.

Optional: Radiopaque skin markers may be used to identify the stricture margins.

4. Insert guidewire

Place a 0.038 in (0.97 mm) guidewire or a 0.035 in (0.89 mm) Jagwire Pulmonary Exchange Guidewire through the stricture.

Optional: The bronchoscope may be removed at this time.

5. Select Stent Size

Choose the stent length to completely bridge the tumor stricture plus an additional 10 mm to 20 mm to overlap normal

mucosa. At least 10 mm should overlap at each end of the stent. Choose the stent diameter to approximate the size of the normal proximal lumen.


Warning: If possible, avoid choosing of a stent which would cross bronchial side branches when placed in order to avoid respiratory complications.

6. Prepare the Ultraflex Tracheobronchial Stent System

Remove the stylet at the distal end of the delivery system. Lubricate the stent with water soluble lubricant. Pass the guidewire into the distal end of the delivery system.

7. Advance the Ultraflex Tracheobronchial Stent System Over the Guidewire into the Stricture

Reintubate the bronchoscope in the airway. Under bronchoscopic visualization, advance the Ultraflex Tracheobronchial Stent System over the guidewire into the stricture. Using fluoroscopy, visualize the delivery system radiopaque (RO) markers, and position the RO markers so that the stricture is centered between the RO markers (Figure 2).


[1] Center Stricture between RO Markers [2] Finger Ring

Figure 2. Delivery Catheter Positioning

8. Deploy Stent

A. Distal Release System

Begin deployment by holding the inner catheter stationary with one hand, and, using the other hand, grasp the finger ring attached to handle end of the stent delivery system. Pull the finger ring to start releasing the suture knots and to unravel the suture. For distal release systems, the suture unravels from the distal (lower) stent end first (Figure 3A). Monitor the stent release under fluoroscopy, keeping the identified stricture margins centered between the delivery system RO markers. If necessary, stop deployment and adjust the stent position.

[1] Begins release closest to distal end of catheter [2] Holding the finger ring, pull back to release suture cord loops and begin stent deployment [3] Continue to release suture cord loops to complete full stent deployment

Figure 3A. Stent Deployment – Distal Release System

For distal release systems, the inner catheter and stent can be pulled towards the operator (outward of the airway) to the correct position.


Warning: Do not attempt to reload a deployed or partially deployed stent in order to avoid damage to the stent/delivery system and/or patient injury.

B. Proximal Release System

Begin deployment by holding the inner catheter stationary with one hand, and, using the other hand, grasp the finger ring attached to handle end of the stent delivery system.

Pull the finger ring to start releasing the suture knots and to unravel the suture. For proximal release systems, the suture unravels from the proximal (upper) stent end first (Figure 3B). Monitor the stent release under fluoroscopy, keeping the identified stricture margins centered between the delivery system RO markers. If necessary, stop deployment and adjust the stent position. For proximal release systems, the inner catheter and stent can be pulled away from the operator (inward of the airway) to the correct position.

Warning: Do not attempt to reload a deployed or partially deployed stent in order to avoid damage to the stent/delivery system and/or patient injury.


[1] Begins release closest to proximal end of catheter [2] Holding the finger ring, pull back to release suture cord loops and begin stent deployment [3] Continue to release suture cord loops to complete full stent deployment

Figure 3B. Stent Deployment-Proximal Release System

9. Assess Deployed Stent and Remove Delivery System and Bronchoscope

Confirm Bronchoscopically and fluoroscopically that the stent has completely deployed from the catheter. Carefully remove the delivery system from within the expanded stent, watching not to dislodge the stent with the catheter tip (Figure 4). Remove bronchoscope and recover patient as indicated.

Optional: If necessary, complete stent expansion using balloon dilators up to the diameter of the stent.

[1] Carefully remove the delivery system from within the expanded stent

Figure 4. Delivery System Removal

Warning: The stent is not intended to be removed once it is properly positioned. However, if it becomes necessary to remove the stent immediately post-deployment, the stent may be removed using forceps with teeth or a retrieval snare.

Ultraflex Tracheobronchial Stent System (Covered) – Distal Release

1. Locate Stricture and Dilate

Pass a bronchoscope into the airway and through the tracheobronchial stricture. If necessary, dilate the stricture using a balloon dilator to the minimum diameter of the preloaded stent delivery catheter diameter, or until a bronchoscope can be passed.

Warning: Physicians should use judgement based on experience in dilating airway strictures. Perforation, bleeding and/or migration could result from dilation of a tumor.

Warning: Placement of the Ultraflex Tracheobronchial Stent System is contraindicated in patients with strictures that cannot be dilated to at least 4 mm or cannot pass a bronchoscope in order to avoid patient injuries such as perforation, hemorrhage and/or obstruction.

2. Measure the Stricture Length and Luminal Diameter

This may be done by visual inspection on bronchoscopy. To determine the stricture length, measure the distance while pulling back on the bronchoscope. To determine the lumen diameter, estimate the diameter of the normal appearing tracheobronchial lumen proximal to the stricture. An open biopsy forceps may be used for a reference guide. Alternatively, the stricture length and luminal diameter may be measured by reviewing a recent CT Scan of the strictured tracheobronchial lumen.

3. Examine Stricture

Bronchoscopically visualize below the stricture, noting the distance to any bronchial branches. Examine the stricture location fluoroscopically. The stricture should be dilated to approximately 75 % of the normal lumen diameter.

Optional: Radiopaque skin markers may be used to identify the stricture margins.

4. Insert guidewire

Place a 0.038 in (0.97 mm) guidewire or a 0.035 in (0.89 mm) Jagwire Pulmonary Exchange Guidewire through the stricture.

Optional: The bronchoscope may be removed at this time.

5. Select the Covered Stent Size

Choose the stent length to completely bridge the tumor stricture plus an additional 10 mm to 20 mm to overlap normal mucosa. At least 10 mm should overlap at each end of the stent. Choose the stent diameter to approximate the size of the normal proximal lumen.


Warning: If possible, avoid choosing of a stent which would cross bronchial side branches when placed in order to avoid respiratory complications.

6. Prepare the Ultraflex Tracheobronchial Stent System (Covered)

Remove the stylet at the distal end of the delivery system. Lubricate the stent with water soluble lubricant. Pass the guidewire into the distal end of the delivery system.

7. Advance the Ultraflex Tracheobronchial Stent System (Covered) Over the Guidewire into the Stricture

Reintubate the bronchoscope in the airway. Under bronchoscopic visualization, advance the Ultraflex Tracheobronchial Stent System (Covered) over the guidewire into the stricture. Using fluoroscopy, visualize the delivery system radiopaque (RO) markers, and position the inner RO markers so that the stricture is centered between them (Figure 5). These indicate the position of the covered portion of the expanded stent. The outer RO markers indicate the expected position of the ends of the expanded stent.

[1] Outer RO Markers Identify Expanded Stent Margins [2] Inner RO Markers Identify Covered Segment Margins After Expansion


Figure 5. Delivery Catheter Positioning

8. Deploy Stent

Distal Release System

Begin deployment by holding the inner catheter stationary with one hand, and, using the other hand, grasp the finger ring attached to handle end of the stent delivery system. Pull the

finger ring to start releasing the suture knots and to unravel the suture. For distal release systems, the suture unravels from the distal (lower) stent end first (Figure 6). Monitor the stent release under fluoroscopy, keeping the identified stricture margins centered between the delivery system RO markers. If necessary, stop deployment and adjust the stent position.

[1] Begins release closest to distal end of catheter [2] Holding the finger ring, pull back to release suture cord loops and begin stent deployment [3] Continue to release suture cord loops to complete full stent deployment

Figure 6. Stent Deployment – Distal Release System


For distal release systems, the inner catheter and stent can be pulled towards the operator (outwards of the airway) to the correct position.

Warning: Do not attempt to reload a deployed or partially deployed stent in order to avoid damage to the stent/delivery system and/or patient injury.

9. Assess Deployed Stent and Remove Delivery System and Bronchoscope

Confirm bronchoscopically and fluoroscopically that the stent has completely deployed from the catheter. Carefully remove the delivery system from within the expanded stent, watching not to dislodge the stent with the catheter tip (Figure 7). Remove bronchoscope and recover patient as indicated.

Optional: If necessary, complete stent expansion using balloon dilators up to the diameter of the stent.

[1] Carefully remove the delivery system from within the expanded stent

Figure 7. Delivery System Removal

Warning: The stent is not intended to be removed once it is properly positioned. However, if it becomes necessary to remove the stent immediately post-deployment, the stent may be removed using forceps with teeth or a retrieval snare.

Disposal

To minimize the risk of infection or microbial hazards after use, dispose device and packaging as follows:

After use, device may contain biohazardous substances. The device and packaging should be treated and disposed of as biohazardous waste or have them treated and disposed of in accordance with any applicable hospital, administrative, and/or local government regulations. Use of a biohazardous container with biological hazard symbol is recommended. Untreated biohazardous waste should not be disposed of in the municipal waste system.

Post-Procedure

Note: It may take 24 - 72 hours for the stent to fully expand.

- Perform routine radiography to assess the stent position and patency post-procedure as warranted.
- Exercise caution during tracheal or tracheostomy suctioning, because the suction catheter may become caught in the stent struts and may dislodge the stent.
- Recurrence of dyspnea may indicate that the stent has become occluded by mucous or migrated out of the stricture. Additional bronchoscopic assessment and treatment may be required.
- Tracheobronchial stent placements have been found to decrease some patients' ability to remove mucous secretion out of the lung.

Any serious incident that occurs in relation to this device should be reported to Boston Scientific and the relevant local regulatory authority.

MRI SAFETY INFORMATION

 MRI Safety Information A person with the Ultraflex Stent(s) may be safely scanned under the following conditions. Failure to follow these conditions may result in injury	
Device Name	Ultraflex Stent
Static Magnetic Field Strength (B ₀)	1.5T or 3.0T
Maximum Spatial Field Gradient	30 T/m (3,000 gauss/cm)
RF Excitation	Circularly Polarized (CP)
RF Transmit Coil Type	Cylindrical Whole-body Coil Cylindrical Head Coil
Operating Mode	Normal Operating Mode
Maximum Whole-Body SAR	2 W/kg (Normal Operating Mode)
Maximum Head SAR	3.2 W/kg (Normal Operating Mode)
Scan Duration/ Temperature Rise	Under scan conditions defined above, Ultraflex Stent is expected to have a temperature rise of less than 4 °C and can be used for 60 minutes of continuous RF
MR Image Artifact	Image artifact caused by device may extend approximately 13 mm from the perimeter and 7 mm from the end of the stent with a spin and gradient echo pulse sequence

If information about a specific parameter is not included, there are no conditions associated with that parameter.

PATIENT COUNSELING INFORMATION

Stent description including the materials contacting the patient, including the warning for the risk of nickel allergy.

Provide the patient with the completed implant card to carry and explain that the Boston Scientific website has additional patient information with the summary of stent safety and clinical performance.

Inform the patient to present the implant card it to their Healthcare professionals (doctors, dentist, technicians) including at MRI scans so they can take the necessary precautions. *Boston Scientific recommends that the patient register the MR conditions disclosed in this Instructions for Use with the MedicAlert Foundation (www.medicalert.org) or equivalent organization.

Inform the patient if any serious incident that occurs in relation to this device should be reported to the manufacturer, and to the relevant local regulatory authority.

Brief the patient on any relevant post-procedure instructions, contraindications, warnings, precautions and/or adverse events found within this Instructions for Use (IFU), pertaining to the patient.

*Provide the patient a follow-up care instruction with dietary information and patient positioning as outline in the Post-Procedure section of the IFU. It may take 24 - 72 hours for the stent to fully expand.

Implantable Device Patient Information

Advise the patient that additional information may be available to them on the Boston Scientific website (www.bostonscientific.com/patientlabeling).

Implant Card Instructions

- Apply the peel-off label from the product onto the supplied patient implant card.
- Fill in the implant date, patient name, healthcare institution and/or physician information.
- Give the completed card to the patient.

WARRANTY

For device warranty information, visit (www.bostonscientific.com/warranty).

The following are trademarks of Boston Scientific Corporation or its affiliates: Ultraflex, Jagwire.

All other trademarks are the property of their respective owners.

SYMBOL DEFINITIONS

Commonly used medical device symbols that appear on the labeling are defined at www.bostonscientific.com/SymbolsGlossary.

Table 1. Symbol Definitions

	Contents
	Recommended Guidewire

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or its affiliates. All rights reserved.
Boston Scientific (Australia) Pty Ltd
Boston Scientific Limited
© 2024 Boston Scientific Corporation