Zero Tip™

Airway Retrieval Basket

REF M00513200, M00513210

en	Instructions for Use	2
es	Instrucciones de uso	6
fr	Mode d'emploi	10
de	Gebrauchsanweisung	14
it	Istruzioni per l'uso	18
nl	Instructies voor gebruik	22
pt-EU	Instruções de Utilização	26

51753276-01

2024-07

Zero Tip™

Airway Retrieval Basket

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

REUSE WARNING

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

DEVICE DESCRIPTION

The Zero Tip Airway Retrieval Basket consists of a wire shaft (or cannula) housed in a flexible sheath with a multi-wire 360°, self-expanding basket at the distal end. A handle is attached proximally via a luer lock connector to open and close the basket. The Zero Tip Airway Retrieval Basket is designed to be used through working channel (2.0 mm or larger) of an endoscope, such as a flexible fiberoptic or rigid bronchoscope during the procedure.

Contents

The package includes:

One (1) Airway retrieval basket

One (1) Introducer

Operating Principle

The Zero Tip Airway Retrieval Basket accomplishes bronchoscopic removal of foreign bodies from the airway tree through the manipulation of the proximal handle of the device which opens and closes the basket. The basket located on the distal end of the device, features a low-profile, flexible, knotted-tip Nitinol wire design, which helps to grasp, manipulate and/or remove objects obstructing the airway tree.

Materials

Contains cobalt: CAS No. 7440-48-4; EN No. 231-158-0. Defined as a 1B carcinogen and reproductive toxicant according to the European Commission in a concentration above 0.1% weight by weight.

Note: Current scientific evidence supports that metal alloys containing cobalt used in medical devices, do not cause an increased risk of cancer or adverse reproductive effects.

User Information

Use of this device should be restricted to use by or under the supervision of physicians trained in bronchoscopic procedures. Care should be exercised to prevent perforation or trauma of the linings and associated tissue.

The Zero Tip Airway Retrieval Basket is for single use and has no other user serviceable parts.

- Do not attempt to repair or to alter the device.
- Do not resterilize the device.

- Any use of this device, other than those indicated in these instructions is not recommended.

INTENDED USE/INDICATIONS FOR USE

The Zero Tip Airway Retrieval Basket is intended to be used to endoscopically remove foreign bodies in the airway.

Clinical Benefit Statement

The intended clinical benefit of Zero Tip Airway Retrieval basket is for the removal of foreign bodies in the airway.

CONTRAINDICATIONS

- This procedure should not be attempted in any patient whose general medical condition and degree of respiratory failure would not allow that patient to tolerate bronchoscopy (rigid or fiberoptic) and/or the manipulation required to perform the procedure.
- Uncooperative patient.

WARNINGS

- Some objects may be too large to be removed endoscopically with a retrieval basket. Do not use the Zero Tip Airway Retrieval Basket if the object is too large to be removed endoscopically or held in the basket to avoid patient injury To avoid object impaction, fluoroscopy or x-ray should be used to determine the size of the object.
 - This device must not come in contact with any electrified instrument, as this may result in patient injury such as electric shock, burn, or tissue damage.
-

ADVERSE EVENTS

Adverse events which may result from the use of a basket in a bronchoscopic procedure include, but are not limited to:

- Airway obstruction
- Allergic reaction
- Bronchospasm
- Displacement of foreign body to a more unfavorable or peripheral location
- Edema (eg.-laryngeal, subglottic, tracheal, etc)
- Entrapment
- Evulsion
- Fracture of an organic foreign body
- Hemorrhage
- Laryngospasm
- Pneumomediastinum
- Pneumothorax
- Respiratory Distress
- Respiratory Failure

HOW SUPPLIED

Device is supplied sterile using an ethylene oxide (EO) process. The packaging and device should be inspected prior to use.

Device Details

Do not use if package is damaged or unintentionally opened before use.

Do not use if labeling is incomplete or illegible.

Handling and Storage

This product has no special handling or storage requirements. Check device label for expiration date.

OPERATIONAL INSTRUCTIONS

Precaution: Before using, inspect the blister pouch for any breach of the package to ensure a sterile product and inspect product for any damage. If seal has been broken or product is damaged DO NOT USE. Immediately return package and product to Boston Scientific.

Prior to use, inspect the basket assembly for any visible damage. Ensure that the handle is securely attached to the wire cable by testing the device as outlined below.

Thumb Slide Style Handle

Place handle in hand with thumb tab up. Push thumb tab forward to close basket and pull back to open the basket by advancing it out of the sheath. Test the basket by opening and closing it several times.

Precaution: Kinks in the sheath will hinder the mechanical operation of the basket.

Procedure

Insertion

The Zero Tip Airway Retrieval Basket devices are introduced into a bronchoscope in either of two ways; i.e., via an introducer (Fig. 1) to facilitate introduction into the scope or directly into the working channel port of the scope.


Figure 1. Introducer

1. First, insert Introducer into the nipple of working channel port of the endoscope.
2. Maintain the basket cage in the closed position, while introducing the (distal) working end of the device into the working channel of the endoscope.
3. Carefully advance basket device through the working channel. Using deliberate short strokes (2 to 3 centimeters), advance the sheath through the scope channel while watching for the sheath as it emerges from the distal scope channel.

Capture

1. Advance the closed basket alongside or beyond the object and open the basket by sliding the thumb tab.
2. Opening basket facilitates foreign body entrapment.
3. Slowly sweep the open basket across the face of the object to engage it into the wire cage.

Removal

Once the foreign body is captured within the basket, carefully maintain pressure on thumb tab in closed position to hold the object. The bronchoscope and foreign body should then be removed as a unit.

Warning: If resistance is encountered during advancement or withdrawal of device, DO NOT EXERT EXCESSIVE FORCE.

Precaution: Objects may be too large once captured in the basket to withdraw the basket fully into the sheath or the scope. Do not

exert excessive force to withdraw basketed object into sheath or into scope. Carefully maintain basket position to hold the object in place. The bronchoscope and foreign body are then removed as a unit.

Disposal

To minimize the risk of infection or microbial hazards after use, dispose device and packaging as follows:

After use, device may contain biohazardous substances. The device and packaging should be treated and disposed of as biohazardous waste or have them treated and disposed of in accordance with any applicable hospital, administrative, and/or local government regulations. Use of a biohazardous container with biological hazard symbol is recommended. Untreated biohazardous waste should not be disposed of in the municipal waste system.

Post-Procedure

Any serious incident that occurs in relation to this device should be reported to the manufacturer, and to the relevant local regulatory authority.

PATIENT COUNSELING INFORMATION

Brief the patient on any relevant post-procedure instructions, contraindications, warnings, precautions and/or adverse events found within this Instructions for Use (IFU), pertaining to the patient.

WARRANTY

For device warranty information, visit (www.bostonscientific.com/warranty).

The following are trademarks of Boston Scientific Corporation or its affiliates: Zero Tip.
All other trademarks are the property of their respective owners.

SYMBOL DEFINITIONS

Commonly used medical device symbols that appear on the labeling are defined at www.bostonscientific.com/SymbolsGlossary

Table 1. Symbol Definitions

	Contents		Includes Introducer
---	----------	---	---------------------

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or its affiliates. All rights reserved.
As seguintes marcas são marcas comerciais da Boston Scientific Corporation ou das respetivas afiliadas: Zero Tip.
Boston Scientific (Australia) Pty Ltd
Boston Scientific Limited
© 2024 Boston Scientific Corporation