



*Stone Cone*TM

Nitinol Urological Retrieval Coil

Boston
Scientific



Stone Cone™

Nitinol Urological Retrieval Coil

The Stone Cone Retrieval Coil
is designed for use with
intracorporeal lithotripters,
including the holmium laser.


The sheath is designed to resist
stray holmium laser energy.

Designed to automatically
release stones when necessary.

The coil is designed to
prevent stone migration
during lithotripsy procedures.

Unique design promotes
sweeping multiple stone
fragments in a single pass.

Quick release handle
facilitates scope exchanges.



Stone Cone Nitinol Retrieval Coil Ordering Information

Product Code	Description
M0063903100	Stone Cone Nitinol Urological Retrieval Coil 3.0Fr, 7mm Coil (Each)
M0063903200	Stone Cone Nitinol Urological Retrieval Coil 3.0Fr, 10mm Coil (Each)

CAUTION: Federal Law (USA) restricts this device to sale by or on the order of a physician.
Refer to package insert provided with the product for complete Instructions for Use, Contraindications for Use,
Potential Adverse Effects, Warnings and Precautions prior to using this product.
Stone Cone is a trademark of Med Source Technologies.

Boston Scientific

Delivering what's next.™

Boston Scientific Corporation
One Boston Scientific Place
Natick, MA 01760-1537
www.bostonscientific.com/urology

Ordering Information
1.888.272.1001

© 2006 Boston Scientific Corporation
or its affiliates. All rights reserved.

BVU940 5M 12/06-12/08